



# Newsletters

Advancing the science of physical testing

Head Office : 401, Tirupati Udyog, I. B. Patel road, Off. Western Express highway, Goregaon (East), Mumbai - 400 063, India.

Tel : +91 - 22 - 4041 3131 • E-mail : [sales@electrolabgroup.com](mailto:sales@electrolabgroup.com)

Factory : EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.


Tel : +91 - 22 - 4161 3122 • Website : [www.electrolabgroup.com](http://www.electrolabgroup.com)



## Frequent Waterbath Maintenance? ⇒ Get ClearView™

### US FDA form 483 excerpt (2017)

...Your analyst recorded the water in water bath for \*\*\* dissolution tester as clear. I observed the water in the water bath was not clear as required in your written procedure for operation and calibration of the \*\*\* dissolution tester. I also observed unidentified whitish \*\*\* floating in the water bath.



Waterbath Temperature



Circulation




Promotes growth of microorganisms in waterbath



Demands frequent waterbath cleaning


### General practices of cleaning of water bath



Cleaning Out of Place



Mechanical Qualification could be affected



Cleaning In Place



Breakdown due to microorganism deposition in tubings and pump

### Alternative methods to maintain clean a water bath for a longer period of time



Chemical Agents

**Limited Effectiveness**  
When used independently




Ultraviolet Sterilization

**Broad Effectiveness**

## ELECTROLAB Solution

### ClearView™ (ECV-08): Retrofittable to most brands of dissolution testers.



Within weeks of routine use



After several months of integration with ClearView™

Description	Part No.
ECV-08 ClearView M.S (230V)	1523600
ECV-08 Lamp	1536A00001

### Note:

- To enhance the efficiency, use ClearView™ in combination with deconex® 12 BASIC.
- Efficiency of ClearView™ depends on the source water microorganism load.

### Advantages:



Unaffected MQ



Enhances Visual Observation



Reduces Breakdown



Reduces Consumables



Reduces Water Consumption



Reduces Labour Time

### ELECTROLAB Cares

- Use of de-mineralized water in water bath is recommended to enhance the efficiency of ClearView™.
- To be more effective, ClearView™ should be used on water that is clean from the beginning as dead microorganisms already present in the water cannot be removed.
- To achieve optimum kill rate, ClearView™ should be kept on continuously.
- Replace the UV lamp after every 12 months of continuous use.
- Ensure that there is no leak or air bubble entrapped in the waterbath tubings.

To know more, please contact our *Specialist* on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : [sales@electrolabgroup.com](mailto:sales@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Rotating Bottle Apparatus

### Introduction

- In 1970s, Professor Beckett and many workers used the rotating bottle apparatus in order to assess the dissolution profiles of extended-release products especially those in the pellet form.
- At present, the rotating bottle apparatus is used in release-rate testing of many dosage forms.
- In this apparatus, bottles are clipped on a rotating shaft which is immersed in a water bath and the shaft is rotated at a desired speed.

**What does USP state about rotating bottle apparatus?**  
**USP 39-NF 34 <1092> The Dissolution Procedure: Development and Validation states:**

A rotating bottle or dialysis tubes may have utility for microspheres and implants...

### Advantages

Rotating bottle apparatus presents important advantages over the basket and paddle apparatuses, especially with regards to hydrodynamics and the possibility of using pH gradient. The advantages of rotating bottle apparatus are as follows:

- Provides a sound hydrodynamic system i.e. the dosage form moves freely through the dissolution medium as the bottle is rotated.
- This is in contrast to the movement of media relative to the dosage form in USP apparatus 1 and 2, wherein various portions of the bulk medium move at different rates.
- Facilitates pH profiling.
- Offers higher reproducibility.
- Small volume dissolution environment.
- Prevents evaporative loss.


### Applications

- Tablets and capsules
- Microspheres
- Implants
- Coated pellets
- Nanoparticulates
- Stents
- Suppository
- Liposomes
- Milk-clotting activity

### Methodology

1. Add the dosage form and extracting fluid to the standard test bottles.
2. Close the caps of the test bottles and snap the bottles on the rack.
3. Carefully place the rack into the water bath and close the lid. Turn on the rotating bottle apparatus.
4. At appropriate sampling intervals, the apparatus stops automatically. Remove the rack assembly and withdraw the desired volume of sample.
5. Replenish the volume of sample withdrawn with fresh media, if required.
6. In case of media changeover, decant the previous media using a suitable tool and replace it with new media.

Comparative studies on drug release testing from microparticulate drug delivery systems using different dissolution techniques.




(Ref: Conti B. et al., 1995.)

**Note:** FDA recommends the use of rotating bottle apparatus for drug release testing of Naltrexone (extended release suspension) and Risperidone (injectable).

## ELECTROLAB Solution

### Rotating Bottle Apparatus (ERB-Wu)



**Note:** Test bottle capacity - 10 mL, 100 mL, 200 mL and 250 mL

Integrity Operating System	Multi-Level Security	10 User Programmable Roles	100 Users	Excursion Report	Non-Volatile Memory	User Defined Report Templates	Report Printout	Technology

### ELECTROLAB Cares

- Ensure the cap of the test bottle is closed tightly to prevent leakage.
- Ensure the bottles are secured firmly onto the rack assembly.
- Weigh the dosage form accurately before loading it into the test bottles.
- As a safety precaution, wear gloves to remove the rack assembly from the water bath.
- Clean the water bath regularly using mild detergent like deconex 12 BASIC to prevent algal growth.

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
 Off Western Express Highway,  
 Goregaon (E) Mumbai - 400 063, India.  
 ☎ : 91 22 4041 3131  
 ✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

**DISCLAIMER :** The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Virtually Eliminate Waste Volume in Autosampling

**USP 38-NF 33 <1092>**

**The Dissolution Procedure: Development and Validation, states ...**

- Comparison of manual and automated procedures should be performed to evaluate the interchangeability of the procedures.
- Autosampling is a useful alternative to manual sampling, especially if the test includes several time points.

**USP 38-NF 33 <1092>**

**The Dissolution Procedure: Development and Validation, states ...**

- Automated methods offer opportunities for increased precision and reproducibility; however, bias may be introduced.
- A recirculated sampling alignment can be operated either by "Discharging the tubing contents into the vessel after each sampling or by allowing the tubing to remain filled with solution in the intervals between sampling points. In the latter case, the dead volume and carryover effects are important considerations".


**Advantages of Auto Sampling**



## ELECTROLAB Solution

**Loop-Back™ technology in Dx sample collector (ESC-8Dx/12Dx)**

### Schematic of Loop-Back™ Technology



### Example

Media Consumption	Manual	With Loop-Back™	Without Loop-Back™
Number of samples	5	5	5
Sample volume	10	10	10
Rinse volume for each sample	0	Rinse volume - Loop-Back volume (3 ml - ~ 3 ml) = ~ 0 ml	3
Total volume withdrawn from vessel	50	~ 50	~ 65 ml

### ELECTROLAB Cares

- Electrolab offers complimentary phone support to enable smooth functioning of Loop-Back™ technology
- With your existing Dx sample collector the Loop-Back™ feature can be activated in the dissolution tester by menu: Menu >> Configure >> Sampling >> Auto >> Syringe >> 10 ml (Syringe Capacity) >> Rinse by Loopback.

### Contact Us

401 Tirupati Udyog, I. B. Patel Road, Off Western Express Highway, Goregaon (E), Mumbai - 400 063, India.  
☎ : 91 22 4041 3131 ✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)



## 150 ml dissolution apparatus now compendial in United States Pharmacopeia (USP)



### What

- When the conditions present at the actual site of application are not suitable for the conventional dissolution method to provide relevant data, a smaller dissolution vessel is preferred



### Why

- Small volume dissolution is a useful technique when dosage contains very small amount of active ingredient, or for ER tablets where small amount of drug is released and may be difficult to detect using HPLC or UV-Vis spectrophotometer
- Smaller media volume increases the concentration of the sample and allows proper detection



### How


- Smaller dissolution vessels with special baskets and paddles have been designed to overcome these challenges
- The lowest volume that can be used in a 150 mL vessel is about 40 mL

### What did USP state about small volume dissolution?


- Till 2014, small volume dissolution was not compendial in USP, although; a 250 mL vessel was detailed in the Chinese Pharmacopoeia, but was not widely used
- As the dimensions for the smaller dissolution vessels, baskets and paddles were not specified in the pharmacopeia, specification variations between different apparatus manufacturers were observed

### Addition in USP about small volume dissolution

- In 2014, a small volume dissolution apparatus using 150 mL or 200 mL vessel has been discussed in USP forum - Chapter 40 (6) In-Process Revision: Doxercalciferol Capsules
- The development and regulatory submission of the vessel and paddle assembly for Doxercalciferol Capsules was a result of joint efforts by Cipla Ltd. and Electrolab.



### Dimensions of mini vessel and paddle



\*All measuring unit in mm

## ELECTROLAB Solutions

### 150 mL Dissolution Apparatus vessel and paddle

#### Features :

- Compatible with current USP 1 and USP 2 dissolution testers
- Easy Snap-fit™ shaft locking mechanism for positive engagement and wobble free operation
- Authenticity of vessels and paddles with ELECTROLAB hologram and laser marking
- Clear and amber colored vessels available

#### Accessories available:

- Standard:** Composite top plate, 150 mL glass vessels, jarlocks and jarlids
- Optional:** Validation kit - depth gauge, centering device, dial gauge etc

### 150 mL Conversion kit for 8 Station Dissolution Tester



Part no. 0104A00303



Conversion Kit for 14 Station Dissolution Tester also available

## ELECTROLAB Cares

- Clean the vessels with mild detergent /soap or cleansing agents
- Specially designed, washable and stackable basket & paddle tray



### Contact Us

401 Tirupati Udyog, I. B. Patel Road, Off Western Express Highway, Goregaon (E), Mumbai - 400 063, India.  
 ☎ +91 22 4041 3131 ✉ [sales@electrolabgroup.com](mailto:sales@electrolabgroup.com)

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Effective Deaeration → Reduces Variability

Recommend your friend or colleague

### What

Removal of dissolved air/gases from the dissolution media

### Why

Non - deaerated media may lead to dissolution failures as dissolved air/gases (bubbles) tend to change:

- Fluid flow characteristics
- pH of an unbuffered solution
- Nature of the active ingredient
- Dispersion characteristics of particles and aggregates
- Wetting pattern and cause the dosage form to become buoyant

< 6 ppm is USP recommended Dissolved Oxygen (DO) level

Effective deaeration reduces variability in results by :


- 52.3% of mean percent dissolved
- 11.8% of standard deviation

(Ref : Eaton J. et al., 2007)

- Different deaeration procedures are permissible, but they should be demonstrated to be equivalent to the USP procedure. Verification can be done using a dissolved gas meter.


(Ref : Dissolution toolkit version 2.0, 2010)

Variability is significantly higher when DO level is > 6 ppm :



(Ref : Nithyanandan P. et al., 2006)

Smaller batches of deaeration may lead to greater variability in DO level



## ELECTROLAB Solutions

### Dissolution Media Preparator (EMP 21-DO)

21 L

Batch Capacity

Rugged-SS Construction

Virtually 0

Dead Volume

Dearates Surfactants like SLS

+0.5°C

Temperature Accuracy

Volume Change Increment


10 Reports

Memory

<4 ppm

DO Level

Single batch deaeration reduces process time and ensures inter-vessel consistency in DO level



\*Data on file

Media Preparator with integrated DO probe

#### ELECTROLAB Media Preparator Main Cycle Report

Set Vacuum level : 500mbar Max vacuum level attained: 503mbar

Pilot Reading: YES Oxy.Level of Indv.Vessel: YES

Oxygen Level of Pilot Reading: 3.80 ppm

S.R.No.	Date	Time (Hrs:min)	Dispense Volume (ml)	Volume Stock (Ltrs)	Oxygen Level (ppm)
1]	28/11/2013	17:43	0500.00	10.64	3.69
2]	28/11/2013	17:44	0500.00	10.14	3.99
3]	28/11/2013	17:45	0500.00	09.64	3.82
4]	28/11/2013	17:46	0500.00	09.14	3.75
5]	28/11/2013	17:47	0500.00	08.64	3.89

Report with individual vessel DO values

### ELECTROLAB Cares

- EMP-21 DO can hold the media under vacuum, without electrical power for several hours. It is recommended that the dissolution test is started promptly after deaeration.
- Media for all the vessels in the dissolution tester should be deaerated in a single batch to ensure consistency.
- Cleaning cycle should be run after deaeration of every media.
- All tubes and pipes should be kept away from sharp edges.

To know more, please contact our Specialist on support@electrolabgroup.com

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com


DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click unsubscribe@electrolabgroup.com

## Disintegration Test - Procedure and Specifications


Recommend your friend or colleague

### USP procedure for Disintegration Test



Ref : USP Chapter <701>

### Changed USP specifications for Disintegration Test



The USP requirements for volume of immersion medium should be followed correctly as the volume plays the most crucial part to obtain correct and reproducible test results

Ref : Schmid K. et al. Influence of the Changed USP Specifications on Disintegration Test Performance. Dissolution Technologies;6-10

**Note :** Beaker should be USP compliant (ID = 97-115 mm) as the measurement of the immersion medium volume is directly related to the size of the beaker

## ELECTROLAB Solutions


### 2 Station Semi Automatic Park Out (EDI-2SA)



- Complies with USP, IP, BP and Ph. Eur. specifications
- Bathless heating saves: • Time • Electricity • Water
- Quick heating of media attains media temperature 30% faster as compared to disintegration tester with bath
- Offers disintegration time registration for individual tablets
- Green illumination from bottom for active arm offers better visibility of disintegration process
- Ergonomically designed and magnetically coupled basket for easy loading




Narrower tolerance for ELECTROLAB beakers



Test Report

### ELECTROLAB Cares

- Thorough cleaning of the basket is recommended for accurate and reliable results
- The entire basket assembly can be disassembled without using any special tool
- During the cleaning process, the rubber retainer will prevent the tubes from falling
- The use of online visual verification of stroke amplitude may assure the user of compliance with regulatory standards



### Contact Us

401 Tirupati Udyog,  
I. B. Patel Road,  
Off Western Express  
Highway, Goregaon (E),  
Mumbai - 400 063, India.  
+ 91 22 4041 3131  
sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Machine Verification replaces human verification in Blue Dye Leak Test method

What	Why	How
<ul style="list-style-type: none"> <li>A 'Leak Test' is a method that detects the presence of a defect in a package; the defect is capable of permitting loss of product contents or critical head space gases, and/or capable of permitting entry of liquids, reactive gases, microorganisms or nonviable particulate (Ref. USP 40 (5) &lt;1207&gt;).</li> </ul>	<ul style="list-style-type: none"> <li>Any sort of physical breach in the packaging structure, such as weak seals, capillaries, pin holes or micro cracks can result in product integrity and safety being compromised. Hence, evaluating the integrity of packaging in the pharmaceutical and food industries is essential.</li> <li>It is not only critical in terms of preserving the product shelf-life and efficacy, but is also a FDA regulatory requirement as specified in cGMP 21 CFR 211.166 for finished pharmaceuticals.</li> </ul>	<ul style="list-style-type: none"> <li>A number of methods exist to detect leaks and various media are used as tracing agents depending upon the method used.</li> <li>A vast majority of pharmaceutical companies are exclusively using the dye ingress test (Tracer Liquid Test) to assess the integrity of the blister packs, with varying degrees of success.</li> </ul>

### USP (USP 40 (5) <1207.2>)

"Tracer liquid test requires the use of either analytical detection through instrumentation or visual inspection aids" that may reduce significant human dependency.  
 "Analytical detection offers the advantage of minimizing the human error that is inherent in visual discernment of low dye concentrations".


#### Tracer Liquid Test

This is the simplest and the oldest method used for integrity testing of blister packs. In this method, the package is placed in a bath of methylene blue dye within a desiccator and a set vacuum is drawn on the package. This technique involves the following phases:

- Evacuation:** During this phase, the instrument attempts to draw air out of the package cavity.
- Retention Time (Vacuum Hold Time):** Following evacuation, the submerged test samples are subjected to vacuum at a predetermined vacuum level for a predetermined time.
- Venting:** During venting, vacuum is released from the test chamber and the next phase begins.
- Penetration:** Following vacuum release, test samples remain submerged in the dye for a predetermined time.

If the package cavity has a leak, it will have a negative pressure during evacuation, drawing the methylene blue dye inside the cavity. The lab operator inspects the package for any blue dye ingress following the completion of the test.

#### Graphical Representation of Tracer Liquid Test Procedure



#### Advantages

- Well established method
- Quick and economical method
- Method provides precise leak detection data

#### Disadvantages

- Significant human dependency may lead to errors
- Lack of verifiable documentation
- Not eco-friendly disposal of dye

## ELECTROLAB Solution

### ELT-201WP (Patented Design)

Electrolab Leak Tester (ELT-201WP) is the **first & only** leak tester which replaces human subjectivity by machine verification in tracer liquid test. Leak detection is automatically verified by weight, reducing human error and providing traceability to the results.



Compliance



Analytical Verification



Programmable Vacuum, Retention & Penetration Time



Multiple Packages

CE

Patented



Technology



Security for Data Integrity



Storage



Printout of Report

### ELECTROLAB Cares

- As leak detection is verified by weight, test can be carried out using water rather than methylene blue dye, thus helping in waste management
- Use lint free cloth to wipe the packages before recording the post-test readings
- Record the readings into the system only when the reading on the weighing balance is stable

### Contact Us

401 Tirupati Udyog, I. B. Patel Road, Off Western Express Highway, Goregaon (E), Mumbai - 400 063, India.  
 Phone : 91 22 4041 3131 Email : sales@electrolabgroup.com

To know more, please contact our **Specialist** on support@electrolabgroup.com

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice


To unsubscribe click unsubscribe@electrolabgroup.com



## STARLITE - 12

### Thermal Paper with 12 years life

- GMP regulations from PIC/S, USFDA and EU regulations require maintenance of accurate back-up files of input data for a minimum period of 5 years which should be secured from alteration, loss or inadvertent erasure.
- If a record is poorly documented, then the manufacturing or QA/QC of a product can be negatively impacted, potentially reducing its reliability.




Impact Printing v/s Thermal Printing		
Features	Impact Printing	Thermal Printing
1. Printing Details		
Print head	Moving	Stationary
Printing speed	Low (8 Lines/sec)	High (52 Lines/sec)
Noise	More (60 dB - 62 dB)	Less (48 dB)
Printing ability	Unable to print complex typographic characters as well as graphs	Prints complex typographic characters as well as graphs efficiently
2. Operational Cost		
Ribbons and cartridges	Yes : Requires ink for printing	No : Inkless printing
Cartridge change	More frequent	Not applicable
Printer outage and repairs	Frequent due to more moving parts	Less
Maintenance	High	Low
3. Type Of Paper		
Non-thermal paper	Fibers of paper dislodge from paper to printer head, eventually causing breakdown of printer head	-
Thermal paper with top coat	No dislodging of fibers on printer head	Extended life

## ELECTROLAB Solution

### STARLITE-12


#### Commonly available thermal paper v/s STARLITE-12



Use STARLITE-12 in weighing scale printers (impact printers) for increased life of printer head.



STARLITE-12 Thermal paper roll (Set of 24)  
Part No. - 0905A00007



### ELECTROLAB Cares

- Avoid exposure to direct sunlight and hot or humid environment for long-term durability
- Avoid prolonged contact with plasticizers (PVC) and solvents (acetone, alcohols etc)
- Avoid rough handling of thermal paper. Friction, scratching, or pressure from fingernails, paper clips and metal objects can cause impressions to develop, obstructing the printed image

#### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)



## Sampling Zone

Recommend your friend or colleague

The USP defines the sampling point as "A zone half way between the top of the media and top of the paddle or basket not less than 1 cm from the vessel wall". Ref: USP <711>

### Advantages of sampling

- Reproducible results
- Higher accuracy in results
- Unattended long duration test with multiple samples
- Reduces cost per test
- Efficient sampling & cleaning system


### Sampling Time

- As per USP,  $\pm 2\%$  variation is allowed in sampling. Ref: USP <711>


### Sampling Volume

- Within  $\pm 1\%$  (Data on file)


Note:  Sampling position




PADDLE WITH 1000ml MEDIUM



PADDLE WITH 500ml MEDIUM



BASKET WITH 1000ml MEDIUM



BASKET WITH 500ml MEDIUM

## ELECTROLAB Solutions

Feature	Manual sampling	Sampling with peristaltic pump	Sampling with syringe pump
Height adjustment	Fixed	Programmable	Programmable
Accuracy	User dependent	Low	High
Time taken in routine calibration	More	Less	—
Contact material	User dependent	Marprene and Tygon	Inert syringes, glass and teflon material
Split rinse	User dependent	Yes	Yes
Dead volume (wastage)	User dependent	More	less
Coarse filtration	Very good	Good	Good
Fine inline filtration	2 ways (cannula and disc filter)	1 way (cannula filter)	2 ways (cannula and disc filter)
Dilution	User dependent	*	**
Mixing	User dependent	*	**
Carryover removal	User dependent	Recirculation method	Recirculation and split rinse
Cleaning	User dependent	-	Auto
Validation	User dependent	Required	Not required

## ELECTROLAB Cares

For the long life of the tubes always check tubings for in automated sampling for -

- Leakage
- Cracks
- Reduced elasticity

Pull the release lever for removing cassettes for stand by situation

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our *Specialist* on support@electrolabgroup.com

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click unsubscribe@electrolabgroup.com

## Topical & Transdermal Drug Delivery

Recommend your friend or colleague

Topical / Transdermal drug delivery system works by diffusion, which is self-contained, discrete dosage form that, when applied to intact skin, is designed to deliver the drug(s) through the skin to the systemic circulation

[Ref : Pharmacopeial forum Vol. 35 (3)]

### Commonly used dosage forms include:

- ▶ Ointments
- ▶ Emulsions
- ▶ Creams
- ▶ Lotions
- ▶ Gels
- ▶ Foams
- ▶ Pastes
- ▶ Powders
- ▶ Sprays
- ▶ Collodion
- ▶ Topical suspensions
- ▶ Ophthalmic ointments
- ▶ Other novel drug delivery systems

### Franz Diffusion Cell

The Franz Diffusion Cell is a static diffusion system used for characterizing drug permeation through a skin model

#### Different membrane can be used :

##### Natural

Stratum corneum plus epidermis , Stratus corneum, Mouse, Hairless mouse, Rat, Guinea pig, Rabbit, Monkey, Pig, Egg -shell membrane, etc

##### Synthetic

Cellulose media, Synthetic polymers, etc

## ELECTROLAB Solutions


### Diffusion Cell Apparatus ( EDC - 07 )

7  
Station

Size of Cell  
7.5 ml  
12.5 ml  
20 ml



Cell Centering Feature for Efficient Stirring



Diffusion Cell

## ELECTROLAB Cares

- To increase the service life, change the water regularly in circulation system of Diffusion Cell Apparatus

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Investigation and evaluation of OOS & OOT

Recommend your friend or colleague

FDA Part 211, subparts I (Laboratory Controls) and J (Records and Reports) recommends to investigate the cause of OOS and OOT to avoid the batch rejection. A meaningful OOS & OOT investigation should be: thorough, timely, unbiased, well-documented and scientifically defensible.

### Some of the obvious errors that could invalidate results

- Spilling of a sample
- Incomplete transfer of sample
- Obvious pipetting errors
- Improper sampling
- Dropping more than one tablet into dissolution vessel

### Practices used in the laboratory phase of an OOS & OOT investigation

- Retesting a portion of the original sample
- Re-sampling test data
- Using outlier testing
- Testing a specimen from the collection of a new Sample from the batch

If the investigation results of OOS & OOT are unsatisfactory, the batch is suspect and must be rejected or held pending for further investigation.

## ELECTROLAB Solutions

### iDisso™ - o6

### Observation and recording of the drug dissolution process



Patented



- Retrofitting system i.e. compatible with most ELECTROLAB Dissolution Testers made over 10 years
- Perfect tool to help in evaluation of OOS & OOT

- Accessory : BlackOut™ with photo -protective monochromatic lights for light sensitive drugs
- Over 1000 hours continuous recording can be stored

To view the video of iDisso™ - o6 [click here](#)

## ELECTROLAB Cares

- To save energy & environment, iDisso™ - o6 has eco-friendly and low power consumption lights

iDisso™ - o6 : 6 W  
Other lights : 40 W

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Commonly Associated Problems with Particle Analysis

Recommend your friend or colleague

### Problems, Probable Causes & Solutions

Problems	Probable Causes	Solutions
Sieving inaccuracy	Resulting from load and power frequency differences, crevices and label rivets which trap sample powder	Sieves of corrosion resistant material to prevent contamination of sample
Sieving inefficiency	Mesh size change due to attached material /particle within the mesh	3D sieving increases efficiency
Change of pore size	Due to using harsh cleaning agent / brush	Use gentle cleaning agent and regular verification of mesh size
Fine particles containment during sieving	Improper fitting of sieves result in the release of fine particles in environment	Snug fitting sieves to prevent release of fine particles in environment
Retention of excess powder on sieves	Excessive mesh binding and blockage, adhesive properties of powder, high ambient or powder humidity	Proper vibration decreases the powder retention on the sieve
Traceability	Label stuck or riveted to sieve body comes apart	Individualized laser engraving allows full traceability throughout life of sieve

## ELECTROLAB Solutions

### Electromagnetic Sieve Shaker (EMS-8 Plus)



Particle size analysis

Programmable amplitude



Larger LCD display

Programmable sieving duration



Low noise

Low maintenance



3D Sifting Motion

- Intermittent and continuous sifting motion with programmable amplitude (0.5 mm to 2 mm) in increments of 0.1 mm
- Continuous amplitude measurement and correction to give reproducible results independent of load
- Programmable sieving duration from 1 min to 99 min
- Suitable for dry and wet sieving



SieveEasy™  
Particle Size Analyzer Software

### Application industries



& Many more....

### ELECTROLAB Cares

- Clean sieves by sonicating them in water or a special detergent (Eg. deconex 12) to ensure accurate and reproducible sieving results
- It is advised not to run the machine without proper locking and tightening both the locking knobs, this may result in an abnormal sound and can damage the machine

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : [sales@electrolabgroup.com](mailto:sales@electrolabgroup.com)

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Filtration in Dissolution Testing

Filtration stops the dissolution process, removes undissolved drug, and allows for sample analysis. The downstream analytical evaluation technique – UV spectrophotometer, HPLC or UPLC – following dissolution usually warrants filtration before analysis because particulate impurities can affect the analysis.

Recommend your friend or colleague

### Commonly Used Types of Syringe Filters


Frits	Syringe filters / Disc type filters
Pore size: 10 $\mu$ – 15 $\mu$	Pore size: 0.22 $\mu$ – 1 $\mu$
<ul style="list-style-type: none"> <li>Stainless steel (Reusable after sonication)</li> <li>Disposable (Single use)</li> </ul>	<ul style="list-style-type: none"> <li>PTFE</li> <li>PVDF</li> <li>Nylon 66</li> <li>PES</li> </ul>

### Properties of Syringe Filters

Features	PTFE	PVDF	Nylon 66	PES
Chemical compatibility	★★★★	★★★	★★★	★★
Dual direction operation	✓	✓	✓	Limited
Cost	₹ ₹ ₹	₹ ₹	₹ ₹	₹

### Filter Pore Sizes Used in Various Analysis

Pore size	0.22 $\mu$	0.45 $\mu$	1 $\mu$	10 $\mu$ – 15 $\mu$
UV				
HPLC				
UPLC				



Get 25 % more filtration area with larger filters (33 mm Dia)

\* Data on file

## ELECTROLAB Solutions

### Most Versatile Solution in the Market with Multiple Filtration Options



- Powerful Syringe Pump enables the filtration process up to 2 or 3 stages to enable sample to be used directly in HPLC / UPLC without further filtration

- Benefits of using Manual Syringe Filters:
  - Ease of method transfer
  - Reduced operating cost
  - Filter re-validation studies not required

ELECTROLAB Pre pump filter manifold & Post pump filter manifold are compatible with 25 mm and 33 mm manual syringe filters.



\*Syringe filters tests done on Merck-Millipore filters  
For more details on Merck Millipore membranes, email techsupport.bioscience@merckgroup.com

### ELECTROLAB Cares

- Verify filter retention volume in addition to chemical compatibility while making a choice of filters
- Reduction in aliquot volume may suggest filter clogging

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice


To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)



## USP Apparatus 1 (Basket)

### Critical Considerations

Recommend your friend or colleague



**Basket Height**  
(USP:  $37.0 \pm 3.0$  mm)



**Basket ID**  
(USP:  $20.2 \pm 1.0$  mm)




**Basket OD**  
(USP:  $22.2 \pm 1.0$  mm)

Cylindrical symmetry



Wobble




Opening




## ELECTROLAB Solutions

### USP 35 Basket Go-NoGo Gauge

Height  
Cylindrical  
symmetry



ID  
OD



Cylindrical asymmetry

USP  
Ph.Eur  
JP

Regulatory Compliant



Reduce Investigation Report



Method Transfer

1min  
Verification

Time Saving



Increases Productivity

Basket Go-NoGo gauge not only controls the parameters set in USP <711> (height, cylinder i.d., and o.d.) but also controls the cylinder symmetry

### ELECTROLAB Cares

- Always hold the basket via the top lip of the basket and avoid putting pressure on the mesh
- Store the baskets in specially designed storage containers



### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Impact of Paddle / Basket Height on Dissolution Testing

Recommend your friend or colleague

"The distance between the inside bottom of the vessel and the basket or paddle blade should be 25 mm  $\pm$  2 mm" Ref : USP <711>

### Effects of inappropriate distance



Shorter distance increases hydrodynamics on dosage



Greater distance reduces hydrodynamics on dosage



Precisely controlled paddle height can reduce variability in results by :

- 3.5% of standard deviation
- 0.5% of mean percent dissolved



As a component of 2-factor analysis when combined with non-manufacturer's vessels, it could contribute to :


- 24.2% variability in mean percent dissolved
- 9.8% variability in standard deviation results

Ref : Eaton J et al. Perturbation Study of Dissolution Apparatus Variables - A Design of Experiment Approach. Dissolution Technologies;14(2):20-26.

USP recommends use of vernier caliper or depth gauge to measure the distance

## ELECTROLAB Solutions

### Calibrated Digital Depth Gauge



Digital value



Easy to use



Accurate value

Serialized laser marked

Digital Depth Gauge for 1 ltr vessel



Available for 2 ltr and 4 ltr vessel



Data transfer



Print



Certified

### ELECTROLAB Cares

- Heights should be periodically checked to ensure reliability of results
- ELECTROLAB recommends to calibrate the depth gauge once a year
- For longer operational life of paddles, baskets and depth gauge, store them in specially designed storage containers



### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice


To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Importance of Rinsing in Dissolution Testing

Recommend your friend or colleague


### What is Contamination?

The sample which is rendered impure by the addition of foreign substance - usually debasing, or devaluing the first



### What is Carryover?

Effect of previous sample on the next sample due to the remainder in the transfer tubings




	Contamination	Carryover
How to reduce?	Cleaning cycle	Rinse cycle
When to run the cycle?	After the test	Before sample withdrawal

### Types of rinse cycle

1) **Simple Rinse** : Purge volume of sampling media to waste


2) **Split rinse** :  
Purge small volume of sampling media by splitting into 2 parts to waste to ensure effective rinsing



Split Rinse

**57%** reduction in rinse volume compared to simple rinse

3) **Loop back** :  
Circulate the initial volume through the transfer tubings back to vessel to eliminate waste volume



Loop Back

**0%** Waste volume

For your complimentary copy of rinse cycle verification SOP, contact ELECTROLAB

## ELECTROLAB Solutions

### Smart Syringe Pump & Sample Collector



Syringe Pump (ESP-64)



Rinsing



Filtration



Maintenance free



Virtually 0 Dead Volume



Sticky Media

All Media Handled



Sample Collector (Dx)

### ELECTROLAB Cares

- Rinsing cycle should be run after every sample to ensure better sampling results
- No. of cleaning cycles should be programmed to run after every test
- The use of dedicated sampling cannulas, for manual sampling minimizes contamination
- Thoroughly clean entire apparatus before start of test

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)


DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Dosage Dispensing in Dissolution Testing


Recommend your friend or colleague

- Timing of drop of one dosage to another
- Timing of test initiation
- Ensuring non rotation of paddles during dosage introduction




- Temperature of media when dosage is introduced
- Position of dosage introduction
- Drop velocity


## ELECTROLAB Solutions




### Working of Auto DropSync™




① Paddle is rotating as set temperature is yet to be attained



② Set temperature is attained and recorded in test report and paddle has stopped rotating



③ Dosage form is introduced in the vessel



④ The test is initiated

## ELECTROLAB Cares

- Regular cleaning ensures smooth operation of the dosage dispenser
- Removal process of the dosage dispenser :




Step 1

Remove gasket from dosage dispenser



Step 2

Press two slider release buttons from below



Step 3

Remove the slider to clean

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : [sales@electrolabgroup.com](mailto:sales@electrolabgroup.com)

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice


To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)

## Friability Testing

Recommend your friend or colleague

Problem	Reason	Solution
<p>The strength of a tablet plays a very important role in its marketing &amp; dissolution. Less strong tablets by losing their shape not only lack elegance &amp; consumer acceptance but also spoil the areas of manufacturing such as coating &amp; packaging.</p>	<p>Measuring the hardness of a tablet is not a reliable indicator for tablet strength as some formulations when compressed into hard tablets tend to 'cap' or lose their crown portion on attrition. Such tablets tend to powder, chip &amp; fragment.</p> <p>Factors affecting friability:-</p> <ol style="list-style-type: none"> <li>1) Poor/worn punches</li> <li>2) Internal factors like moisture</li> </ol>	<p>The Friability test is closely related to tablet hardness and is designed to evaluate the ability of the tablet to withstand abrasion, Friction &amp; shock in packaging, handling &amp; shipping. The loss in weight due to abrasion is a measure of the tablet friability.</p>

### Procedure



### Formula

$$\% \text{ Weight Loss} = \frac{\text{Initial Weight} - \text{Final Weight}}{\text{Initial Weight}} \times 100$$

## ELECTROLAB Solutions

 <p><b>Friability Tester (EF – 2W)</b> with weighing scale interface</p>	 <p><b>Abrasion Drum</b> Abrasion drum for more rigorous informational testing of tablet and coating integrity.</p>	 <p><b>10° Tilt</b> If tablet size or shape causes irregular tumbling, adjust the drum base so that the base forms an angle of about 10° with the bench top and the tablets no longer bind together.</p>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------	------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------	---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

## ELECTROLAB Cares

### Validation parameters\*

Parameters	A	B	C	D	E	F	G
Description	Center hole dia.	Center disc dia.	Inside dia.	Outer dia.	Depth	Tablet drop	Inside curve
USP Specification (mm)	10 ± 0.1	25 ± 0.5	287 ± 4.0	302.5 ± 4.0	38 ± 2.0	156.0 ± 2.0	80.5 ± 5.0

\* Contact ELECTROLAB for detailed procedure for verification of above parameters.

### Maintenance

- Routinely check drum(s) for signs of cracks or chips
- Clean drum thoroughly using de-ionized water and dry with a soft cloth

### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our **Specialist** on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)




## Is your Tap Density Tester compliant to USP 39-NF 34?

### USP Chapter <616> Bulk Density and Tapped Density of Powders

Features	METHOD 1		METHOD 2	
Taps/min and Fixed drop height	Taps/min	Height	Taps/min	Height
	300 ± 15	14 mm ± 2 mm	250 ± 15	3 mm ± 0.2 mm
	250 ± 15	3 mm ± 0.2 mm		
Number of tap counts	10, 500, 1250			
Specifications for graduated cylinder	250 ml (readable 2 ml)		100 ml (readable 1 ml)	
Weight of graduated cylinder	220 g ± 44 g		130 g ± 16 g	
Weight of support for the graduated cylinder and holder	450 g ± 10 g		240 g ± 12 g	
Difference between succeeding measurements	≤ 2 ml		≤ 1 ml	

#### Methodology



#### Note:

- Weight of sample ≤ 100 g.
- Tapped density, Compressibility index and Hausner ratio can be calculated using the tap density tester.

### ELECTROLAB Solution



(ETD – 1020x)



USP Method 1 and USP Method 2



Simultaneous Rotation and Tapping Motion



Snap-Lock™ Mechanism



Date and Time Lock



SS Construction (optional)



Alphanumeric Display



Results



Print



LAN (optional)



User Friendly Operation

- Test can be performed at two different modes :
  - **USP mode** : Test will run at set number of taps as per USP requirements
  - **User mode** : Test will run at set number of taps as per user requirements
- Facilitates calculations of Tapped density, Compressibility index and Hausner ratio.

For firmware upgradation, contact ELECTROLAB

#### Contact Us

401 Tirupati Udyog, I. B. Patel Road,  
Off Western Express Highway,  
Goregaon (E) Mumbai - 400 063, India.  
☎ : 91 22 4041 3131  
✉ : sales@electrolabgroup.com

To know more, please contact our *Specialist* on [support@electrolabgroup.com](mailto:support@electrolabgroup.com)

DISCLAIMER : The information contained in this document is believed to be correct, But ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice

To unsubscribe click [unsubscribe@electrolabgroup.com](mailto:unsubscribe@electrolabgroup.com)