

CLEARVIEW



ClearView™

Features:

The ClearView system uses UV technology as a water purification system serving as a germicide against algae, bacteria and other microorganisms helping to maintain a consistently clear water bath.

ClearView

- ❖ Increases instrument uptime and throughput,
- ❖ Reduces labour costs due to bath cleaning and maintenance.
- ❖ Lowers microorganism count
- ❖ Enhances visual observation of the dissolution process

Life of the bulb : 8000 hours

ClearView™ can be retrofitted to almost any ELECTROLAB dissolution tester.

Ordering Information:

Products		
Product code	Description	Quantity
0221300	ClearView™(Body made of M.S.) -230V	1
0211300	ClearView™ (Body made of M.S.) -110V	1
0104A00432	EDT Clearview (230V) Version-2	1
0104A00491	EDT Clearview (110V) Version-2.1	1
0104A00477	EDT Clearview (230V) Version-2.1	1

Standard accessories		
Part number	Description	Quantity
1536A00002	ECV-16 Filter Cartridge (Set of 6)	1

Optional accessories		
Product code	Description	Quantity
1536A00004	Clearview Filter Assembly Nozzle	1
1536A00012	ECV Flow Pure Spun 10" long 100 micron (Set of 3)	1
1536A00011	ECV Super Pure Spun 10" long 50 micron Filter (Set of 3)	1
1536A00010	ECV coupling (Uv to flow sensor)	1
1536A00009	ECV filter Mounting Bracket assembly	1
1536A00008	ECV-16 Ultra Filter Cartridge (set of 3)	1
1536A00007	ECV-16 50micron Filter Cartridge (Set of 6)	1
1536A00006	ECV-16 20micron Filter Cartridge - 10" (Set of 6)	1
1536A00005	ECV Y connector	1
1536A00003	ECV-16 Filter	1
1536A00001	ECV-16 UV Lamp	1

Breakdown accessories		
Product code	Description	Quantity
1536B00014	ECV UV Lamp Adapter	1
1536B00013	ECV Front Panel	1
1536B00012	ECV Wiring set with water flow sensor	1
1536B00011	ECV Flow Indicator PCB Ver.1.0 (1803001)	1
1536B00010	ECV-16 Hours Monitor	1
1536B00009	ECV-16 Fuse Slow Blow (250mA)	1
1536B00008	ECV-16 Mains Cord – 5Amp	1
1536B00007	ECV-16 Wago Connector (2 Way)	1
1536B00006	ECV-16 Fuse Holder	1
1536B00005	ECV-16 2 Pin Screw Terminal	1
1536B00004	ECV-16 BLDC Pum	1
1536B00003	ECV-16 SMPS	1
1536B00002	ECV-16 BLDC ECV-16 LED Lamp (Panel Mount)	1
1536B00001	ECV-16 Ballast	1

Assurance Plan		
Product code	Description	Quantity
01A491S01	Platinum Assurance Plan for EDT Clearview	1
01A491S02	Gold Assurance Plan for EDT Clearview	1
01A491S03	Silver Assurance Plan for EDT Clearview	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.