


sentina⁰

Packaging Data Digitization



ANY Brand
Devices



MetaData
acquisition

Can Be Connected !

FEATURES

- Sentinai provides up to 4 port connectivity between legacy instruments from any brand and the server For seamless data transfer.
- Sentinai having inbuilt storage of 2,00,000 reports for each port.
- It is designed to work with lims and mes for Transferring data over the network.
- Generate results reports through a web Browser and print them on any local or Network printer.
- Enhanced data security, precision, and Traceability.
- Visual & led indication for data capture and Transfer for each port.
- User Managment with three level configurable roles.



ADVANTAGES

- Inbuilt storage eliminates chances of data loss during transfer.
- Automatic PDF/CSV transfer to FTP/SMB server.
- There is no requirement to install any software.

BENEFITS

- Inbuilt storage enhances data accuracy to minimize out-of-specification instances.
- Replaces up to 4 physical printers to save maintenance & paper cost & space.
- Data accessible in digital format facilitating Readability and searchability.
- IQOQ + Validation.
- Personalized Company Logo can be added.

Sentinai Licenses	Basic	Machine Learning	Machine Learning Advance
PDF Transfer	✓	✓	✓
Report Capture	Batch	Batch/Online	Batch + Online
CSV & other file format			✓

SPECIFICATIONS		
Model	Sentinai	
No. of Instrument Connection	2	4
Display	OLED Screen Display	
Network connectivity	LAN + Wifi	
WIFI connectivity	Available	
IP	DHCP/Static	
Utility required for configuration	No web access	
Instrument Wise Storage	Yes (Auto)	
Report Storage	Auto on device + Network	
Report Format	PDF , CSV	
Report Print	Network Printer	
Internal Memory	8 lakh+ Test report Storage	
Connectivity	FTP/SMB	
LIMS,MES Connectivity (OPC Module)	Yes (Optional)	
Data Out	LAN / Auto Print / PDF Over FTP/ SMB	
Dimension (LxWxH)	150 mm X 105 mm X 53 mm	
Weight	1 kg	