

# ELECTROLAB

Your Quality, Our Assurance

## XTRACTR

LCD  
40X4

Display

D'light™  
● ● ● ●

Technology

02


Station



Xtainer



Report Printing



### Applications


- ❖ Fast and efficient sample preparator for content uniformity, assay and impurity testing samples
- ❖ According to USP <905> chapter, Content Uniformity is required for dosage forms containing 25 mg or less of a drug substance

### Robust Design

- ❖ 02 station individually powered closed system
- ❖ Autodetection of Xtainer type
- ❖ Variable Xtainer capacities to accommodate all dosage forms and label claims
- ❖ Amber Xtainer to protect light sensitive actives and inactives
- ❖ Threaded lids secure Xtainer contents
- ❖ Motor capable of reversing direction in case of stalling
- ❖ Audio visual alerts to indicate instrument status


### User friendly Interface

- ❖ Simple and easy operable LCD display
- ❖ 10 method storage capacity
- ❖ Report printing via serial printer



## Accessories

Part no	Description (Quantity - set of 10)
0104A00536	Xtractr 50 mL clear Xtainer with stainless steel blade
0104A00537	Xtractr 50 mL amber Xtainer with stainless steel blade
0104A00540	Xtractr 250 mL clear Xtainer with stainless steel blade
0104A00541	Xtractr 250 mL amber Xtainer with stainless steel blade
0104A00544	Xtractr 500 mL clear Xtainer with stainless steel blade
0104A00545	Xtractr 500 mL amber Xtainer with stainless steel blade



## Specifications

Model	Xtractr- 2T
Number of stations	02
Stirring range	1500 – 6000 RPM
Jar capacities	50 mL, 250 mL and 500 mL
Jar type	Clear, Amber
Blade type	Stainless steel
Test duration	10 sec – 99 min 59 sec
Display	Alphanumeric, 40 x 4 LCD
Methods	10
Storage	Last report
Printer Interface	RS232

ELECTROLAB

## Our Products

- Complete range of Dissolution Testers • Dissolution Media Preparator • Friability Testers • Tablet Hardness Testers
- Electromagnetic Sieve Shakers • Tap Density Tester • Bulk Density Tester • Powder Flow Tester • Leak Testers • Peristaltic Pumps

The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice  
December, 2020

**ELECTROLAB (India) PVT. LTD.** (CIN NO : U33110MH2002PTC136735)

**Head Office:** 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

**Factory:** EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

