

Dissolve Oxygen Meter EzDO with Ezintegrity



Features:

- 21 CFR part 11 compliant with audit trail
- Corrosion resistant HIPS plastic body
- Easy to use, built-in, menu- driven instructions on 7" touch screen TFT display
- Mouse and keyboard connectivity
- Direct connectivity (without PC) for wi-fi, network or USB printer (printer having Linux support)
- Pen drive connectivity
- GLP Compliance reports with Time, Date and sample identification number facility
- 1000 reports can be stored
- Automatic Temperature Compensation
- Data Logging and Reporting :
The authenticity of the reports generated can be ensured with the Date / Time indication in the report.
The reports can be seen on the display
- Automatic Data Logging :
To store the DO value at specified time intervals. (Selectable from 1 to 999 minutes)
Upto 50 readings can be taken. This facility is particularly useful for stability study of samples.
- Manual Data Logging :
Large memory to store 1000 reports

Specifications:

Model	EzDO with Ezintegrity
Parameter	DO
DO Range	0 ppm to 20 ppm
DO Resolution	0.001 ppm
DO Accuracy	± 2 % FS
DO Units	ppm
Dimensions (L x W x H)	230 mm x 190 mm x 80 mm
Weight	1 Kg (approx)
Power Supply	24V DC Mains Adaptor

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com