

Compact Peristaltic Pump

3 years warranty

EL- Handy Pump




Product Features

- 304 stainless steel shell, anti-corrosion, resistant to rust, conform to GMP request.
- OLED display motor speed and working status, digital knob control speed.
- EL- HandyPump external control adopts standard wiring terminal, which is more suitable for supporting industrial use.
- Support RS232, RS485 communication, Modbus protocol, more convenient to connect with PLC.
- With power-down memory function, cost-effective.
- The external control function interface is rich, and the external control isolation signal is 5-24V.
- Compact structure, can be stacked.

Model number - EL-HandyPump

Dimension Drawing (Unit:mm)



Pump Head



EL-HandyPump01



EL-HandyPump02

Technical Specifications

Flow rate range	0.0033~365.69mL/min	Power supply	DC 24V/1A with power adapter
Speed range	0.1~300rpm	External control	Start/stop control(switch signal)
Speed resolution	0.1rpm		Speed: 0-5V, 4-20mA(standard), 0-10V (optional)
Control method	Digital knob control/Mechanical keypad	Communication interface	RS485, RS232
Motor type	Stepper motor	Condition temperature	0-40℃
Display	OLED display(0.96")	Relative humidity	<80%
Power consumption	15W	Output interface	Output motor working status
Weight	1.75Kg	IP rate	IP31
Keypad lifetime	300,000 times	Drive dimension	198*113*118mm(L*W*H)

Product Composition and Flow Rate Range

Model	Pump Head	Speed (rpm)	Tubing Size	Flow Rate Range (mL/min)
EL- Handy Pump	EL- HandyPump01	0.1-300	13",14", 19", 16", 25"	0.0033-365.69
	EL- HandyPump02		13",14", 19", 16"	0.0033-190.76

Our Products

- Complete range of Dissolution Testers • Dissolution Media Preparator • Friability Testers • Tablet Hardness Testers
- Electromagnetic Sieve Shakers • Tap Density Tester • Bulk Density Tester • Powder Flow Tester • Leak Testers • Peristaltic Pumps

The information contained in this document is believed to be correct but **ELECTROLAB** accepts no liability for any errors and reserves the right to alert specifications without notice.

JUNE, 2021

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory : EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 71

+91 - 22 - 4161 3122 www.electrolabgroup.com

