

Content Uniformity Sample Preparator Xtractr



Xtractr-10

Features:

- Fast and efficient sample preparator for content uniformity, assay and impurity testing samples
- According to USP <905> chapter, Content Uniformity is required for dosage forms containing 25 mg or less of a drug substance
- 10 station individually powered closed system
- Autodetection of Xtainer presence
- Variable Xtainer capacities to accommodate all dosage forms and label claims
- Amber Xtainers to protect light sensitive actives and inactives
- Threaded lids secure Xtainer contents
- Motor capable of reversing direction in case of stalling
- Audio visual alerts to indicate instrument status
- Non-contact temperature detection of individual Xtainers
- Motor auto shutdowns on reaching 'off' temperature and auto restarts on reaching 'on' temperature
- Simple and easy operable LCD display
- Sufficient method storage capacity
- Report printing via serial printer

Specifications:

Model	Xtractr-10
No. of stations	10
Stirring range	1500 – 6000 rpm
Temp. range	20 - 80°C
Jar capacities	50 mL, 250 mL and 500 mL
Jar type	Clear, amber
Blade type	Stainless steel
Test duration	10 sec – 99 min 59 sec
Display	Alphanumeric, 20 x 4 LCD
Methods	10
Storage	Last report
Printer Interface	RS232
Dimensions (L X W X H)	698 X 396 X 401 mm
Weight	15 kg

Ordering information:

Products		
Product code	Description	Quantity
534020	Xtractr-10 (Body made of GMP Compliance) -230V	1
533020	Xtractr-10 (Body made of GMP Compliance) -110V	1

Optional accessories		
Part number	Description	Quantity
0104A00537	Xtractr 50 mL amber Xtainer with stainless steel blade	1
0104A00541	Xtractr 250 mL amber Xtainer with stainless steel blade	1
0104A00545	Xtractr 500 mL amber Xtainer with stainless steel blade	1
0104A00544	Xtractr 500 mL clear Xtainer with stainless steel blade (Set of 10)	1
0104A00540	Xtractr 250 mL clear Xtainer with stainless steel blade (Set of 10)	1
0104A00536	Xtractr 50 mL clear Xtainer with stainless steel blade (Set of 10)	1
0104A00666	Rpm Validation Tool For Xtractr	1

Breakdown Parts		
Part number	Description	Quantity
0104B00102	Xtractr-10 Top Cover	1

Assurance Plan		
Part number	Description	Quantity
5501S01	Platinum Assurance Plan for Content Uniformity Sample . Preparator Model Xtractr-10	1
5501S02	Gold Assurance Plan for Content Uniformity Sample Preparator Model Xtractr-10	1
5501S03	Silver Assurance Plan for Content Uniformity Sample Preparator Model Xtractr-10	1
5501S20	Bronze Assurance Plan for Content Uniformity Sample Preparator Model Xtractr-10	1

IQ/OQ Document		
Part number	Description	Quantity
5501S15	IQ/OQ Document for Content Uniformity Sample Preparator Model Xtractr-10	1

Inspection & Service Charges		
Part number	Description	Quantity
5501S15	Inspection & Service Charges for Content Uniformity Sample Preparator Model Xtractr-10	1

Service Charges		
Part number	Description	Quantity
5501S05	Service Charges for Content Uniformity Sample Preparator Model Xtractr-10	1

Upgradation		
Part number	Description	Quantity
5501S13	Upgradation Charges For Content Uniformity Sample Preparator Model Xtractr-10	1

Reinstallation Charges		
Part number	Description	Quantity
5501S11	Reinstallation Charges with IQ/OQ Document for Content Uniformity Sample Preparator Model Xtractr-10	1
5501S10	Reinstallation Charges without IQ/OQ Document for Content Uniformity Sample Preparator Model Xtractr-10	1

Re-calibration Charges		
Part number	Description	Quantity
5501S09	Re-calibration Charges for Content Uniformity Sample Preparator Model Xtractr-10	1

Inspection Charges		
Part number	Description	Quantity
5501S04	Inspection Charges for Content Uniformity Sample Preparator Model Xtractr-10	1

In-House Repairing		
Part number	Description	Quantity
5501S07	In-House Repairing for Content Uniformity Sample Preparator Model Xtractr-10	1

Dismantling Charges		
Part number	Description	Quantity
5501S12	Dismantling Charges for Content Uniformity Sample Preparator Model Xtractr-10	1