


DermaX - 07

Manual Dry Heat Diffusion Cell



DermaX - 07

Features:

- Biometric login
- Interactive touch screen
- Cells available in variable size
- Electronic Logbook and Audit Trail
- Power saving heaters for each cell
- Infrared based temperature monitoring
- Variable volumes between each cell and the heating block
- >60% contact between each cell and the heating block
- Provided 21 CFR Part 11 Compliance
- Increased protocol storage space, number of users and access levels

Specifications:

Model	DermaX - 07
Number of stations	7
Speed range	200 RPM to 1000 RPM
Speed accuracy	± 2 RPM
Temperature range	20°C to 45°C
Temperature accuracy	± 1.0°C
Display	Multi Touch Capacitive 7" TFT Touch Screen
Test Report Storage	10,000
Heater	200 W
Report Printout	Ethernet, RS-232
No. of Sampling Interval	18
Dosage Volume	0.25 - 6.2 ml
Cell Type	Small, Medium and Large (Clear and Amber)
Water bath	One piece moulded clear
Cell Volume	8 - 30 ml
Protocol Storage	1000
Dimensions (L x W x H)	640 mm x 540 mm x 660 mm
Electrical power	220/230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/10 Amp
Illumination	D'Light
Security	Multilevel with password for 250 users with 10 user programmable role
Weight	11 Kg

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

DermaX - 07

Manual Dry Heat Diffusion Cell

Ordering information:

Products		
Product code	Description	Quantity
1020300	Manual Dry Heat Diffusion Cell Model DermaX-07 (MS)	1

Standard Accessories		
Product code	Description	Quantity
1001A00094	Dermax Bubble Detection Mirror	1

Optional Accessories		
Product code	Description	Quantity
1001A00083	Dermax Small Cell, Clear	1
1001A00105	Small Cell Sleeve	1
1001A00092	Cell Cap Kit, Small, 9mm X 4mm	1
1001A00093	Cell Cap Kit, Small, 9mm X 20mm	1
1001A00106	Cell Cap Kit, Small, 11.3mm X 4mm	1
1001A00107	Cell Cap Kit, Small, 11.3mm X 20mm	1
1001A00087	Small Cell Stirrer, Small	1
1001A00086	Small Cell Stirrer, Big	1
1001A00084	Dermax Medium Cell, Clear	1
1001A00104	Medium Cell Sleeve	1
1001A00096	Cell Cap Kit, Large, 15mm X 4mm	1
1001A00097	Cell Cap Kit, Large, 15mm X 20mm	1
1001A00106	Cell Cap Kit, Small, 11.3mm X 4mm	1
1001A00107	Cell Cap Kit, Small, 11.3mm X 20mm	1
1001A00089	Medium Cell Stirrer, Small	1
1001A00088	Medium Cell Stirrer, Big	1
1001A00085	Dermax Large Cell, Clear	1
1001A000102	Cell Cap Kit, Large, 20mm X 4mm	1
1001A000103	Cell Cap Kit, Large, 20mm X 20mm	1
1001A00090	Large Cell Stirrer, Big	1
1001A00091	Large Cell Stirrer, Small	1
1001A00102	Dermax Replenishment Bottle (Set of 7)	1
1001A00100	Dermax Cell Evaporation Cap (Set of 15)	1
0104A00439	IntegrOS lite version	1

IQ/OQ Document		
Product code	Description	Quantity
1003S15	IQ/OQ Document for Manual Dry Heat Diffusion Cell Model DermaX-07	1

Extended IQ/OQ Document		
Product code	Description	Quantity
E1003S15	Extended IQ/OQ Document for Manual Dry Heat Diffusion Cell Model DermaX-07	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)
Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.
 ☎ +91 - 22 - 4041 3131 ✉ sales@electrolabgroup.com
Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.
 ☎ +91 - 22 - 4161 3122 🌐 www.electrolabgroup.com