

Tablet Dissolution Tester Inspire-08 (Basic)



Inspire-08 (Basic)

Features:

- LCD display of 40 characters and 4 lines for user friendly operation
- Multi level security with password
- Alphanumeric data entry of product and user details
- Unique TimeAction™ function which allows media changeover and infinity test run, ideal for sustained and controlled release products
- Complies with USP, IP, Ph. Eur. and JP specifications
- Automatic stirrer height positioning for USP 1, 2, 5, 6, intrinsic test methods
- Easy Snap-Fit™ shaft locking mechanism for positive engagement and wobble free operation
- Clear water bath with metal top-plate
- Precise individual vessel centering system
- Adjustable leveling of instrument
- Vibration free smooth motorized lift movement
- Isolated water circulating pump for precise temperature control of the water bath
- 20 programmable protocols
- 24 programmable sampling intervals
- Easy to empty and clean water bath with bottom drainage
- Designed to minimize routine validation
- Report for online validation of test, protocol, error log during the test
- Programmable wake-up for temperature control
- Programmable sample volume from 1 mL to 20 mL
- Audio-visual indicators for system status
- Power failure detection facility

Specifications:

| Model | Inspire-08 (Basic) |
|-------------------------|--|
| Speed range | Variable speed from 20 to 300 RPM |
| Speed accuracy | +1 RPM |
| Temperature range | 20° C to 40° C (Extended range 10°C to 40°C & 20°C to 55°C optional) |
| Temperature accuracy | (±0.1° C) |
| Display | 40 x 4 character backlighted LCD |
| No. of protocols | 20 maximum |
| Sample interval | Maximum 24 of 1 min to 23 hrs, 59 min |
| Printer interface | Serial port |
| Printer supported | 80 column Dot Matrix printer |
| Bath | Clear molded water bath |
| Heater | 1 kW/230 VAC, material S.S. 304 with sensor |
| Vessel capacity | 150 ml, 250 ml, 1 L |
| Electrical power | 220/ 230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/10 Amp |
| Dimensions (W x L x H) | 517 x 712 x 665 mm |
| Weight | 68 KG |

Ordering Information:

| Products | | |
|--------------|--|----------|
| Product code | Description | Quantity |
| 0128204 | Tablet Dissolution Tester Model Inspire-08 -230V (Basic) | 1 |
| 0118204 | Tablet Dissolution Tester Model Inspire-08 -110V (Basic) | 1 |

| Standard accessories | | |
|----------------------|--|----------|
| Part number | Description | Quantity |
| 0104BV4155 | EDT V-5 External Probe Without Trim Pot | 1 |
| 0104A00075 | EDT laser marked basket - 40 mesh (Set of 6) | 1 |
| 0104A00066 | EDT basket rod | 6 |
| 0104A00068 | EDT Basket Holder with Clip (Set of 6) | 1 |
| 0104A00592 | Dissolution Inspire Detachable paddle blade (Set of 6) | 1 |
| 0104A00610 | Stirrer for replenishment and blank vessels | 2 |
| 0104A00041 | EDT moulded jar lid - plastic (Set of 6) | 1 |
| 0104A00242 | EDT merlon vessel -1000ml clear (N) | 8 |
| 0104BV4114 | EDT Temperature Controller (ETC-15) | 1 |

| Assurance Plan | | |
|----------------|---|----------|
| Part number | Description | Quantity |
| 0182S01 | Platinum Assurance Plan for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| 0182S02 | Gold Assurance Plan for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| 0182S03 | Silver Assurance Plan for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| 0182S20 | Bronze Assurance Plan for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |

| Reinstallation Charges | | |
|------------------------|--|----------|
| Part number | Description | Quantity |
| 0182S11 | Reinstallation Charges with IQ/OQ Document for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |

| IQ/OQ Document | | |
|--|--|----------|
| Part number | Description | Quantity |
| 0104BV4225 | IQ/OQ Document for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Inspection & Service Charges | | |
| Part number | Description | Quantity |
| 0182S06 | Inspection & Service Charges for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Upgradation | | |
| Part number | Description | Quantity |
| 0182S13 | Upgradation for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| ASTM Validation Charges | | |
| Part number | Description | Quantity |
| 0182S19 | ASTM validation charges for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Performance Qualification (PQ Documents) | | |
| Part number | Description | Quantity |
| 0182S17 | Performance Qualification (PQ Documents) for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Performance Verification Test | | |
| Part number | Description | Quantity |
| 0182S14 | Performance Verification Test (PVT Assist) for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Single Validation Charges | | |
| Part number | Description | Quantity |
| 0182S08 | Single Validation Charges for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |
| Dismantling Charges | | |
| Part number | Description | Quantity |
| 0182S12 | Dismantling Charges for Tablet Dissolution Tester Model Inspire-8 Basic | 1 |

| Sr.No. | Part code | Breakdown accessories |
|---------------|------------------|--|
| 1 | 0104BV4154 | EDT-Lx V5 Heater Card for ETC-15 Without Trim Pot |
| 2 | 0104BV4153 | EDT-Lx V5 Compatch Card Without Trim Pot |
| 3 | 0104BV4151 | EDT-08Lx V-5 Individual Vessel Temp. Sensor PCB with Cable & Connector |
| 4 | 0104BV4147 | EDT 08L V-1234 Diving Manifold Bearing's Housing (Set of 2) |
| 5 | 0104BV4145 | EDT-Lx V-4 Driver Card (110V) |
| 6 | 0104BV4141 | EDT-Lx V-45 Mains Supply Filter for ETC-15 (230V) |
| 7 | 0104BV4140 | EDT-Lx V-4 SMPS - 12V |
| 8 | 0104BV4137 | EDT-LX V-4 Mains Cord |
| 9 | 0104BV4135 | EDT-Lx V-4 Inlet Connector for Side Circulation Waterbath |
| 10 | 0104BV4115 | EDT-08Lx V-4 - 10 Pin FRC Cable (Processor Card to Stirrer Opto Card) |
| 11 | 0104BV4114 | EDT-Lx V-4 Temperature Controller ETC-15 |
| 12 | 0104BV4091 | EDT Lx V-4 Communication Cable (Syringe Pump to Sample Collector) |
| 13 | 0104BV4090 | EDT Lx V4 Communication Cable (Disso to Syringe Pump) |
| 14 | 0104BV4087 | EDT-Lx V-4 Serial Printer Cable |
| 15 | 0104BV4086 | EDT-Lx Opto PCB For Encoder |
| 16 | 0104BV4071 | EDT-Lx Encoder With Flywheel |
| 17 | 0104BV4054 | EDT-08Lx V-4 – 20 Pin FRC Cable (Processor Card to Keypad) |
| 18 | 0104BV4048 | EDT-Lx V-4 Lift Opto Bottom |
| 19 | 0104BV4030 | EDT-Lx V4 Main Cable 15 AMP |
| 20 | 0104BV4016 | EDT-Lx V-4 LCD Display |
| 21 | 0104BV4005 | EDT-Lx V-4 RPM Stirrer Motor |
| 22 | 0104B04138 | EDT V-1234 Opto Disc for Lift Motor |
| 23 | 0104B00168 | EDT-08L/14L V-234- 3V DC Battery for Processor Card |
| 24 | 0104B00099 | EDT V-1234 Link Cord |
| 25 | 0104B00096 | EDT V-1234 Bath Coupling Mounting Plate (Set of 2) |
| 26 | 0104B00087 | EDT V-1234 Inlet & Outlet Connector for |
| 27 | 0104B00065 | EDT-08L V-1234 Bath Top Plate Knob with Washer (Set of 4) |
| 28 | 0104B00052 | EDT-08L V-1234 Small Rubber Leg for Water Bath 22mm (Set of 4) |
| 29 | 0104B00050 | EDT V-1234 Drain Rubber for Water Bath |
| 30 | 0104B00048 | EDT Dispensor Plate Holder (Set of 2) |

| | | |
|----|------------|--|
| 31 | 0104B00047 | EDT V-1234 Inline Fuse Holders |
| 32 | 0104B00046 | EDT V1234 Mains Switch (Yellow) |
| 33 | 0104B00045 | EDT V-1234 Rocker Switch Mains Red |
| 34 | 0104B00041 | EDT V-1234 Lift Motor 230V-2A |
| 35 | 0104B00037 | EDT V-1234 Lift Motor 110V-3A |
| 36 | 0104B00032 | EDT V-1234 Links Cable 5 Amp |
| 37 | 0104B00029 | EDT V-1234 Mains Fuse Holder |
| 38 | 0104B00028 | EDT V-1234 Thermal Cut-Off |
| 39 | 0104B00026 | EDT V-1234 Timing Pulley Set |
| 40 | 0104B00025 | EDT V-1234 Drainage Valve |
| 41 | 0104B00009 | EDT V-1234 Manifold Opto |
| 42 | 0104B00008 | EDT-08L V-1234 Manifold Bearing (Set of 2) |
| 43 | 0104B00003 | EDT-08L/14L V-34 Manifold Bush Highspeed |
| 44 | 0104BV4173 | EDT V -1234 LED with Cable (Set of 2) |
| 45 | 0104BV4188 | EDT-Inspire BLDC Drive 8 Pole |
| 46 | 0104BV4177 | EDT V-1234 Cut spacer 0.2mm |
| 47 | 0104BV4176 | EDT V-1234 Cut Spacer 0.5mm |
| 48 | 0104B00012 | EDT V-1234 Individual Vessel Temperature Sensor (Set of 8) |
| 49 | 0104B00014 | EDT V-1234 Single Pulley Spindle |
| 50 | 0104B00015 | EDT V-1234 Double Pulley Spindle |
| 51 | 0104BV4027 | EDT-Lx V-4 Fuse Set |
| 52 | 0104BV4163 | EDT-8/14Lx V-4-Inspire Driver Card |
| 53 | 0104BV4098 | EDT-Lx V-4 Heater- 110V (For ETC-15) |
| 54 | 0104B00101 | EDT-Lx V-4 Heater- 110V (For ETC-15) Old Straight Nozzle |
| 55 | 0104BV4128 | EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 230V |
| 56 | 0104BV4155 | EDT V-5 External Probe Without Trim Pot |
| 57 | 0104BV4183 | EDT-LX V4 Inspire Processor Card |
| 58 | 0104BV4099 | EDT-08L Bath Support Frame |
| 59 | 0104BV4080 | EDT-Lx V4 Compatch Card |
| 60 | 0104BV4085 | EDT-Lx V4 Gear for Lift Motor |
| 61 | 0104BV4158 | EDT-Lx V-45 Advance Filter for ETC-15 (230V) |
| 62 | 0171B00001 | EDT-08/14Tx LCD Module (Touch Screen) |
| 63 | 0104B00039 | EDTV-1234 Heater 230V-2A |
| 64 | 0104BV4185 | EDT- V-4-Inspire Stirrer opto card |
| 65 | 0104B00022 | EDT-08L V-1234 Timing Belt Set |
| 66 | 0104BV4024 | EDT-Lx V-4 Main Socket with Power Supply Switch & Filter |
| 67 | 0104B00134 | EDT-08L V-1234- 16 Pin FRC Cable (Processor Card to Driver Card) |

| | | |
|-----|------------|--|
| 68 | 0104BV4156 | EDT-Lx V-4 Mains Cord For ETC-11Lx |
| 69 | 0104BV4059 | EDT-Lx V-4 Lift Disc Opto Card |
| 70 | 0104BV4223 | EDT-08LX Offline With Fraction Collector Platinum Kit Electronic Spares |
| 71 | 0104B00083 | EDT V-1234 MF Connector for Manifold |
| 72 | 0104B00018 | EDT V-1234 External Probe |
| 73 | 0104B00013 | EDTV-1234 Bath Temperature Sensor |
| 74 | 0171B00009 | EDT-08/14Tx Lift motor drivr card |
| 75 | 0104BV5159 | EDT-08Lx V-5 Individual Vessel Temperature Sensor PCB with Cable & Connector with out trim pot |
| 76 | 0104BV4001 | EDT-Lx V4 Stirrer Opto Card |
| 77 | 0104B00049 | EDT V-1234 Individual Vessel Temperature Sensor |
| 78 | 0104BV4081 | EDT-Lx V-4 Temperature Controller ETC-11 |
| 79 | 0104BV4084 | EDT-Lx V-4 Lead Screw for Lift |
| 80 | 0105BV4164 | EDT- Wear & Tear Kit (For Dx) |
| 81 | 0104BV4178 | EDT V-1234 Height adjustment spacer 0.5mm |
| 82 | 0104BV4074 | EDT-Lx V-4 Processor Card |
| 83 | 0168B00002 | EDT-08Tx Processor Card |
| 84 | 0104BV4159 | EDT-08Lx V-45 Wear & Tear Kit (For Dx) |
| 85 | 0104BV4149 | EDT-08Lx V-4 Individual Vessel Temp. Sensor PCB (Black) with Cable & Connector |
| 86 | 0104BV4221 | EDT-08LX Offline With D Sample Collector Platinum Kit Electronic Spares |
| 87 | 0104B00058 | EDT-08L V-1234 Water Circulation Tygon Tube Set |
| 88 | 0104BV4127 | EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 110V |
| 89 | 0104BV4092 | EDT-Lx V-4 Circulation Pump – 24V (For ETC-15) |
| 90 | 0104BV4095 | EDT-Lx V-4 Heater Card for ETC-15 |
| 91 | 0104BV4175 | EDT V-1234 Height adjustment spacer 0.2mm |
| 92 | 0104B00143 | EDT V-1234 LED For ETC-15 |
| 93 | 0104BV4095 | EDT-Lx V-4 Heater Card for ETC-15 |
| 94 | 0104BV4239 | Inspire 8 Base Frame Assembly |
| 95 | 0104BV4237 | Inspire Processor card for 10 to 55 degC |
| 96 | 0104BV4201 | Inspire Lift Park Position Opto Sensor |
| 97 | 0104BV4238 | Inspire Program Patch |
| 98 | 0104BV4225 | Inspire-08 Platinum Kit Mechanical Spares. |
| 99 | 0104BV4226 | Inspire-08 Platinum Kit Electronic Spares. |
| 100 | 0104BV4225 | Inspire-08 Platinum Kit Mechanical Spares. |
| 101 | 0104BV4226 | Inspire-08 Platinum Kit Electronic Spares. |
| 102 | 0104BV4239 | Inspire 8 Base Frame Assembly |

| sr.no. | Part code | Optional accessories |
|--------|-------------|---|
| 1 | 0104A00583 | EDT 20 Mesh Stationary Basket (Set of 6) |
| 2 | 0104A00706 | EDT-Inspire Syringe manifold Tygon Tube with SS Sleeve (Set Of 6) |
| 3 | 0104A00693 | EDT Rubber Holder For External Probe.(Set Of 6) |
| 4 | 0104A00609 | EDT 62mm Paddle Over Disc Assembly With 16 Mesh(set of 6) |
| 5 | 0104A00625 | EDT Moulded Jar LID - Plastic |
| 6 | 0104A00631 | EDT – Adjustable Legs Long (Set of 4) |
| 7 | 0104A005893 | EDT Teflon Paddles (set of 8) |
| 8 | 0104A00596 | EDT Paddle/Basket Rod |
| 9 | 0104A00362 | EDT 85mm Paddle Over Disc Assembly |
| 10 | 0104A00142 | EDT-08L/08Lx Water Bath with Plastic QRC - Clear |
| 11 | 0104A00780 | EDT-08 Bath Top Plate For Hexa Vessel |
| 12 | 0104A00048 | EDT Adjustable Sampling Cannula for low evaporation jar lid (Set of 8) |
| 13 | 0104A00272 | EDT Gasket For Tablet Dispensor Plate |
| 14 | 0104A00112 | EDT Helical Coil Sinkers |
| 15 | 0104A00046 | EDT Adjustable Sampling Cannulas For Tablet Input Device (Set of 6) Length 144 mm |
| 16 | 0104A00205 | EDT Lid for Basket 40 Mesh 24.5 mm (Set of 6) |
| 17 | 0104A00182 | EDT Inside Caliper 6" |
| 18 | 0104A00782 | EDT 90mm Watch Glass Disc (Set of 6) |
| 19 | 0104A00400 | EDT Laser Marked Glass Vessel 1000 ML - Clear With Sigma Delta Certificate |
| 20 | 0122A00010 | EDT Depth Gauge for 2 Ltr Vessels |
| 21 | 0104A00204 | EDT Basket 40 Mesh 24.5 mm (Set of 6) |
| 22 | 0104A00128 | EDT Intrinsic Attachment |
| 23 | 0104A00146 | EDT Plastic Quick Release Coupling with Shut Off Valve (Set of 2) |
| 24 | 0104A00399 | EDT Glass Vessel 1000 ML - Amber With Sigma Delta Certificate |
| 25 | 0104A00398 | EDT Glass Vessel 1000 ML – Clear With Sigma Delta Certificate |
| 26 | 0104A00138 | EDT Jar Lock for Glass Vessel Suitable for Composite Top Plate (Set of 16) |
| 27 | 0104A00262 | EDT Bath Rod |
| 28 | 0104A00394 | EDT-08L/08Lx Aluminum Top Plate for 250ml |
| 29 | 0104A00155 | EDT Digital Depth Gauge |
| 30 | 0104A00771 | EDT Customized Cage Sinkers 9.70mm Height X 13.50mm ID (Set of 6) |
| 31 | 0104A00421 | EDT Customized Sinkers 14mm Internal Diameter X 34 mm Internal Length (Set of 6) |
| 32 | 0104A00358 | EDT 33mm Paddle Over Disc Assembly |

| | | |
|----|-------------|--|
| 33 | 0104A00359 | EDT 50mm Paddle Over Disc Assembly |
| 34 | 0104A00360 | EDT 55mm Paddle Over Disc Assembly |
| 35 | 0104A00361 | EDT 65mm Paddle Over Disc Assembly |
| 36 | 0104A00353 | EDT Laser Marked Paddle for 33mm Disc (Set of 6) |
| 37 | 0104A00354 | EDT Laser Marked Paddle for 49mm Disc (Set of 6) |
| 38 | 0104A00355 | EDT Laser Marked Paddle for 55mm Disc (Set of 6) |
| 39 | 0104A00356 | EDT Laser Marked Paddle for 65mm Disc (Set of 6) |
| 40 | 0104A00357 | EDT Laser Marked Paddle for 85mm Disc (Set of 6) |
| 41 | 0104A00093 | EDT Moulded Tray For Paddles |
| 42 | 0104A00168 | EDT Paddle Over Disc 33 |
| 43 | 0104A00169 | EDT Paddle Over Disc 50 (Set of 6) |
| 44 | 0104A00170 | EDT Paddle Over Disc 55 (Set of 6) |
| 45 | 0104A00171 | EDT Paddle Over Disc 65 (Set of 6) |
| 46 | 0104A00172 | EDT Paddle Over Disc 85 (Set of 6) |
| 47 | 0104A00317 | EDT Paddles – 1000ml For Stationary Intrinsic Attachment (Set of 6) |
| 48 | 0104A00085 | EDT Paddles For 1000 MI Glass Vessel (Set of 6) |
| 49 | 0104A00088 | EDT Paddles For Enhancer Cell (Set of 6) |
| 50 | 0104A00526 | EDT Paddles For 1000 ml Vessel (Set of 8) |
| 51 | 0104A00568 | EDT Paddles For 1000 MI Glass Vessel |
| 52 | 0104A00576 | EDT 90mm Watch Glass Disc Assembly With 120 Mesh Screen Without Paddle |
| 53 | 0104A005894 | EDT Teflon Paddle |
| 54 | 0104A00630 | EDT Laser Marked Paddles For 1000 MI Glass Vessel |
| 55 | 0104A00637 | EDT 55mm Paddle Over Disc Assembly With 120 Mesh(set of 6) |
| 56 | 0104A00699 | EDT 33mm Paddle Over Disc Assembly With 20 Mesh (set of 6) |
| 57 | 0104A00702 | EDT 33mm Paddle Over Disc Assembly With 10 Mesh (set of 6) |
| 58 | 0104A00704 | EDT 55mm Paddle Over Disc Assembly With 20 Mesh(set of 6) |
| 59 | 0104A00703 | EDT 55mm Paddle Over Disc Assembly With 10 Mesh(set of 6) |
| 60 | 0104A00765 | EDT 76mm Paddle Over Disc Assembly With Standard Mesh(Set Of 6) |
| 61 | 0104A00124 | EDT Digital Thermometer |

| | | |
|----|---------------|--|
| 62 | 0164A00001 | EDT I-Disso 6U |
| 63 | 0104A00010S19 | Laser Marking Charges For EDT Kone Vessel 1000ml |
| 64 | 0104A00452 | EDT Hybrid Cannula for Sampling Manifold (Set of 6) |
| 65 | 0104A00070 | EDT Basket LID 10 Mesh (Set of 6) |
| 66 | 0104A00016 | EDT Glass Vessel 'V' Type 1000ml - Clear |
| 67 | 0104A00310 | EDT-08L/08Lx Bath Top Plate – Aluminum |
| 68 | 0104A00068 | EDT Basket Holder with Clip (Set of 6) |
| 69 | 0104A00014 | EDT Glass Vessel 'D' Type 1000ml - Clear |
| 70 | 0104A00162 | EDT- 08L/EDT- 08Lx Water Bath With Quick Release Coupling With Shut Off Valve Connector |
| 71 | 0104A00622 | EDT Japanese Sinker |
| 72 | 0104A00463 | EDT Sampling Cannula Rubber 'O' Ring for Dispenser plate (Set of 10) |
| 73 | 0104A00331 | EDT-08L/08Lx Blackout Curtain for I-Disso 6U |
| 74 | 0104A00097 | EDT Sinker 11 mm ID x 31 mm OD (Set of 6) |
| 75 | 0104A00382 | EDT-08L/08Lx Aluminum Top Plate for 150ml Conversion Kit |
| 76 | 0104A00071 | EDT Basket LID 40 Mesh (Set of 6) |
| 77 | 0104A00239 | EDT Digital Validation Kit |
| 78 | 0104A00330 | EDT Cooling Instrument (Temp. From 10° to 40° C) |
| 79 | 0104A00139 | EDT Jar Lock For Merlon Vessel For Composite Top Plate |
| 80 | 0104A00163 | EDT Teflon Coated Sinker 22 mm Long & 13 cm Length of Gauge Wire |
| 81 | 0104A00256 | EDT Adjustable Sampling Cannulas For 500 ML |
| 82 | 0122A00006 | EDT Sampling Cannulas for 2 Ltr Vessel |
| 83 | 0104A00260 | EDT-08L/08Lx Tablet Dispenser Plate Without Tubing |
| 84 | 0104A00285 | EDT-08L/08Lx – Enhancer Cell Kit- 4 CM2 |
| 85 | 0104A00109 | EDT Butterfly Sinkers(Set of 6) |
| 86 | 0104A00129 | EDT KBR Bench Press |
| 87 | 0104A00234 | EDT Low Evaporation Jar Lid for Adjustable Cannulas Suitable for Glass Vessel (Set of 6) |
| 88 | 0104A00242 | EDT Merlon Vessel 1000ml - Clear (N) |
| 89 | 0104A00318 | EDT-08L/08Lx -Enhancer Cell Kit 2CM |
| 90 | 0104B00013 | EDTV-1234 Bath Temperature Sensor |
| 91 | 0104A00159 | EDT-08L/08Lx Tablet Dispenser Plate With Tubing |
| 92 | 0104A00133 | EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 16) |

| | | |
|-----|------------|--|
| 93 | 0104A00483 | EDT Silicon Tube Set For ETC-11 |
| 94 | 0104A00482 | EDT Plunger Dial Gauge |
| 95 | 0104A00669 | EDT-08L/08Lx Light with supply Adaptor for I-Disso 6U |
| 96 | 0104A00672 | EDT Pin Sinker Type - B (Set of 12) |
| 97 | 0104A00678 | EDT Customized Japanese sinker ID 15mm x IL 32mm (Set of 6) |
| 98 | 0104A00346 | EDT-08L/08Lx 150 ML Low Volume Conversion Kit (Aluminum) |
| 99 | 0104A00419 | EDT-08L/08LX IMMERSION (ENHANCER) CELL KIT- 0.5 CM2 For Acrylic Bath |
| 100 | 0104A00667 | EDT Basket Sinker 40 Mesh (19mm Height) with Lid (set of 6 |
| 101 | 0104A00029 | EDT Laser Marked Kone Vessel 1000 ML - Clear |
| 102 | 0104A00030 | EDT Laser Marked Kone Vessel 1000 ML - Amber |
| 103 | 0104A00252 | EDT Laser Marked Basket 20 Mesh |
| 104 | 0104A00280 | EDT Basket 20 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate |
| 105 | 0104A00329 | EDT Luer Lock FTLLC-6 Male (Set of 8) |
| 106 | 0104A00278 | EDT – Intrinsic Plate |
| 107 | 0104A00225 | EDT Merlon Vessel 1000 ML - Clear with Jar Lock Suitable for S.S. Top Plate |
| 108 | 0104A00258 | EDT Merlon Vessel 1 Ltr Amber With Jar Centring Ring For Composite Top Plate |
| 109 | 0104A00157 | EDT Sampling Cannulas Fixed Height 750 ml (Set of 6) Length 108 mm |
| 110 | 0104A00156 | EDT Sampling Cannulas Fixed Height 900 ML (Set of 6) Length 98 mm |
| 111 | 0104A00049 | EDT Sampling Cannulas - Peak (Set of 6) |
| 112 | 0104A00050 | EDT Disposable Filters (PKT of 100) |
| 113 | 0104A00094 | EDT Sinker 5 mm ID x 15.5 mm OD (Set of 6) |
| 114 | 0104A00050 | EDT Disposable Filters (PKT of 100) |
| 115 | 0104A00095 | EDT Sinker 6 mm ID x 18 mm OD (Set of 6) |
| 116 | 0104A00096 | EDT Sinker 8 mm ID x 23 mm OD (Set of 6) |
| 117 | 0104A00098 | EDT Japanese Sinker (Set of 6) |
| 118 | 0104A00250 | EDT V-1234 Tygon Tube With SS Sleeve (Set of 6) |
| 119 | 0104A00373 | EDT Luer Lock BSFTL6-210-6 Female (Set of 8) |
| 120 | 0104A00435 | EDT Sampling Cannula Fixed Height 250 ML (Set of 6) Length 65mm |
| 121 | 0104A00432 | EDT Clearview (230V) Version-2 |
| 122 | 0104A00010 | EDT Kone Vessel 1000ml - Clear |

| | | |
|-----|-------------|---|
| 123 | 0104A00056 | EDT Basket 30 Mesh (Set of 6) |
| 124 | 0104A00718 | Re-Coating For EDT/08LX/08L |
| 125 | 0104A00715 | Re-Coating For EDT/14LX/14L |
| 126 | 0104A00720 | EDT-08LX Lift Cover Set Top & Bottom. |
| 127 | 0104A00723 | EDT Deconex 12 basic (7ltr.) |
| 128 | 0104A00705 | EDT Rubber Holder For Disposable Filter (Set Of 6) |
| 129 | 0104A00245 | EDT Customized Japanese Sinkers 15 mm ID X 25-26 mm Height |
| 130 | 0104A00705 | EDT Rubber Holder For Disposable Filter (Set Of 6) |
| 131 | 0104A00309 | EDT Tygon Tube for Drainage |
| 132 | 0104A00674 | EDT Basket 40 Mesh Teflon Coated (Set of 6) |
| 133 | 0104A00381 | EDT Teflon Coated Basket 20 Mesh (Set of 6) |
| 134 | 0104A00047 | EDT Adjustable Sampling Cannula For Low Evaporation Jar LID (Set of 6) |
| 135 | 0104A00257 | EDT Adjustable Sampling Cannulas For 900 ML (Set of 6) |
| 136 | 0104A00024 | EDT Laser Marked Glass 1000 ML - Amber |
| 137 | 0104A00282 | TDT-08L/EDT-08Lx – Belt Tightener |
| 138 | 0104A00671 | EDT Pin Type Sinkers - A (Set of 12) |
| 139 | 0104A00737 | EDT Deconex 12Basic (Pack of 10 Bottle each of 0.9L/1.3Kg) |
| 140 | 0104A00150 | EDT Tygon Tube for ETC-11 |
| 141 | 0104A00249 | EDT Jar Lid for 100ml Vessels |
| 142 | 0104A00570 | EDT Cooling Instrument (Temp. From 10° to 40° C) - 3 litre |
| 143 | 0104A00741 | Up-gradation from EDT- 08LX Basic to Offline Sampling Accessories With DC Sample collector |
| 144 | 0104A00389 | EDT Jar Lid for 250ml Vessels |
| 145 | 0104A00385 | EDT Jar Lock for Glass Vessel Suitable for Aluminum Top Plate (Set of 16) |
| 146 | 0104A00101 | EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big (Set of 6) |
| 147 | 0104A00681 | Tygon Tube Replacement And Performance Enhancement Kit For EDT-08LX With D Sample Collector |
| 148 | 0104A00401 | EDT Laser Marked Glass Vessel 1000 ML - Amber With Sigma Delta Certificate |
| 149 | 0104A00005 | EDT Glass Vessel 1000ml - Amber |
| 150 | 0104XA00004 | EDT Glass Vessel 1000ml - Clear |
| 151 | 0104A00023 | EDT Laser Marked Glass Vessel 1000 ML - Clear |

| | | |
|-----|------------|--|
| 152 | 0104A00004 | EDT Glass Vessel 1000ml - Clear |
| 153 | 0104A00629 | EDT Laser Marked Basket 40 Mesh |
| 154 | 0104A00567 | EDT Basket 40 Mesh |
| 155 | 0104A00300 | EDT Basket 40 Mesh for 150ml Vessel (Set of 6) |
| 156 | 0104A00443 | EDT -8 Mesh Felodipine Basket with Spacer (Set of 6) |
| 157 | 0104A00077 | EDT Laser Marked Basket 80 Mesh (Set of 6) |
| 158 | 0104A00299 | EDT Basket Holder with 'O' Ring (Set of 6) for 150ml Vessel |
| 159 | 0104A00365 | EDT Digital Vernier Caliper for 1000 MI Vessel (150mm) |
| 160 | 0104A00212 | EDT-08L/08Lx 100 MI Low Volume Conversion Kit |
| 161 | 0122A00019 | EDT Sampling Manifold Cannula S.S for 2Ltr (Set of 6) |
| 162 | 0104A00118 | EDT Validation Kit - Analog |
| 163 | 0104A00069 | EDT Basket Holder with 'O' Ring (Set of 6) |
| 164 | 0104A00451 | EDT Digital Dial Gauge |
| 165 | 0104A00515 | EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 8) |
| 166 | 0104A00709 | Offline Sampling Accessories for EDT-8Lx Dissolution Tester with Dx Sample Collector |
| 167 | 0104A00064 | EDT Felodipine Basket with Spacer (Set of 6) |
| 168 | 0104A00076 | EDT Laser Marked Basket 60 Mesh (Set of 6) |
| 169 | 0104A00078 | EDT Laser Marked Basket 100 Mesh (Set of 6) |
| 170 | 0104A00283 | EDT-08L/08Lx – Adjustable Legs (Set of 4) |
| 171 | 0104A00231 | EDT Sampling Cannula Fixed Height 500 MI (Set of 6) Length 108 mm |
| 172 | 0104A00042 | EDT Low Evaporation Jar LID (Set of 6) |
| 173 | 0104A00180 | EDT Tube Connecting Coupling With Nozzle (Set of 2)(Replacement Of QRC) |
| 174 | 0104A00307 | EDT-08Lx Side Circulation Water Bath - Clear with Assembly |
| 175 | 0104A00766 | EDT Transdermal Rod For Cylinder - Height 107.9 mm (Set of 8) |
| 176 | 0104A00181 | EDT Tube Connecting Coupling (Set of 2)(Replacement Of QRC) |
| 177 | 0104A00404 | EDT 33mm Watchglass Disc Assembly With 16 Mesh Screen |
| 178 | 0104A00405 | EDT 55mm Watchglass Disc Assembly With 16 Mesh Screen |
| 179 | 0104A00404 | EDT 33mm Watchglass Disc Assembly With 16 Mesh Screen |

| | | |
|-----|------------|--|
| 180 | 0104A00406 | EDT 90mm Watchglass Disc Assembly With 16 Mesh Screen |
| 181 | 0105A00133 | EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 28) |
| 182 | 0104A00579 | EDT-08L/08Lx 1 Ltr. Conversion Kit for 250 ML Vessels |
| 183 | 0104A00577 | EDT-08L/08Lx 1 Ltr. Conversion Kit for 150 ML Vessels |
| 184 | 0104A00571 | EDT Cooling Instrument (Temp. From 10° to 40° C) - 6 litre |
| 185 | 0104A00566 | EDT Basket Holder with 'O' Ring |
| 186 | 0104A00534 | EDT-08L/08Lx -Motorized Sampling Manifold without Accessories |
| 187 | 0104A00525 | EDT 33mm Watch glass Disc.Assembly with 120 Mesh Screen |
| 188 | 0104A00524 | EDT 90mm Watch glass Disc.Assembly with 120 Mesh Screen |
| 189 | 0104A00523 | EDT 55mm Watch glass Disc.Assembly with 120 Mesh Screen |
| 190 | 0104A00522 | EDT Plastic 100ml Beakers For cleaning(Set of 8) |
| 191 | 0104A00521 | EDT-8LX Cleaning Tray |
| 192 | 0104A00520 | EDT-8/14LX Tablet Locater (set of 8) |
| 193 | 0104A00519 | EDT-8/14LX Sampling port seal (set of 8) |
| 194 | 0104A00506 | EDT-V 1234 Merlon Vessels set of 6 OLD type |
| 195 | 0104A00476 | EDT Customized Cage Type Sinkers Size 16mm ID x 13 mm HT With Spacer (Set of 6) |
| 196 | 0104A00475 | EDT Customized Cage Type Sinkers Size 12mm ID X 11mm Length (Set of 6) |
| 197 | 0104A00473 | EDT Replenishment Cannula (Set of 6) |
| 198 | 0104A00472 | EDT Laser Marked Kone Vessel 1000 ml - Clear With Sigma Delta Certificate |
| 199 | 0104A00471 | EDT Kone Vessel 1000 ml – Clear With Sigma Delta Certificate |
| 200 | 0104A00465 | EDT Customized Cage Type Sinkers Size 10mm ID Length x 10mm HT With Lid (Set of 6) |
| 201 | 0104A00464 | EDT Tablet Dispensing Derlin 'O' Ring for Dispenser plate (Set of 6) |
| 202 | 0104A00462 | EDT Manifold Rubber 'O' Ring for Dispenser plate (Set of 2) |
| 203 | 0104A00461 | EDT Customized Cage Type Sinkers Size 14mm ID x 13 mm HT With Spacer (Set of 6) |
| 204 | 0104A00460 | EDT Cage Type Sinkers - Small With Spacer (Set of 6) |

| | | |
|-----|------------|---|
| 205 | 0104A00459 | EDT Cage Type Sinkers - Big With Spacer (Set of 6) |
| 206 | 0104A00458 | EDT Customized Granules Dispensing Adapter For Dispenser Plate (Set of 12) |
| 207 | 0104A00457 | EDT Spacer For Cage Type Sinkers Size 22ID X 22mm HT with lid (Set of 6) |
| 208 | 0104A00456 | EDT Spacer For Small Cage Type Sinkers (Set of 6) |
| 209 | 0104A00455 | EDT Spacer For Big Cage Type Sinkers (Set of 6) |
| 210 | 0104A00453 | EDT-08L/08Lx – Auto Dispenser Plate |
| 211 | 0104A00447 | EDT 40 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6) |
| 212 | 0104A00446 | EDT Customized Cage Type Sinkers – Medium With Lid (Set of 6) |
| 213 | 0104A00444 | EDT Teflon Coated Basket 10 Mesh (Set of 6) |
| 214 | 0104A00441 | EDT Basket LID 20 Mesh (Set Of 6) |
| 215 | 0104A00437 | EDT Bath Nozzle Rubber Gasket (Set of 2) |
| 216 | 0104A00434 | EDT-08L/08Lx Bath Guard (Set of 2) |
| 217 | 0104A00433 | EDT Silicon Tube Set For ETC-15 |
| 218 | 0104A00429 | EDT-08Lx Conversion Kit for Enhancer Cell to 1Ltr Vessels |
| 219 | 0104A00428 | EDT Basket Rod for 250ml Vessel |
| 220 | 0104A00425 | EDT Basket Sinkers 40 Mesh (9.00mm) with Lid (Set of 6) |
| 221 | 0104A00424 | EDT Disposable Filter 3.00mm (Pkts of 100) |
| 222 | 0104A00423 | EDT-08L/08Lx Support Strip for Bath (Set of 2) |
| 223 | 0104A00422 | EDT-08Lx Top Cover |
| 224 | 0104A00415 | EDT Big Retaining Ring for for Enhancer Cell Kit – 4cm (Set of 8) |
| 225 | 0104A00414 | EDT Top Cover with Base Plate for ETC-11LX |
| 226 | 0104A00413 | EDT Sampling Cannula Fixed Height 250 ML (Set of 6) Length 71 mm |
| 227 | 0104A00412 | EDT Luer Lock PMLP-21 (Male) (Set of 8) |
| 228 | 0104A00411 | EDT Silicon Tube for Drainage |
| 229 | 0104A00409 | EDT 120 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6) |
| 230 | 0104A00408 | EDT 120 Mesh Screen For 55mm Watchglass Disc Assembly (Set of 6) |
| 231 | 0104A00407 | EDT 120 Mesh Screen For 33mm Watchglass Disc Assembly (Set of 6) |
| 232 | 0104A00402 | EDT Luer Lock FLLR-1 (Female) (Set of 8) |
| 233 | 0104A00397 | EDT-08L/08Lx – Water Circulation Silicon Tube Set |

| | | |
|-----|------------|---|
| 234 | 0104A00393 | EDT-08L/08Lx PTFE Carrier Tube From UV to Dissolution (OD -1.6mm / ID-0.8mm) - (Set Of 8) |
| 235 | 0104A00392 | EDT-08L/08Lx PTFE Carrier Tube From Agilent Pump To UV (OD-1.6mm / ID-0.8mm) - (Set Of 8) |
| 236 | 0104A00391 | EDT-08L/08Lx PTFE Carrier Tube From Disso to Agilent Pump (OD-1.6mm / ID-0.8mm) -(Set of 8) |
| 237 | 0104A00390 | EDT Washer for Stationary Basket (Set of 2) |
| 238 | 0104A00388 | EDT Analog Pressure Gauge For KBR Press |
| 239 | 0104A00388 | EDT Analog Pressure Gauge For KBR Press |
| 240 | 0104A00386 | EDT-08L/08Lx Aluminum Top Plate for Enhancer Cell Kit |
| 241 | 0104A00384 | EDT Luer Lock FTLL-6 Female (Set of 8) |
| 242 | 0104A00379 | EDT Glass Vessel Trolley |
| 243 | 0104A00378 | EDT Customized Transparent Moulded Jar Lid - Plastic (Set of 6) |
| 244 | 0104A00377 | EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 6) |
| 245 | 0104A00376 | EDT Sampling Cannula Fixed Height 500 ML (Set of 6) Length 127 mm |
| 246 | 0104A00375 | EDT-08L/08Lx – Tablet Dispenser Plate with Provision for Felodipine Basket |
| 247 | 0104A00374 | EDT Customized Helix Coil Sinkers 9.5mm ID X 20mm Height (Set of 6) |
| 248 | 0104A00372 | EDT Hardened Polished Surface Plate (for Intrinsic) |
| 249 | 0104A00371 | EDT Bolts for Fixing the Die and Surface Plate (for Intrinsic)(Set of 3) |
| 250 | 0104A00370 | EDT Stand For 1000ml Vessel (Composite) |
| 251 | 0104A00369 | EDT Deconex CIP Power-X (14 Kg) |
| 252 | 0104A00368 | EDT Online Connector (Set of 8) |
| 253 | 0104A00364 | EDT Glass Vessel 'D' Type 1000 ML - Amber With Ring |
| 254 | 0104A00351 | EDT 15ml Replenishment Bottle With Rubber Cap |
| 255 | 0104A00350 | EDT Luer Lock MLRL210-6 Male (Set of 8) |
| 256 | 0104A00348 | EDT Customized 10 Mesh Basket Sinkers with Lid (12.7mm ID – 24.5mm Height) (Set of 6) |
| 257 | 0104A00347 | EDT Customized 10 Mesh Basket Sinkers with Lid (9.1mm ID – 18.8mm Height) (Set of 6) |
| 258 | 0104A00345 | EDT Prednisone Tablet (Pkt of 30) |

| | | |
|-----|------------|---|
| 259 | 0104A00344 | EDT Prednisone Standard USP-1559006 Reference 250 Mg |
| 260 | 0104A00343 | EDT Disc-Filter-Millex-HN, 33mm, 0.45um, Nylon Membrane, Non-Sterile (Pkt of 250) |
| 261 | 0104A00342 | EDT Disc Filter-Millex -HN, 33mm, 0.45um PES Membrane, Non-Sterile (Pkt of 250) |
| 262 | 0104A00339 | EDT-08L/08Lx PTFE Carrier Tube From Disso to Watson Marlow Pump (Set of 8) |
| 263 | 0104A00337 | EDT-08L/08Lx PTFE Carrier Tube From Disso to Ismatec Pump (Set of 8) |
| 264 | 0104A00336 | EDT Probe for Vibration Meter |
| 265 | 0104A00335 | EDT Sampling Cannula Fixed Height For Enhancer Cell Kit -200ml (Set of 6) |
| 266 | 0104A00333 | EDT Cylindrical Flat Bottom Glass Vessel 200ml – Clear for Enhancer Cell (Set of 6) |
| 267 | 0104A00327 | EDT Enhancer Cell for Enhancer Cell Kit – 2 CM (Set of 6) |
| 268 | 0104A00325 | EDT-08L/08Lx -Motorised Sampling Manifold (For Peristaltic Offline) |
| 269 | 0104A00324 | EDT Sampling Manifold Cannula S.S. (SET Of 6) |
| 270 | 0104A00323 | EDT-08L V-34 Syringe Manifold (for Offline Syringe) |
| 271 | 0104A00322 | EDT-08L/08Lx PTFE Carrier Tube For Syringe Manifold (From Sampling Probe To Syringe Manifold) (Set Of 6) |
| 272 | 0104A00321 | EDT Luer Lock BSMLRL-6 Male (Set of 8) |
| 273 | 0104A00320 | EDT Enhancer Cell for Enhancer Cell Kit - 4 CM (Set of 6) |
| 274 | 0104A00316 | EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6) |
| 275 | 0104A00314 | EDT Sinkers 20 Mesh with Lid (Set of 6) |
| 276 | 0104A00316 | EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6) |
| 277 | 0104A00312 | EDT-08L/08Lx S.S. Strip for Composite Top Plate (Set of 5) |
| 278 | 0104A00288 | EDT Basket Sinkers 40 Mesh (11.7mm Height) with Lid (Set of 6) |
| 279 | 0104A00284 | EDT Sinkers 10 Mesh with Lid (Set of 6) |
| 280 | 0104A00279 | EDT Basket 10 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate |
| 281 | 0104A00277 | EDT-08L/08Lx -Motorised Sampling Manifold (For Syringe Offline) |
| 282 | 0104A00276 | EDT Deconex-CIP Fresh |
| 283 | 0104A00275 | EDT Deconex-CIP Alpha- X |

| | | |
|-----|------------|---|
| 284 | 0104A00274 | EDT U Glass Tube |
| 285 | 0104A00273 | EDT Teflon Coated Basket Shaft With O Ring & Holder |
| 286 | 0104A00271 | EDT Rubber Ring For Low Evaporation Jar Lid |
| 287 | 0104A00270 | EDT-08L/EDT-08Lx Composite Top Plate for 100ml Low Volume Conversion Kit |
| 288 | 0104A00269 | EDT Spirit Level |
| 289 | 0104A00268 | EDT Top Cover with Base Plate for ETC-11 |
| 290 | 0104A00266 | EDT Basket 100 Mesh "V" Type |
| 291 | 0104A00264 | EDT Basket 40 Mesh "H" Type |
| 292 | 0104A00263 | EDT Basket 100 Mesh "H" Type |
| 293 | 0104A00261 | EDT S.S Knob for Bath Rod |
| 294 | 0104A00259 | EDT Merlon Vessel 1000ml Amber (N) |
| 295 | 0104A00254 | EDT SS Filters with Holders – 40 Micron (Set of 6) |
| 296 | 0104A00255 | EDT SS Filters with Holders – 80 Micron (Set of 6) |
| 297 | 0104A00253 | EDT Rubber Holder For SS Filter |
| 298 | 0104A00248 | EDT Peak Glass Lense to Mesure Basket Mesh |
| 299 | 0104A00246 | EDT Polished Punch |
| 300 | 0104A00241 | EDT Deconex 24 Liq (7 Kg) |
| 301 | 0104A00240 | EDT Deconex 25 Organ Acid (5 Kg) |
| 302 | 0104A00237 | EDT Vernier Caliper with Stopper for 1000ml Vessels |
| 303 | 0104A00235 | EDT Low Evaporation Jar Lid for Sampling Cannulas Fixed Height Suitable for Glass Vessel (Set of 6) |
| 304 | 0104A00232 | EDT Sampling Cannula Fixed Height 100 ML (Set of 6) Length 73.5 mm |
| 305 | 0104A00229 | EDT Dispenser Jar Ring Black (Set of 8) |
| 306 | 0104A00228 | EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate |
| 307 | 0104A00228 | EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate |
| 308 | 0104A00237 | EDT Vernier Caliper with Stopper for 1000ml Vessels |
| 309 | 0104A00229 | EDT Dispenser Jar Ring Black (Set of 8) |
| 310 | 0104A00228 | EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate |
| 311 | 0104A00227 | EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for S.S. Top Plate |
| 312 | 0104A00226 | EDT Merlon 1000 ML - Clear with Jar Centering Ring Suitable for Composite Top Plate |
| 313 | 0104A00224 | EDT Spindle Latch with Spring (Set of 8) |

| | | |
|-----|------------|---|
| 314 | 0104A00220 | EDT Japanese Sinkers with 40 Mesh Cover (Set of 6) |
| 315 | 0104A00219 | EDT Sinkers 40 Mesh with Lid (Set of 6) |
| 316 | 0104A00218 | EDT Customized Japanese Sinkers 15 mm ID X 40 mm Height (Set of 6) |
| 317 | 0104A00206 | EDT Holder for Basket 40 Mesh 24.5 mm (Set of 6) |
| 318 | 0104A00203 | EDT Deconex 12 Basic (500ml) |
| 319 | 0104A00185 | EDT Stop Watch |
| 320 | 0104A00184 | EDT Transdermal Cylinder - Height 50.8 mm (Set of 6) |
| 321 | 0104A00177 | EDT 'D' Type Basket 40 Mesh (Set of 6) |
| 322 | 0104A00174 | EDT 'D' Type Basket 10 Mesh (Set of 6) |
| 323 | 0104A00167 | EDT Luer Lock FTLL210-6 Female (Set of 8) |
| 324 | 0104A00166 | EDT Luer Lock MTLL220-6 Male (Set of 8) |
| 325 | 0104A00160 | EDT Gauge Block 25 mm |
| 326 | 0104A00154 | EDT Basket Validation Gauge with Certificate |
| 327 | 0104A00152 | EDT Moulded Tray For Basket, Basket Holder & Rod |
| 328 | 0104A00145 | EDT Plastic Quick Release Coupling with Nozzle (Set of 2) |
| 329 | 0104A00141 | EDT Vessel Centering Ring for Merlon Vessel Suitable for Composite Top Plate (Set of 8) |
| 330 | 0104A00140 | EDT Vessel Centering Ring for Glass Vessel Suitable for Composite Top Plate (Set of 8) |
| 331 | 0104A00136 | EDT-08L/08Lx Composite Top Plate for Glass Vessel |
| 332 | 0104A00135 | EDT Vessel Centering Plate for Merlon Vessel Suitable for S.S. Top Plate |
| 333 | 0104A00134 | EDT TNT Washer for Merlon Vessel Suitable for S.S. Top Plate (Set of 16) |
| 334 | 0104A00132 | EDT Jar Lock for Glass Vessel Suitable for S.S. Top Plate (Set of 16) |
| 335 | 0104A00131 | EDT-08L/08Lx S.S. Top Plate for Merlon Vessel |
| 336 | 0104A00130 | EDT-08L/08Lx S.S. Top Plate for Glass Vessel |
| 337 | 0104A00125 | EDT Digital Verticality Meter |
| 338 | 0104A00123 | EDT Thermometer - Mercury |
| 339 | 0104A00122 | EDT Digital Tachometer |
| 340 | 0104A00121 | EDT Depth Gauge for 1000 ML Vessel |
| 341 | 0104A00120 | EDT Centering Gauge |
| 342 | 0104A00119 | EDT Centering Device for 1000 ML Vessel |
| 343 | 0104A00116 | EDT Helix Coil Sinkers - CAPWHT-2S (Set of 6) |
| 344 | 0104A00119 | EDT Centering Device for 1000 ML Vessel |
| 345 | 0104A00115 | EDT Helix Coil Sinkers - CAPWHT-02 (Set of 6) |

| | | |
|-----|------------|---|
| 346 | 0104A00114 | EDT Metformin Sinkers - Small (Set of 6) |
| 347 | 0104A00113 | EDT Metformin Sinkers - Big (Set of 6) |
| 348 | 0104A00110 | EDT Sinkers 'U' Type (Set of 6) |
| 349 | 0104A00107 | EDT Alfuzosin Sinkers 37.5 mm ID x 20 mm OD (Set of 6) |
| 350 | 0104A00105 | EDT 3 Prong Sinkers - SS (Set of 6) |
| 351 | 0104A00104 | EDT 3 Prong Sinkers - Plastic (Set of 6) |
| 352 | 0104A00102 | EDT Spring Type Sinkers 9 mm ID x 30 mm OD - Small (Set of 6) |
| 353 | 0104A00100 | EDT Cage Type Sinkers - Small (Set of 6) |
| 354 | 0104A00081 | EDT Laser Marked Felodipine Basket (Set of 6) |
| 355 | 0104A00080 | EDT Laser Marked Stationary Basket (Set of 6) |
| 356 | 0104A00075 | EDT Laser Marked Basket 40 Mesh (Set of 6) |
| 357 | 0104A00072 | EDT Laser Marked Basket 10 Mesh (Set of 6) |
| 358 | 0104A00067 | EDT Basket 'O' Ring (Set of 6) |
| 359 | 0104A00066 | EDT Basket Rod |
| 360 | 0104A00065 | EDT Moulded Tray For Basket |
| 361 | 0104A00062 | EDT Stationary Basket (Set of 6) |
| 362 | 0104A00060 | EDT Basket 100 Mesh (Set of 6) |
| 363 | 0104A00059 | EDT Basket 80 Mesh (Set of 6) |
| 364 | 0104A00058 | EDT Basket 60 Mesh (Set of 6) |
| 365 | 0104A00055 | EDT Basket 20 Mesh (Set of 6) |
| 366 | 0104A00055 | EDT Basket 20 Mesh (Set of 6) |
| 367 | 0104A00052 | EDT Inline Filters - 0.45 Micron (Set of 6) |
| 368 | 0104A00051 | EDT SS Filters with Holders (Set of 6) |
| 369 | 0104A00045 | EDT Sampling Cannulas Fixed Height (Set of 6) |
| 370 | 0104A00044 | EDT Jar LID For Staggered Start - White (Set of 6) |
| 371 | 0104A00041 | EDT Moulded Jar LID - Plastic (Set of 6) |
| 372 | 0104A00040 | EDT Rubber Stand For Glass Vessel |
| 373 | 0104A00037 | EDT Laser Marked Glass Vessel With Ring 1000ml - Clear |
| 374 | 0104A00035 | EDT Laser Marked Glass Vessel 'V' Type 1000ml - Clear |
| 375 | 0104A00033 | EDT Laser Marked Glass Vessel 'D' Type 1000 ml - Clear |
| 376 | 0104A00032 | EDT Laser Marked Flat Bottom 1000ml - Amber |
| 377 | 0104A00032 | EDT Laser Marked Flat Bottom 1000ml - Amber |
| 378 | 0104A00022 | EDT Laser Marked Glass Vessel 250 ml - Amber |
| 379 | 0104A00022 | EDT Laser Marked Glass Vessel 250 ml - Amber |
| 380 | 0104A00018 | EDT Glass Vessel With Ring 1000ml - Clear |
| 381 | 0104A00017 | EDT Glass Vessel V Type 1000 ml - Amber |

| | | |
|-----|-------------|---|
| 382 | 0104A00012 | EDT Flat Bottom 1000 ml - Clear |
| 383 | 0104A00127 | EDT Vibration Meter |
| 384 | 0104A00484 | EDT Spring Assembly for Tablet Dispenser Plate (Set of 2) |
| 385 | 0104A00442 | EDT Deconex 12 Basic (10 Bottles of 1000 ml) |
| 386 | 0104A00623 | EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big |
| 387 | 0104A00063 | EDT Felodipine Basket (Set of 6) |
| 388 | 0104A00584 | EDT 40 Mesh Stationary Basket (Set of 6) |
| 389 | 0104A00057 | EDT Basket 40 Mesh (Set of 6) |
| 390 | 0104A00410 | EDT Clearview (230V) |
| 391 | 0104A00267 | EDT Basket 20 Mesh "V" Type |
| 392 | 0104A00175 | EDT 'D' Type Basket 20 Mesh (Set of 6) |
| 393 | 0104A00265 | EDT Basket 20 Mesh "H" Type |
| 394 | 0104A00613 | EDT-08Lx PTFE Bath Top Plate |
| 395 | 0104A00581 | EDT 16 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6) |
| 396 | 0104A00692 | Inspire-08 Manual Syringe Sampling Attachment |
| 397 | 0104A00763 | Inspire-8 Tablet dispenser plate |
| 398 | 0104A005891 | Embedded Sliding drawer for Paddle and Basket Inspire-8/TrustE-8 all Configuration |
| 399 | 0104A00785 | Connexion+ |