

8 Station Dissolution Tester

Manual Syringe Sampling



Inspire-8 Manual Syringe Sampling

Features:

- LCD display of 40 characters and 4 lines for user friendly operation
- Multilevel security with password
- Alphanumeric data entry of product and user details
- Unique TimeAction™ function which allows media changeover and infinity test run, ideal for sustained and controlled release products
- Complies with USP, IP, Ph. Eur. and JP specifications
- Automatic stirrer height positioning for USP 1, 2, 5, 6, intrinsic test methods
- Easy Snap-Fit™ shaft locking mechanism for positive engagement and wobble free operation
- Clear water bath with metal top-plate
- Precise individual vessel centering system
- Adjustable leveling of instrument
- Vibration free smooth motorized lift movement
- Isolated water circulating pump for precise temperature control of the water bath
- 20 programmable protocols
- 24 programmable sampling intervals
- Easy to empty and clean water bath with bottom drainage
- Designed to minimize routine validation
- Report for online validation of test, monograph, error log during the test
- Programmable wake-up for temperature control
- Programmable sample volume from 1 mL to 20 mL
- Audio-visual indicators for system status
- Power failure detection facility
- Individual temperature probe for monitoring temperature of each vessel

Specifications:

Model	Inspire-8 Manual Syringe Sampling
Speed range	Variable speed from 20 to 300 RPM
Speed accuracy	±1 RPM
Temperature range	20°C to 40°C (Extended range 10°C to 40°C & 20°C to 55°C optional)
Temperature accuracy	±0.1°C
Display	40 x 4 character backlit LCD
No. of protocols	20 maximum
Sample interval	Maximum 24 of 1 min to 23 hrs, 59 min
Printer interface	Serial port
Printer supported	80 column Dot Matrix printer
Bath	Clear molded water bath
Dimensions (W x L x H)	517 x 712 x 665 mm
Heater	1 kW/230 VAC, material S.S. 304 with sensor
Vessel capacity	150 ml, 250 ml, 500 ml, 1 L
Electrical power	220/ 230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/10 Amp
Weight	68 KG

Ordering Information:

Products		
Product code	Description	Quantity
0129104	Inspire-08 with Manual Syringe Sampling - 230V	1
0119104	Inspire-08 with Manual Syringe Sampling - 110V	1

Standard accessories		
Part number	Description	Quantity
0104A00075	EDT laser marked basket - 40 mesh (Set of 6)	1
0104A00066	EDT basket rod	6
0104A00068	EDT Basket Holder with clip (Set of 6)	1
0104A00592	Dissolution Inspire Detachable Paddle Blade (Set of 6)	1
0104A00610	Stirrer for replenishment	2
0104A00242	EDT merlon vessel -1000ml clear (N)	8
0104BV4155	EDT V-5 External Probe Without Trim Pot	1
0104BV4114	EDT-Lx V-4 Temperature Controller ETC-15	1

Manual syringe sampling attachment		
Part number	Description	Quantity
0104A00277	Motorised sampling manifold	1
0104BV4050	Individual vessel temperature sensors with cable and connector	1
0104A00334	Syringe filter manifold	1
0104A00763	Inspire-8 Tablet dispenser plate	1

Assurance Plan		
Part number	Description	Quantity
0191S03	Silver Assurance Plan for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1
0191S01	Platinum Assurance Plan for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1
0191S02	Gold Assurance Plan for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1
0191S20	Bronze Assurance Plan for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

IQ/OQ Document		
Part number	Description	Quantity
0191S15	IQ/OQ Document for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

Single Validation Charges

Part number	Description	Quantity
0191S04	Single Validation Charges for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

Performance Qualification (PQ Documents)

Part number	Description	Quantity
0191S17	Performance Qualification (PQ Documents) for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Samplig	1

Reinstallation Charges

Part number	Description	Quantity
0191S11	Reinstallation Charges with IQ/OQ Document for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

Upgradation

Part number	Description	Quantity
0191S13	Upgradation for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

Inspection & Service Charges

Part number	Description	Quantity
0191S06	Inspection & Service Charges for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

ASTM Validation Charges

Part number	Description	Quantity
0191S16	ASTM Validation Charges for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

Dismantling Charges

Part number	Description	Quantity
0191S12	Dismantling Charges for Tablet Dissolution Tester Model Inspire-8 with Manual Syringe Sampling	1

*The information contained in this document is believed to be correct but **ELECTROLAB** accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB(India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

sr.no.	Part code	Optional accessories
1	0104A00583	EDT 20 Mesh Stationary Basket (Set of 6)
2	0104A00706	EDT-Inspire Syringe manifold Tygon Tube with SS Sleeve (Set Of 6)
3	0104A00693	EDT Rubber Holder For External Probe.(Set Of 6)
4	0104A00609	EDT 62mm Paddle Over Disc Assembly With 16 Mesh(set of 6)
5	0104A00625	EDT Moulded Jar LID - Plastic
6	0104A00631	EDT – Adjustable Legs Long (Set of 4)
7	0104A005893	EDT Teflon Paddles (set of 8)
8	0104A00596	EDT Paddle/Basket Rod
9	0104A00362	EDT 85mm Paddle Over Disc Assembly
10	0104A00142	EDT-08L/08Lx Water Bath with Plastic QRC - Clear
11	0104A00780	EDT-08 Bath Top Plate For Hexa Vessel
12	0104A00048	EDT Adjustable Sampling Cannula for low evaporation jar lid (Set of 8)
13	0104A00272	EDT Gasket For Tablet Dispensor Plate
14	0104A00112	EDT Helical Coil Sinkers
15	0104A00046	EDT Adjustable Sampling Cannulas For Tablet Input Device (Set of 6) Length 144 mm
16	0104A00205	EDT Lid for Basket 40 Mesh 24.5 mm (Set of 6)
17	0104A00182	EDT Inside Caliper 6"
18	0104A00782	EDT 90mm Watch Glass Disc (Set of 6)
19	0104A00400	EDT Laser Marked Glass Vessel 1000 ML - Clear With Sigma Delta Certificate
20	0122A00010	EDT Depth Gauge for 2 Ltr Vessels
21	0104A00204	EDT Basket 40 Mesh 24.5 mm (Set of 6)
22	0104A00128	EDT Intrinsic Attachment
23	0104A00146	EDT Plastic Quick Release Coupling with Shut Off Valve (Set of 2)
24	0104A00399	EDT Glass Vessel 1000 ML - Amber With Sigma Delta Certificate
25	0104A00398	EDT Glass Vessel 1000 ML – Clear With Sigma Delta Certificate
26	0104A00138	EDT Jar Lock for Glass Vessel Suitable for Composite Top Plate (Set of 16)
27	0104A00262	EDT Bath Rod
28	0104A00394	EDT-08L/08Lx Aluminum Top Plate for 250ml
29	0104A00155	EDT Digital Depth Gauge
30	0104A00771	EDT Customized Cage Sinkers 9.70mm Height X 13.50mm ID (Set of 6)
31	0104A00421	EDT Customized Sinkers 14mm Internal Diameter X 34 mm Internal Length (Set of 6)
32	0104A00358	EDT 33mm Paddle Over Disc Assembly

33	0104A00359	EDT 50mm Paddle Over Disc Assembly
34	0104A00360	EDT 55mm Paddle Over Disc Assembly
35	0104A00361	EDT 65mm Paddle Over Disc Assembly
36	0104A00353	EDT Laser Marked Paddle for 33mm Disc (Set of 6)
37	0104A00354	EDT Laser Marked Paddle for 49mm Disc (Set of 6)
38	0104A00355	EDT Laser Marked Paddle for 55mm Disc (Set of 6)
39	0104A00356	EDT Laser Marked Paddle for 65mm Disc (Set of 6)
40	0104A00357	EDT Laser Marked Paddle for 85mm Disc (Set of 6)
41	0104A00093	EDT Moulded Tray For Paddles
42	0104A00168	EDT Paddle Over Disc 33
43	0104A00169	EDT Paddle Over Disc 50 (Set of 6)
44	0104A00170	EDT Paddle Over Disc 55 (Set of 6)
45	0104A00171	EDT Paddle Over Disc 65 (Set of 6)
46	0104A00172	EDT Paddle Over Disc 85 (Set of 6)
47	0104A00317	EDT Paddles – 1000ml For Stationary Intrinsic Attachment (Set of 6)
48	0104A00085	EDT Paddles For 1000 ML Glass Vessel (Set of 6)
49	0104A00088	EDT Paddles For Enhancer Cell (Set of 6)
50	0104A00526	EDT Paddles For 1000 ml Vessel (Set of 8)
51	0104A00568	EDT Paddles For 1000 ML Glass Vessel
52	0104A00576	EDT 90mm Watch Glass Disc Assembly With 120 Mesh Screen Without Paddle
53	0104A005894	EDT Teflon Paddle
54	0104A00630	EDT Laser Marked Paddles For 1000 ML Glass Vessel
55	0104A00637	EDT 55mm Paddle Over Disc Assembly With 120 Mesh(set of 6)
56	0104A00699	EDT 33mm Paddle Over Disc Assembly With 20 Mesh (set of 6)
57	0104A00702	EDT 33mm Paddle Over Disc Assembly With 10 Mesh (set of 6)
58	0104A00704	EDT 55mm Paddle Over Disc Assembly With 20 Mesh(set of 6)
59	0104A00703	EDT 55mm Paddle Over Disc Assembly With 10 Mesh(set of 6)
60	0104A00765	EDT 76mm Paddle Over Disc Assembly With Standard Mesh(Set Of 6)
61	0104A00124	EDT Digital Thermometer
62	0164A00001	EDT I-Disso 6U
63	0104A00010S19	Laser Marking Charges For EDT Kone Vessel 1000ml
64	0104A00452	EDT Hybrid Cannula for Sampling Manifold (Set of
65	0104A00070	EDT Basket LID 10 Mesh (Set of 6)
66	0104A00016	EDT Glass Vessel 'V' Type 1000ml - Clear
67	0104A00310	EDT-08L/08Lx Bath Top Plate – Aluminum
68	0104A00068	EDT Basket Holder with Clip (Set of 6)

69	0104A00014	EDT Glass Vessel 'D' Type 1000ml - Clear
70	0104A00162	EDT- 08L/EDT- 08Lx Water Bath With Quick Release Coupling With Shut Off Valve Connector
71	0104A00622	EDT Japanese Sinkers
72	0104A00463	EDT Sampling Cannula Rubber 'O' Ring for Dispenser plate (Set of 10)
73	0104A00331	EDT-08L/08Lx Blackout Curtain for I-Disso 6U
74	0104A00097	EDT Sinkers 11 mm ID x 31 mm OD (Set of 6)
75	0104A00382	EDT-08L/08Lx Aluminum Top Plate for 150ml
76	0104A00071	EDT Basket LID 40 Mesh (Set of 6)
77	0104A00239	EDT Digital Validation Kit
78	0104A00330	EDT Cooling Instrument (Temp. From 10° to 40° C)
79	0104A00139	EDT Jar Lock For Merlon Vessel For Composite Top
80	0104A00163	EDT Teflon Coated Sinkers 22 mm Long & 13 cm Length of Gauge Wire
81	0104A00256	EDT Adjustable Sampling Cannulas For 500 ML
82	0122A00006	EDT Sampling Cannulas for 2 Ltr Vessel
83	0104A00260	EDT-08L/08Lx Tablet Dispenser Plate Without
84	0104A00285	EDT-08L/08Lx – Enhancer Cell Kit- 4 CM2
85	0104A00109	EDT Butterfly Sinkers (Set of 6)
86	0104A00129	EDT KBR Bench Press
87	0104A00234	EDT Low Evaporation Jar Lid for Adjustable Cannulas Suitable for Glass Vessel (Set of 6)
88	0104A00242	EDT Merlon Vessel 1000ml - Clear (N)
89	0104A00318	EDT-08L/08Lx -Enhancer Cell Kit 2CM
90	0104B00013	EDTV-1234 Bath Temperature Sensor
91	0104A00159	EDT-08L/08Lx Tablet Dispenser Plate With Tubing
92	0104A00133	EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 16)
93	0104A00483	EDT Silicon Tube Set For ETC-11
94	0104A00482	EDT Plunger Dial Gauge
95	0104A00669	EDT-08L/08Lx Light with supply Adaptor for I-Disso
96	0104A00672	EDT Pin Sinkers Type - B (Set of 12)
97	0104A00678	EDT Customized Japanese sinkers ID 15mm x IL 32mm (Set of 6)
98	0104A00346	EDT-08L/08Lx 150 ML Low Volume Conversion Kit (Aluminum)
99	0104A00419	EDT-08L/08Lx IMMERSION (ENHANCER) CELL KIT- 0.5 CM2 For Acrylic Bath
100	0104A00667	EDT Basket Sinkers 40 Mesh (19mm Height) with Lid
101	0104A00029	EDT Laser Marked Kone Vessel 1000 ML - Clear
102	0104A00030	EDT Laser Marked Kone Vessel 1000 ML - Amber
103	0104A00252	EDT Laser Marked Basket 20 Mesh

104	0104A00280	EDT Basket 20 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate
105	0104A00329	EDT Luer Lock FTLLC-6 Male (Set of 8)
106	0104A00278	EDT – Intrinsic Plate
107	0104A00225	EDT Merlon Vessel 1000 ML - Clear with Jar Lock Suitable for S.S. Top Plate
108	0104A00258	EDT Merlon Vessel 1 Ltr Amber With Jar Centring Ring For Composite Top Plate
109	0104A00157	EDT Sampling Cannulas Fixed Height 750 ml (Set of 6) Length 108 mm
110	0104A00156	EDT Sampling Cannulas Fixed Height 900 ML (Set of 6) Length 98 mm
111	0104A00049	EDT Sampling Cannulas - Peak (Set of 6)
112	0104A00050	EDT Disposable Filters (PKT of 100)
113	0104A00094	EDT Sinkers 5 mm ID x 15.5 mm OD (Set of 6)
114	0104A00050	EDT Disposable Filters (PKT of 100)
115	0104A00095	EDT Sinkers 6 mm ID x 18 mm OD (Set of 6)
116	0104A00096	EDT Sinkers 8 mm ID x 23 mm OD (Set of 6)
117	0104A00098	EDT Japanese Sinkers (Set of 6)
118	0104A00250	EDT V-1234 Tygon Tube With SS Sleeve (Set of 6)
119	0104A00373	EDT Luer Lock BSFTL6-210-6 Female (Set of 8)
120	0104A00435	EDT Sampling Cannula Fixed Height 250 ML (Set of 6) Length 65mm
121	0104A00432	EDT Clearview (230V) Version-2
122	0104A00010	EDT Kone Vessel 1000ml - Clear
123	0104A00056	EDT Basket 30 Mesh (Set of 6)
124	0104A00718	Re-Coating For EDT/08LX/08L
125	0104A00715	Re-Coating For EDT/14LX/14L
126	0104A00720	EDT-08LX Lift Cover Set Top & Bottom.
127	0104A00723	EDT Deconex 12 basic (7ltr.)
128	0104A00705	EDT Rubber Holder For Disposable Filter (Set Of 6)
129	0104A00245	EDT Customized Japanese Sinkers 15 mm ID X 25-26 mm Height
130	0104A00705	EDT Rubber Holder For Disposable Filter (Set Of 6)
131	0104A00309	EDT Tygon Tube for Drainage
132	0104A00674	EDT Basket 40 Mesh Teflon Coated (Set of 6)
133	0104A00381	EDT Teflon Coated Basket 20 Mesh (Set of 6)
134	0104A00047	EDT Adjustable Sampling Cannula For Low Evaporation Jar LID (Set of 6)
135	0104A00257	EDT Adjustable Sampling Cannulas For 900 ML (Set of 6)
136	0104A00024	EDT Laser Marked Glass 1000 ML - Amber
137	0104A00282	EDT-08L/EDT-08Lx – Belt Tightener
138	0104A00671	EDT Pin Type Sinkers - A (Set of 12)
139	0104A00737	EDT Deconex 12Basic (Pack of 10 Bottle each of 100ml)

140	0104A00150	EDT Tygon Tube for ETC-11
141	0104A00249	EDT Jar Lid for 100ml Vessels
142	0104A00570	EDT Cooling Instrument (Temp. From 10° to 40° C) -
143	0104A00741	Up-gradation from EDT- 08LX Basic to Offline Sampling Accessories With DC Sample collector
144	0104A00389	EDT Jar Lid for 250ml Vessels
145	0104A00385	EDT Jar Lock for Glass Vessel Suitable for Aluminum Top Plate (Set of 16)
146	0104A00101	EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big (Set of 6)
147	0104A00681	Tygon Tube Replacement And Performance Enhancement Kit For EDT-08LX With D Sample
148	0104A00401	EDT Laser Marked Glass Vessel 1000 ML - Amber With Sigma Delta Certificate
149	0104A00005	EDT Glass Vessel 1000ml - Amber
150	0104XA00004	EDT Glass Vessel 1000ml - Clear
151	0104A00023	EDT Laser Marked Glass Vessel 1000 ML - Clear
152	0104A00004	EDT Glass Vessel 1000ml - Clear
153	0104A00629	EDT Laser Marked Basket 40 Mesh
154	0104A00567	EDT Basket 40 Mesh
155	0104A00300	EDT Basket 40 Mesh for 150ml Vessel (Set of 6)
156	0104A00443	EDT -8 Mesh Felodipine Basket with Spacer (Set of 6)
157	0104A00077	EDT Laser Marked Basket 80 Mesh (Set of 6)
158	0104A00299	EDT Basket Holder with 'O' Ring (Set of 6) for 150ml
159	0104A00365	EDT Digital Vernier Caliper for 1000 ML Vessel
160	0104A00212	EDT-08L/08Lx 100 ML Low Volume Conversion Kit
161	0122A00019	EDT Sampling Manifold Cannula S.S for 2Ltr (Set of 6)
162	0104A00118	EDT Validation Kit - Analog
163	0104A00069	EDT Basket Holder with 'O' Ring (Set of 6)
164	0104A00451	EDT Digital Dial Gauge
165	0104A00515	EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 8)
166	0104A00709	Offline Sampling Accessories for EDT-8Lx Dissolution Tester with Dx Sample Collector
167	0104A00064	EDT Felodipine Basket with Spacer (Set of 6)
168	0104A00076	EDT Laser Marked Basket 60 Mesh (Set of 6)
169	0104A00078	EDT Laser Marked Basket 100 Mesh (Set of 6)
170	0104A00283	EDT-08L/08Lx – Adjustable Legs (Set of 4)
171	0104A00231	EDT Sampling Cannula Fixed Height 500 ML (Set of 6) Length 108 mm
172	0104A00042	EDT Low Evaporation Jar LID (Set of 6)
173	0104A00180	EDT Tube Connecting Coupling With Nozzle (Set of 2)(Replacement Of QRC)
174	0104A00307	EDT-08Lx Side Circulation Water Bath - Clear with

175	0104A00766	EDT Transdermal Rod For Cylinder - Height 107.9 mm (Set of 8)
176	0104A00181	EDT Tube Connecting Coupling (Set of 2)(Replacement Of QRC)
177	0104A00404	EDT 33mm Watchglass Disc Assembly With 16
178	0104A00405	EDT 55mm Watchglass Disc Assembly With 16
179	0104A00404	EDT 33mm Watchglass Disc Assembly With 16
180	0104A00406	EDT 90mm Watchglass Disc Assembly With 16
181	0105A00133	EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 28)
182	0104A00579	EDT-08L/08Lx 1 Ltr. Conversion Kit for 250 ML
183	0104A00577	EDT-08L/08Lx 1 Ltr. Conversion Kit for 150 ML
184	0104A00571	EDT Cooling Instrument (Temp. From 10° to 40° C) -
185	0104A00566	EDT Basket Holder with 'O' Ring
186	0104A00534	EDT-08L/08Lx -Motorized Sampling Manifold without Accessories
187	0104A00525	EDT 33mm Watch glass Disc.Assembly with 120 Mesh Screen
188	0104A00524	EDT 90mm Watch glass Disc.Assembly with 120 Mesh Screen
189	0104A00523	EDT 55mm Watch glass Disc.Assembly with 120 Mesh Screen
190	0104A00522	EDT Plastic 100ml Beakers For cleaning(Set of 8)
191	0104A00521	EDT-8LX Cleaning Tray
192	0104A00520	EDT-8/14LX Tablet Locator (set of 8)
193	0104A00519	EDT-8/14LX Sampling port seal (set of 8)
194	0104A00506	EDT-V 1234 Merlon Vessels set of 6 OLD type
195	0104A00476	EDT Customized Cage Type Sinkers Size 16mm ID x 13 mm HT With Spacer (Set of 6)
196	0104A00475	EDT Customized Cage Type Sinkers Size 12mm ID X 11mm Length (Set of 6)
197	0104A00473	EDT Replenishment Cannula (Set of 6)
198	0104A00472	EDT Laser Marked Kone Vessel 1000 ml - Clear With Sigma Delta Certificate
199	0104A00471	EDT Kone Vessel 1000 ml – Clear With Sigma Delta Certificate
200	0104A00465	EDT Customized Cage Type Sinkers Size 10mm ID Length x 10mm HT With Lid (Set of 6)
201	0104A00464	EDT Tablet Dispensing Derlin 'O' Ring for Dispensor plate (Set of 6)
202	0104A00462	EDT Manifold Rubber 'O' Ring for Dispensor plate
203	0104A00461	EDT Customized Cage Type Sinkers Size 14mm ID x 13 mm HT With Spacer (Set of 6)
204	0104A00460	EDT Cage Type Sinkers - Small With Spacer (Set of 6)

205	0104A00459	EDT Cage Type Sinkers - Big With Spacer (Set of 6)
206	0104A00458	EDT Customized Granules Dispensing Adapter For Dispenser Plate (Set of 12)
207	0104A00457	EDT Spacer For Cage Type Sinkers Size 22ID X 22mm HT with lid (Set of 6)
208	0104A00456	EDT Spacer For Small Cage Type Sinkers (Set of 6)
209	0104A00455	EDT Spacer For Big Cage Type Sinkers (Set of 6)
210	0104A00453	EDT-08L/08Lx – Auto Dispenser Plate
211	0104A00447	EDT 40 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)
212	0104A00446	EDT Customized Cage Type Sinkers – Medium With Lid (Set of 6)
213	0104A00444	EDT Teflon Coated Basket 10 Mesh (Set of 6)
214	0104A00441	EDT Basket LID 20 Mesh (Set Of 6)
215	0104A00437	EDT Bath Nozzle Rubber Gasket (Set of 2)
216	0104A00434	EDT-08L/08Lx Bath Guard (Set of 2)
217	0104A00433	EDT Silicon Tube Set For ETC-15
218	0104A00429	EDT-08Lx Conversion Kit for Enhancer Cell to 1Ltr
219	0104A00428	EDT Basket Rod for 250ml Vessel
220	0104A00425	EDT Basket Sinkers 40 Mesh (9.00mm) with Lid (Set of 6)
221	0104A00424	EDT Disposable Filter 3.00mm (Pkts of 100)
222	0104A00423	EDT-08L/08Lx Support Strip for Bath (Set of 2)
223	0104A00422	EDT-08Lx Top Cover
224	0104A00415	EDT Big Retaining Ring for Enhancer Cell Kit – 4cm (Set of 8)
225	0104A00414	EDT Top Cover with Base Plate for ETC-11LX
226	0104A00413	EDT Sampling Cannula Fixed Height 250 ML (Set of 6) Length 71 mm
227	0104A00412	EDT Luer Lock PMLP-21 (Male) (Set of 8)
228	0104A00411	EDT Silicon Tube for Drainage
229	0104A00409	EDT 120 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)
230	0104A00408	EDT 120 Mesh Screen For 55mm Watchglass Disc Assembly (Set of 6)
231	0104A00407	EDT 120 Mesh Screen For 33mm Watchglass Disc Assembly (Set of 6)
232	0104A00402	EDT Luer Lock FSLLR-1 (Female) (Set of 8)
233	0104A00397	EDT-08L/08Lx – Water Circulation Silicon Tube Set
234	0104A00393	EDT-08L/08Lx PTFE Carrier Tube From UV to Dissolution (OD -1.6mm / ID-0.8mm) - (Set Of 8)
235	0104A00392	EDT-08L/08Lx PTFE Carrier Tube From Agilent Pump To UV (OD-1.6mm / ID-0.8mm) - (Set Of 8)
236	0104A00391	EDT-08L/08Lx PTFE Carrier Tube From Disso to Agilent Pump (OD-1.6mm / ID-0.8mm) -(Set of 8)

237	0104A00390	EDT Washer for Stationary Basket (Set of 2)
238	0104A00388	EDT Analog Pressure Gauge For KBR Press
239	0104A00388	EDT Analog Pressure Gauge For KBR Press
240	0104A00386	EDT-08L/08Lx Aluminum Top Plate for Enhancer
241	0104A00384	EDT Luer Lock FTLL-6 Female (Set of 8)
242	0104A00379	EDT Glass Vessel Trolley
243	0104A00378	EDT Customized Transparent Moulded Jar Lid - Plastic (Set of 6)
244	0104A00377	EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 6)
245	0104A00376	EDT Sampling Cannula Fixed Height 500 ML (Set of 6) Length 127 mm
246	0104A00375	EDT-08L/08Lx – Tablet Dispenser Plate with Provision for Felodipine Basket
247	0104A00374	EDT Customized Helix Coil Sinkers 9.5mm ID X 20mm Height (Set of 6)
248	0104A00372	EDT Hardened Polished Surface Plate (for Intrinsic)
249	0104A00371	EDT Bolts for Fixing the Die and Surface Plate (for Intrinsic)(Set of 3)
250	0104A00370	EDT Stand For 1000ml Vessel (Composite)
251	0104A00369	EDT Deconex CIP Power-X (14 Kg)
252	0104A00368	EDT Online Connector (Set of 8)
253	0104A00364	EDT Glass Vessel 'D' Type 1000 ML - Amber With
254	0104A00351	EDT 15ml Replenishment Bottle With Rubber Cap
255	0104A00350	EDT Luer Lock MLRL210-6 Male (Set of 8)
256	0104A00348	EDT Customized 10 Mesh Basket Sinkers with Lid (12.7mm ID – 24.5mm Height) (Set of 6)
257	0104A00347	EDT Customized 10 Mesh Basket Sinkers with Lid (9.1mm ID – 18.8mm Height) (Set of 6)
258	0104A00345	EDT Prednisone Tablet (Pkt of 30)
259	0104A00344	EDT Prednisone Standard USP-1559006 Reference
260	0104A00343	EDT Disc-Filter-Millex-HN, 33mm, 0.45um, Nylon Membrane, Non-Sterile (Pkt of 250)
261	0104A00342	EDT Disc Filter-Millex -HN, 33mm, 0.45um PES Membrane, Non-Sterile (Pkt of 250)
262	0104A00339	EDT-08L/08Lx PTFE Carrier Tube From Disso to Watson Marlow Pump (Set of 8)
263	0104A00337	EDT-08L/08Lx PTFE Carrier Tube From Disso to Ismatec Pump (Set of 8)
264	0104A00336	EDT Probe for Vibration Meter
265	0104A00335	EDT Sampling Cannula Fixed Height For Enhancer Cell Kit -200ml (Set of 6)
266	0104A00333	EDT Cylindrical Flat Bottom Glass Vessel 200ml – Clear for Enhancer Cell (Set of 6)

267	0104A00327	EDT Enhancer Cell for Enhancer Cell Kit – 2 CM (Set
268	0104A00325	EDT-08L/08Lx -Motorised Sampling Manifold (For Peristaltic Offline)
269	0104A00324	EDT Sampling Manifold Cannula S.S. (SET Of 6)
270	0104A00323	EDT-08L V-34 Syringe Manifold (for Offline Syringe)
271	0104A00322	EDT-08L/08LX PTFE Carrier Tube For Syringe Manifold (From Sampling Probe To Syringe
272	0104A00321	EDT Luer Lock BSMLRL-6 Male (Set of 8)
273	0104A00320	EDT Enhancer Cell for Enhancer Cell Kit - 4 CM (Set
274	0104A00316	EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6)
275	0104A00314	EDT Sinkers 20 Mesh with Lid (Set of 6)
276	0104A00316	EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6)
277	0104A00312	EDT-08L/08Lx S.S. Strip for Composite Top Plate
278	0104A00288	EDT Basket Sinkers 40 Mesh (11.7mm Height) with Lid (Set of 6)
279	0104A00284	EDT Sinkers 10 Mesh with Lid (Set of 6)
280	0104A00279	EDT Basket 10 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate
281	0104A00277	EDT-08L/08Lx -Motorised Sampling Manifold (For Syringe Offline)
282	0104A00276	EDT Deconex-CIP Fresh
283	0104A00275	EDT Deconex-CIP Alpha- X
284	0104A00274	EDT U Glass Tube
285	0104A00273	EDT Teflon Coated Basket Shaft With O Ring &
286	0104A00271	EDT Rubber Ring For Low Evaporation Jar Lid
287	0104A00270	EDT-08L/EDT-08Lx Composite Top Plate for 100ml Low Volume Conversion Kit
288	0104A00269	EDT Spirit Level
289	0104A00268	EDT Top Cover with Base Plate for ETC-11
290	0104A00266	EDT Basket 100 Mesh "V" Type
291	0104A00264	EDT Basket 40 Mesh "H" Type
292	0104A00263	EDT Basket 100 Mesh "H" Type
293	0104A00261	EDT S.S Knob for Bath Rod
294	0104A00259	EDT Merlon Vessel 1000ml Amber (N)
295	0104A00254	EDT SS Filters with Holders – 40 Micron (Set of 6)
296	0104A00255	EDT SS Filters with Holders – 80 Micron (Set of 6)
297	0104A00253	EDT Rubber Holder For SS Filter
298	0104A00248	EDT Peak Glass Lense to Measure Basket Mesh
299	0104A00246	EDT Polished Punch
300	0104A00241	EDT Deconex 24 Liq (7 Kg)
301	0104A00240	EDT Deconex 25 Organ Acid (5 Kg)
302	0104A00237	EDT Vernier Caliper with Stopper for 1000ml

303	0104A00235	EDT Low Evaporation Jar Lid for Sampling Cannulas Fixed Height Suitable for Glass Vessel (Set of 6)
304	0104A00232	EDT Sampling Cannula Fixed Height 100 ML (Set of 6) Length 73.5 mm
305	0104A00229	EDT Dispenser Jar Ring Black (Set of 8)
306	0104A00228	EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate
307	0104A00228	EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate
308	0104A00237	EDT Vernier Caliper with Stopper for 1000ml
309	0104A00229	EDT Dispenser Jar Ring Black (Set of 8)
310	0104A00228	EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for Composite Top Plate
311	0104A00227	EDT Merlon Vessel 1000 ML - Amber with Jar Lock Suitable for S.S. Top Plate
312	0104A00226	EDT Merlon 1000 ML - Clear with Jar Centering Ring Suitable for Composite Top Plate
313	0104A00224	EDT Spindle Latch with Spring (Set of 8)
314	0104A00220	EDT Japanese Sinkers with 40 Mesh Cover (Set of 6)
315	0104A00219	EDT Sinkers 40 Mesh with Lid (Set of 6)
316	0104A00218	EDT Customized Japanese Sinkers 15 mm ID X 40 mm Height (Set of 6)
317	0104A00206	EDT Holder for Basket 40 Mesh 24.5 mm (Set of 6)
318	0104A00203	EDT Deconex 12 Basic (500ml)
319	0104A00185	EDT Stop Watch
320	0104A00184	EDT Transdermal Cylinder - Height 50.8 mm (Set of 6)
321	0104A00177	EDT 'D' Type Basket 40 Mesh (Set of 6)
322	0104A00174	EDT 'D' Type Basket 10 Mesh (Set of 6)
323	0104A00167	EDT Luer Lock FTLL210-6 Female (Set of 8)
324	0104A00166	EDT Luer Lock MTLL220-6 Male (Set of 8)
325	0104A00160	EDT Gauge Block 25 mm
326	0104A00154	EDT Basket Validation Gauge with Certificate
327	0104A00152	EDT Moulded Tray For Basket, Basket Holder & Rod
328	0104A00145	EDT Plastic Quick Release Coupling with Nozzle (Set of 8)
329	0104A00141	EDT Vessel Centering Ring for Merlon Vessel Suitable for Composite Top Plate (Set of 8)
330	0104A00140	EDT Vessel Centering Ring for Glass Vessel Suitable for Composite Top Plate (Set of 8)
331	0104A00136	EDT-08L/08Lx Composite Top Plate for Glass
332	0104A00135	EDT Vessel Centering Plate for Merlon Vessel Suitable for S.S. Top Plate
333	0104A00134	EDT TNT Washer for Merlon Vessel Suitable for S.S. & Top Plate (Set of 16)

334	0104A00132	EDT Jar Lock for Glass Vessel Suitable for S.S. Top Plate (Set of 16)
335	0104A00131	EDT-08L/08Lx S.S. Top Plate for Merlon Vessel
336	0104A00130	EDT-08L/08Lx S.S. Top Plate for Glass Vessel
337	0104A00125	EDT Digital Verticality Meter
338	0104A00123	EDT Thermometer - Mercury
339	0104A00122	EDT Digital Tachometer
340	0104A00121	EDT Depth Gauge for 1000 MI Vessel
341	0104A00120	EDT Centering Gauge
342	0104A00119	EDT Centering Device for 1000 MI Vessel
343	0104A00116	EDT Helix Coil Sinkers - CAPWHT-2S (Set of 6)
344	0104A00119	EDT Centering Device for 1000 MI Vessel
345	0104A00115	EDT Helix Coil Sinkers - CAPWHT-02 (Set of 6)
346	0104A00114	EDT Metformin Sinkers - Small (Set of 6)
347	0104A00113	EDT Metformin Sinkers - Big (Set of 6)
348	0104A00110	EDT Sinkers 'U' Type (Set of 6)
349	0104A00107	EDT Alfuzosin Sinkers Sinkers 37.5 mm ID x 20 mm OD (Set of 6)
350	0104A00105	EDT 3 Prong Sinkers - SS (Set of 6)
351	0104A00104	EDT 3 Prong Sinkers - Plastic (Set of 6)
352	0104A00102	EDT Spring Type Sinkers 9 mm ID x 30 mm OD - Small (Set of 6)
353	0104A00100	EDT Cage Type Sinkers - Small (Set of 6)
354	0104A00081	EDT Laser Marked Felodipine Basket (Set of 6)
355	0104A00080	EDT Laser Marked Stationary Basket (Set of 6)
356	0104A00075	EDT Laser Marked Basket 40 Mesh (Set of 6)
357	0104A00072	EDT Laser Marked Basket 10 Mesh (Set of 6)
358	0104A00067	EDT Basket 'O' Ring (Set of 6)
359	0104A00066	EDT Basket Rod
360	0104A00065	EDT Moulded Tray For Basket
361	0104A00062	EDT Stationary Basket (Set of 6)
362	0104A00060	EDT Basket 100 Mesh (Set of 6)
363	0104A00059	EDT Basket 80 Mesh (Set of 6)
364	0104A00058	EDT Basket 60 Mesh (Set of 6)
365	0104A00055	EDT Basket 20 Mesh (Set of 6)
366	0104A00055	EDT Basket 20 Mesh (Set of 6)
367	0104A00052	EDT Inline Filters - 0.45 Micron (Set of 6)
368	0104A00051	EDT SS Filters with Holders (Set of 6)
369	0104A00045	EDT Sampling Cannulas Fixed Height (Set of 6)
370	0104A00044	EDT Jar LID For Staggered Start - White (Set of 6)
371	0104A00041	EDT Moulded Jar LID - Plastic (Set of 6)
372	0104A00040	EDT Rubber Stand For Glass Vessel
373	0104A00037	EDT Laser Marked Glass Vessel With Ring 1000ml -
374	0104A00035	EDT Laser Marked Glass Vessel 'V' Type 1000ml -

375	0104A00033	EDT Laser Marked Glass Vessel 'D' Type 1000 MI -
376	0104A00032	EDT Laser Marked Flat Bottom 1000ml - Amber
377	0104A00032	EDT Laser Marked Flat Bottom 1000ml - Amber
378	0104A00022	EDT Laser Marked Glass Vessel 250 MI - Amber
379	0104A00022	EDT Laser Marked Glass Vessel 250 MI - Amber
380	0104A00018	EDT Glass Vessel With Ring 1000ml - Clear
381	0104A00017	EDT Glass Vessel V Type 1000 ml – Amber
382	0104A00012	EDT Flat Bottom 1000 MI - Clear
383	0104A00127	EDT Vibration Meter
384	0104A00484	EDT Spring Assembly for Tablet Dispenser Plate
385	0104A00442	EDT Deconex 12 Basic (10 Bottles of 1000 ml)
386	0104A00623	EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big
387	0104A00063	EDT Felodipine Basket (Set of 6)
388	0104A00584	EDT 40 Mesh Stationary Basket (Set of 6)
389	0104A00057	EDT Basket 40 Mesh (Set of 6)
390	0104A00410	EDT Clearview (230V)
391	0104A00267	EDT Basket 20 Mesh "V" Type
392	0104A00175	EDT 'D' Type Basket 20 Mesh (Set of 6)
393	0104A00265	EDT Basket 20 Mesh "H" Type
394	0104A00613	EDT-08Lx PTFE Bath Top Plate
395	0104A00581	EDT 16 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)

Sr.No.	Part code	Breakdown accessories
1	0104BV4154	EDT-Lx V5 Heater Card for ETC-15 Without Trim Pot
2	0104BV4153	EDT-Lx V5 Compatch Card Without Trim Pot
3	0104BV4151	EDT-08Lx V-5 Individual Vessel Temp. Sensor PCB with Cable & Connector
4	0104BV4147	EDT 08L V-1234 Diving Manifold Bearing's Housing (Set of 2)
5	0104BV4145	EDT-Lx V-4 Driver Card (110V)
6	0104BV4141	EDT-Lx V-45 Mains Supply Filter for ETC-15 (230V)
7	0104BV4140	EDT - Lx V-4 SMPS - 12V
8	0104BV4137	EDT-LX V-4 Mains Cord
9	0104BV4135	EDT-Lx V-4 Inlet Connector for Side Circulation Waterbath
10	0104BV4115	EDT-08Lx V-4 – 10 Pin FRC Cable (Processor Card to Stirrer Opto Card)
11	0104BV4114	EDT-Lx V-4 Temperature Controller ETC-15
12	0104BV4091	EDT Lx V-4 Communication Cable (Syringe Pump to Sample Collector)
13	0104BV4090	EDT Lx V4 Communication Cable (Disso to Syringe Pump)
14	0104BV4087	EDT-Lx V-4 Serial Printer Cable
15	0104BV4086	EDT-Lx Opto PCB For Encoder
16	0104BV4071	EDT-Lx Encoder With Flywheel
17	0104BV4054	EDT-08Lx V-4 – 20 Pin FRC Cable (Processor Card to Keypad)
18	0104BV4048	EDT-Lx V-4 Lift Opto Bottom
19	0104BV4030	EDT-Lx V4 Main Cable 15 AMP
20	0104BV4016	EDT-Lx V-4 LCD Display
21	0104BV4005	EDT-Lx V-4 RPM Stirrer Motor
22	0104B04138	EDT V-1234 Opto Disc for Lift Motor
23	0104B00168	EDT-08L/14L V-234- 3V DC Battery for Processor Card
24	0104B00099	EDT V-1234 Link Cord
25	0104B00096	EDT V-1234 Bath Coupling Mounting Plate (Set of 2)
26	0104B00087	EDT V-1234 Inlet & Outlet Connector for
27	0104B00065	EDT-08L V-1234 Bath Top Plate Knob with Washer (Set of 4)
28	0104B00052	EDT-08L V-1234 Small Rubber Leg for Water Bath 22mm (Set of 4)
29	0104B00050	EDT V-1234 Drain Rubber for Water Bath
30	0104B00048	EDT Dispensor Plate Holder (Set of 2)

31	0104B00047	EDT V-1234 Inline Fuse Holders
32	0104B00046	EDT V1234 Mains Switch (Yellow)
33	0104B00045	EDT V-1234 Rocker Switch Mains Red
34	0104B00041	EDT V-1234 Lift Motor 230V-2A
35	0104B00037	EDT V-1234 Lift Motor 110V-3A
36	0104B00032	EDT V-1234 Links Cable 5 Amp
37	0104B00029	EDT V-1234 Mains Fuse Holder
38	0104B00028	EDT V-1234 Thermal Cut-Off
39	0104B00026	EDT V-1234 Timing Pulley Set
40	0104B00025	EDT V-1234 Drainage Valve
41	0104B00009	EDT V-1234 Manifold Opto
42	0104B00008	EDT-08L V-1234 Manifold Bearing (Set of 2)
43	0104B00003	EDT-08L/14L V-34 Manifold Bush Highspeed
44	0104BV4173	EDT V -1234 LED with Cable (Set of 2)
45	0104BV4188	EDT-Inspire BLDC Drive 8 Pole
46	0104BV4177	EDT V-1234 Cut spacer 0.2mm
47	0104BV4176	EDT V-1234 Cut Spacer 0.5mm
48	0104B00012	EDT V-1234 Individual Vessel Temperature Sensor (Set of 8)
49	0104B00014	EDT V-1234 Single Pulley Spindle
50	0104B00015	EDT V-1234 Double Pulley Spindle
51	0104BV4027	EDT-Lx V-4 Fuse Set
52	0104BV4163	EDT-8/14Lx V-4-Inspire Driver Card
53	0104BV4098	EDT-Lx V-4 Heater- 110V (For ETC-15)
54	0104B00101	EDT-Lx V-4 Heater- 110V (For ETC-15) Old Straight Nozzle
55	0104BV4128	EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 230V
56	0104BV4155	EDT V-5 External Probe Without Trim Pot
57	0104BV4183	EDT-LX V4 Inspire Processor Card
58	0104BV4099	EDT-08L Bath Support Frame
59	0104BV4080	EDT-Lx V4 Compatch Card
60	0104BV4085	EDT-Lx V4 Gear for Lift Motor
61	0104BV4158	EDT-Lx V-45 Advance Filter for ETC-15 (230V)
62	0171B00001	EDT-08/14Tx LCD Module (Touch Screen)
63	0104B00039	EDTV-1234 Heater 230V-2A
64	0104BV4185	EDT- V-4-Inspire Stirrer opto card
65	0104B00022	EDT-08L V-1234 Timing Belt Set
66	0104BV4024	EDT-Lx V-4 Main Socket with Power Supply Switch & Filter
67	0104B00134	EDT-08L V-1234- 16 Pin FRC Cable (Processor Card to Driver Card)
68	0104BV4156	EDT-Lx V-4 Mains Cord For ETC-11Lx

69	0104BV4059	EDT-Lx V-4 Lift Disc Opto Card
70	0104BV4223	EDT-08LX Offline With Fraction Collector Platinum Kit Electronic Spares
71	0104B00083	EDT V-1234 MF Connector for Manifold
72	0104B00018	EDT V-1234 External Probe
73	0104B00013	EDTV-1234 Bath Temperature Sensor
74	0171B00009	EDT-08/14Tx Lift motor drivr card
75	0104BV5159	EDT-08Lx V-5 Individual Vessel Temperature Sensor PCB with Cable & Connector with out trim pot
76	0104BV4001	EDT-Lx V4 Stirrer Opto Card
77	0104B00049	EDT V-1234 Individual Vessel Temperature Sensor
78	0104BV4081	EDT-Lx V-4 Temperature Controller ETC-11
79	0104BV4084	EDT-Lx V-4 Lead Screw for Lift
80	0105BV4164	EDT- Wear & Tear Kit (For Dx)
81	0104BV4178	EDT V-1234 Height adjustment spacer 0.5mm
82	0104BV4074	EDT-Lx V-4 Processor Card
83	0168B00002	EDT-08Tx Processor Card
84	0104BV4159	EDT-08Lx V-45 Wear & Tear Kit (For Dx)
85	0104BV4149	EDT-08Lx V-4 Individual Vessel Temp. Sensor PCB (Black) with Cable & Connector
86	0104BV4221	EDT-08LX Offline With D Sample Collector Platinum Kit Electronic Spares
87	0104B00058	EDT-08L V-1234 Water Circulation Tygon Tube Set
88	0104BV4127	EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 110V
89	0104BV4092	EDT-Lx V-4 Circulation Pump – 24V (For ETC-15)
90	0104BV4095	EDT-Lx V-4 Heater Card for ETC-15
91	0104BV4175	EDT V-1234 Height adjustment spacer 0.2mm
92	0104B00143	EDT V-1234 LED For ETC-15
93	0104BV4095	EDT-Lx V-4 Heater Card for ETC-15
94	0104BV4239	Inspire 8 Base Frame Assembly
95	0104BV4237	Inspire Processor card for 10 to 55 degC
96	0104BV4201	Inspire Lift Park Position Opto Sensor
97	0104BV4238	Inspire Program Patch
98	0104BV4225	Inspire-08 Platinum Kit Mechanical Spares.
99	0104BV4226	Inspire-08 Platinum Kit Electronic Spares.