

EHT (Gen2) Tablet Tester



EHT-5P (ET) / EHT-3P (ET) (CFR License)

Features:

- Biometric Login
- Robust Design
- Report print serial / LAN port
- 21 CFR part 11 Compliance
- Automatic calculation
- User friendly interface
- Audit Trail
- Based on measurement parameter requirement model can be selected (Refer model table)
- Up to 200 samples per test
- Jaw with groove
- Adjustable break between the samples
- Constant speed measuring principles available as standard
- Programmable calibration and verification modes
- Auto/Manual run mode with programmable delay time
- User selectable templates for printing data
- Basic security
- Plausibility Checking
- Maximum load cell force 0-800N
- Multi Point Verification
- Calibration and verification reminders
- Body made of S.S to comply with CGMP regulations
- Perfect companion for IPQC cubicals
- Compatible with Isolators
- Protective cover to protect operator from Inhaling tablet dust from tablets that 'explode'
- IP 41 equivalent
- Machine to machine data transfer
- DHCP/Static LAN option
- Network time sync option
- Advance security
- Logbook
- Report template
- Backup/restore
- Weighing parameter through external balance

Specifications:

Model	EHT-5P (ET) / EHT-3P (ET)
Tablet Shape	Support all shapes
CFR License	Users : 250 Roles : 10 Audit Trails : 200000 Report template : 10 Reports : 50000 Products : 1000 Methods : 1000 Logbook : 200000
Display	7" full color intuitive touch screen display
Hardness Unit	Newton, Kilo pond, Strong Cobb, Custom factor
Hardness Range	0 to 800 N
Accuracy	better than ± 1 N
Hardness calibration	Static with reference weights
Illimination	D'Light
Dia / Width Tolerance	50 microns
Thickness	up to 15 mm
Thickness Tolerance	50 microns
Weighing balance	Mettler, Sartorius, Contech, Ohaus
Weighing unit	mg / gm
Sample per test	up to 200 configurable
Print Templates	User defined template for printing data
Dimensions (L x W x H)	640 mm x 540 mm x 660 mm
Electrical power	220/230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/10 Amp
Weight	11 Kg

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

EHT (Gen2) Tablet Tester

EHT (Gen2) Models							
				*E - External Thickness		*T - Thermal Printer	
EHT Models	Hardness	Weight	Thickness	Length / Diameter	Width	External Thickness	Inbuilt Printer
2P	✓	✓					
3P	✓	✓	✓				
3P (E)	✓	✓	✓			✓	
5P	✓	✓	✓	✓	✓		
5P (E)	✓	✓	✓	✓	✓	✓	
2P (T)	✓	✓					✓
3P (T)	✓	✓	✓				✓
3P (ET)	✓	✓	✓			✓	✓
5P (T)	✓	✓	✓	✓	✓		✓
5P (ET)	✓	✓	✓	✓	✓	✓	✓

Ordering information:

Products		
Product code	Description	Quantity
0942000	Tablet Hardness Tester Gen 2 Model - EHT 5P (ET) (230V)	1
0932000	Tablet Hardness Tester Gen 2 Model - EHT 5P (ET) (110V)	1
0942700	Tablet Hardness Tester Gen 2 Model - EHT 3P (ET) (230V)	1
0932700	Tablet Hardness Tester Gen 2 Model - EHT 3P (ET) (110V)	1

CFR License		
Product code	Description	Quantity
09420SL05	EHT-5P (ET) Gen 2 CFR License 7 year	1
09420SL03	EHT-5P (ET) Gen 2 CFR License 1 year	1
09427SL05	EHT-3P (ET) Gen 2 CFR License 7 year	1
09427SL03	EHT-3P (ET) Gen 2 CFR License 1 year	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Ordering information:

Standard Accessories		
Product code	Description	Quantity
0906A00023	EHT Gen.2 Hood Assembly	1
0906A00024	EHT Gen.2 Waste Collection Tray	1
0906A00025	EHT Gen.2 Calibration Bracket, Pin & Plate	1
0906A00027	EHT Gen.2 Base Platform 41mm Big Groove	1
0901A00002	Brush	1

Optional Accessories		
Product code	Description	Quantity
0906A00026	EHT Gen.2 Stationary Jaw 66mm	1
0906A00028	EHT Gen.2 Base Platform 41mm Plain	1
0906A00029	EHT Gen.2 Base Platform 41mm Small Groove	1
0906A00030	EHT Gen.2 Base Platform 65mm Big Groove	1
0906A00031	EHT Gen.2 Base Platform 65mm Plain	1
0906A00032	EHT Gen.2 Base Platform 65mm Small Groove	1
0906A00033	25 mm Slip Gauge	1
0905A00013	EBT/EHT Validation Kit (10mm Gauge Block + 5 Kg Calibrated Weight)	1
0104A00448	Biometrics (Gen.2)	1

IQ/OQ Document		
Product code	Description	Quantity
0420S15	IQ/OQ Document For Electrolab Tablet Hardness Tester Gen 2 Model - EHT 5P (ET)	1
0927S15	IQ/OQ Document For Electrolab Tablet Hardness Tester Gen 2 Model - EHT 3P (ET)	1

Upgradation		
Product code	Description	Quantity
	Upgradation For Electrolab Tablet Hardness Tester Gen 2 Model - EHT 5P (ET)	1

One time development for Report Customization		
Product code	Description	Quantity
	One time development for Report Customization For Electrolab Tablet Hardness Tester Gen 2 Model - EHT 5P (ET)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com