

# EHT-01MBP Tablet Hardness Tester



## EHT-01MBP

### Features:

- Compliance with the USP <1217> Pharmacopoeia
- High accuracy, high resolution
- Date and Time are included on all reports
- Strap for easy carrying from one place to another  
Instrument identification and calibration sticker is placed on the front panel for easy reference
- Three units N, Kgf, lbf for measurement provided and can convert between the three units automatically
- Report upto 50 reading print with instrument details, test details and Statistics (min, max and average)
- Calibration and Verification available with reports
- Thermal Printer

### Specifications:

Model	EHT-01MBP
Hardness Range	Up to 500 Newton
Measuring Units	N, Kgf, lbf
Display	OLED 1.3"
Test Mode	Peak, Auto
Platform Length	30 mm
Accuracy	± 1 N
Calibration	Yes with Report
Validation	Yes with Report
Tablet Size	2 mm to 28 mm
Number of Samples	50 (Per test)
Battery	3.7V Battery 5000mAh
Power Supply	5VDC - C type USB adapter
Dimensions (L x W x H)	295 mm x 90 mm x 100 mm
Weight (Approx)	1.7 Kg

### Ordering information:

Product		
Product code	Description	Quantity
0943400	Electrolab Tablet Hardness Tester Model EHT-01MBP	1
Standard Accessories		
Product code	Description	Quantity
0905A00023	EHT-01M Calibration Platform	1
0905A00028	EHT-01M Knob	1
0905A00027	EHT-01M Sliding Cover	1
0908B00016	EHT-01M Charging Adapter with USB C Type cable, 5V/2.4A	1
0905A00024	EHT-01M Waste tray Acrylic	1
0901A00002	Brush	1
0905A00026	EHT-01MB/MBS 3.7V 5200mAh Lilon Battery Pack	1
0905A00025	EHT-01MP 40 Column Printer with cable	1
Optional Accessories		
Product code	Description	Quantity
0901A00012	S.S. Weight 5kg (Calibrated) With Certificate	1
IQOQ Document		
Product code	Description	Quantity
0934515	IQ/OQ Document For Electrolab Tablet Hardness Tester Model EHT-01MBP	1

\*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupa Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com