

Dissolution Tester

VERITAS
Basic




ELECTROLAB
Your Quality, Our Assurance

www.electrolabgroup.com

Dissolution Testers



8 Station Bathless Dissolution Tester VERITAS-08



Fingerprint Authentication



Continuous Temperature Monitoring



Time Saving



Maintenance Free



Identifies Position



Touch Screen



Ethernet Connectivity



Infinity Test



Report

14 Station Bathless Dissolution Tester VERITAS-14





Key Features

- 7" Tft Touch Screen Display
- Lan Connectivity Allows Central Printing
- Multi Level Security With Password For 250 Users With User Programmable Role
- Bathless Heating (saves Water, Time, Electricity)
- Complies With Current Usp , Ip, Ph. Eur. And Jp Specifications
- Unique Time Action Tm Function Which Allows Media Changeover And Infinity Test Run, Ideal For Sustained And Controlled Release Products.
- Programmable Different Aliquot Volume For Each Sample Interval Allows Variable Dilution
- Automatic Stirrer Height Positioning For Usp 1, 2, 5, 6 Intrinsic Test Methods
- Easy Snap-fit™ Shaft Locking Mechanism For Positive Engagement And Wobble Free Operation
- Bathless Dissolution Tester Is Maintenance-free
- Specially Designed Efficient Heater Jackets For Vessel, Heats The Media Quickly Attains Temperatures From 25°C To 37°C Within Approximately 12 Minutes
- V-stir™ (Patented) Ensures Continuous Wireless Temperature Monitoring Without Affecting The Hydrodynamics In The Vessel.
- V-stir™ (patented) Identifies Position And Apparatus Type And Apparatus Presence To Reduce Manual Errors
- Intuitive 7" Tft Touch-screen Graphic Display
- Visual Indication Along With Sound Alert Of Dissolution.
- Ethernet Connectivity Allows Central Printing Of Data.
- In-built Accessory Drawer For Storing The Baskets And Paddles
- Instrument Status From All Sides Of The Instrument:
 - Idle
 - Ready
 - Run
 - Error