

### Kloudface EKF-24 for EDI-2SA & ED-2SAPOx model



**EKF-24**

#### Features:

- 21 CFR Part 11 compliance
- Visual Indication using D'light
- Multi Level Security
- Electronic Instrument Log Book
- Biometric (Optional)
- Control the parent instrument and all operation through EKF only

#### Specifications:

Model	EKF-24
Display	Multi touch capacitive 7" screen
Firmware	Integrity operating System Integros 2.0
Firmware Type	Lite, Standard and Pro
Printer Interface	Serial / LAN
Printer	80 column
Backup	Automatic backup with koudbase server / LAN (Reports and Audit Trail events)
Memory	200000 Audit Trail
	50000 Reports
	1000 Products
	250 Users
	20 Roles
	10 Report Templates
Dimension (L x W x H)	207 mm x 231 mm x 83 mm
Weight	1 Kg (approx)

\*The information contained in this document is believed to be correct but **ELECTROLAB** accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

**Head Office:** 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

**Factory:** EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

### Kloudface EKF-24 for EDI-2SA & ED-2SAPOx model

## Ordering Information

Products		
Product code	Description	Quantity
1542600_1	Kloudface Backup Model EKF-24	1

  

IQOQ Document		
Product code	Description	Quantity
1526S15	IQOQ Document for Kloudface Backup Model EKF-24	1

  

Extended IQOQ Document		
Product code	Description	Quantity
E1526S15	Extended IQOQ Document for Kloudface Backup Model EKF-24	1

  

Software		
Product code	Description	Quantity
15426S1	Software EKF-24 IntegrOs2 Pro License - 1 Year	1

  

Breakdown Accessories		
Product code	Description	Quantity
0301B00055	EKF LCD Display (7 inch display)	1
0201B00096	EKF Power Adapter ( DC 12V / 2Amp)	1
0201B00097	EKF Power Supply Module	1
0150B00001	EKF-Gen.2 SMPS	1
0150B00002	EKF-Gen.2 Compatch Kit	1
0150B00003	EKF-Gen.2 Touch Screen	1
0150B00010	EKF-Gen.2 Y Cable With Connector	1
0150B00004	EKF-Gen.2 Interface Card Ver.1.2	1
0150B00005	EKF-Gen.2 Micro SD 32GB Card	1
B000150012	EKF-Gen.2 Base And Top Cover	1
0150B00007	EKF-Gen.2 Delight Ver.2	1
0150B00008	EKF-Gen.2 Biometric Sensor	1
0150B00009	EKF-Gen.2 Data Converter Module	1
0150B00011	EKF-Gen.2 Pen drive 32GB	1

\*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

**Head Office:** 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131    sales@electrolabgroup.com

**Factory:** EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122    www.electrolabgroup.com