

Leak Tester



ELT-201

Features:

- Meets USP and package specifications
- Ideal for checking leakage of packed strips, blisters, small sachets and ampoules
- Programmable vacuum level and hold time
- Slow venting for better results
- Printout of test report via Serial port & Ethernet
- Oil-less vacuum pump
- 5.7 inch graphical touch screen display
- Tolerance range setting to give pass / fail results
- Programmable 20 methods and storage of 20 product data

Specifications:

Model	ELT-201
Vacuum level	2 mmHg to 650 mmHg
Timer	Maximum can be set up to 00.01 minutes Minimum can be set up to 99.59 min
Vacuum setting level	Vacuum can stabilize to the set value ± 5 mmHg of Hg within 90 secs.
Power requirement	230V. 50 Hz
Power	300 VA

Ordering Information:

Products		
Product code	Description	Quantity
5040400	Leak Tester Model ELT-201 (230V)	1
5030400	Leak Tester Model ELT-201 (110V)	1

Optional accessories		
Part number	Description	Quantity
5001A00022	12 mm Round Plexiglass Plate for 250 mm Dia Dessicator	1
5001A00023	12 mm Round Plexiglass Plate for 300 mm Dia Dessicator	1
5001A00009	Vacuum desiccators of 300 mm dia	1
5001A00008	Vacuum desiccators of 250 mm dia	1
5001A00014	Vacuum Chamber: 300 mm dia x 450 mm Height	1
5001A00016	Vacuum Chamber: 400 mm dia x 290 mm Height	1
5001A00017	Vacuum Chamber: 400 mm dia x 450 mm Height	1
5001A00013	Vacuum chamber: 300 mm dia x 290 mm height	1
5001A00021	Vacuum Chamber: 300 mm Dia x 145 mm Height	1
5001A00018	Vacuum Chamber: 400 mm dia x 500 mm Height	1
5001A00015	Vacuum Chamber: 300 mm dia x 500 mm Height	1
0201A00020	80 Column serial printer	1
0104A00785	connexion+	1
1544006	Sentin Ai (2-Port)	1
1544004	Sentin Ai (2-Port) Machine Learning	1
1544003	Sentin Ai (2-Port) Machine Learning Advance	1
1544007	Sentin Ai (4-Port)	1
1544005	Sentin Ai (4-Port) Machine Learning	1
1544002	Sentin Ai (4-Port) Machine Learning Advance	1

Breakdown Parts		
Product code	Description	Quantity
5004B00001	ELT-201 Display Card (V-1.2)	1
5004B00002	ELT-201 LCD Module (Touch Screen)	1
5004B00003	Communication Cable 9 pin to 25 pin For Sartorius Weighing	1
5004B00004	ELT Solenoid Valve	1

IQ/OQ Document		
Product code	Description	Quantity
5004S15	IQ/OQ Document for Leak Tester Model ELT-201	1

Extended IQ/OQ Document		
Product code	Description	Quantity
E5004S15	Extended IQ/OQ Document for Leak Tester Model ELT-201	1

In-House Repairing		
Product code	Description	Quantity
5004S07	In-House Repairing for Leak Tester Model ELT-201	1

Single Validation Charges		
Product code	Description	Quantity
5004S08	Single Validation Charges for Leak Tester Model ELT-201	1

Inspection & Service Charges		
Product code	Description	Quantity
5004S06	Inspection & Service Charges for Leak Tester Model ELT-20	1

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

+91 - 22 - 4161 3122 www.electrolabgroup.com