

Tablet Dissolution Tester Trust€-08 with Sample Collector (Dc)



Trust€-08 with Sample Collector (Dc)

Features:

- 7" TFT touch screen display
- LAN connectivity allows central printing
- Multilevel security with password for 250 users with user programmable role
- Alphanumeric data entry of product and user details
- Complies with USP, IP, Ph. Eur. and JP specifications
- Unique TimeAction™ function which allows media changeover and infinity test run, ideal for sustained and controlled release products
- Programmable different aliquot volume for each sample Interval allows variable dilution
- Automatic stirrer height positioning for USP 1, 2, 5, 6, intrinsic test methods
- Easy Snap-Fit™ shaft locking mechanism for positive engagement and wobble free operation
- Molded clearwater bath
- Easy to empty and clean water bath with bottom drainage
- Precise individual vessel centering system
- Adjustable leveling of instrument
- Vibration free smooth motorized lift movement
- Isolated temperature controller and water circulating pump (vibration-free & precise) with wake up and sleep modes which saves time and energy
- 1000 programmable products
- 24 programmable sampling intervals
- 10000 test report storage
- Instrument accessories like paddles, baskets, vessels etc. are laser marked with serial numbers
- Designed to minimize routine validation
- Programmable report template
- Report for online validation of test, monograph, error log during the test
- Programmable sample volume from 1 ml to 20 ml
- Audio-visual indicators for system status
- Power failure detection and auto recovery facility
- D'light technology enabling audio visual indication of instrument status
- Finger print module
- 250 users
- Auto drop sync
- Auto dilution with stirring
- Intense rinse
- Rinse by purge
- Rinse by loopback
- Dynamic Dive to lower the sampling cannula to the regulatory complaint depth depending on the aliquots withdrawn

Specifications:

Model	Trust€-08 with Sample Collector (Dc)
No. of stations	08
Speed range	20 RPM to 300 RPM
Speed accuracy	± 0.5 RPM
Temperature range	20°C to 40°C (extended 10°C to 40°C and 20°C to 55°C optional)
Temperature accuracy	± 0.1°C
Display resolution	0.1°C
Product Storage	1000
Display	7" TFT touch screen display
Test Report Storage	10,000
Sample interval	1 min to 999 hrs, 59 min
Interface Port	Ethernet, RS-232, RS-485
No. of syringes	Standard: 12
Stirrer Drive	High performance BLDC drive (maintenance-free)
Vessel capacity	150ml, 250ml, 1L
Water bath	One piece moulded clear
Bath capacity	22 L
Temperature Controller	Heater: 1kW, SS 316 with thermal protection
Dimensions (W x L x H)	517 mm x 1100 mm x 690 mm
Electrical power	220/230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/6 Amp
Weight	141 Kg

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Tablet Dissolution Tester TrustE-08 with Sample Collector (Dc)

Ordering information:

Products		
Product code	Description	Quantity
01210404	TrustE-08 with Dc - 230V	1
01110404	TrustE-08 with Dc - 110V	1

Standard accessories		
Product code	Description	Quantity
0105A00075	EDT laser marked basket mesh 40 (Set of 6)	1
0104A00066	EDT Basket rod	8
0104A00818	Secure Clip Basket Holder	6
0104A00152	EDT Moulded Tray For Basket, Basket Holder & Rod	1
0104A00091	EDT Laser Marked Paddles For 1000 ML Glass Vessel (Set of 6)	1
0104A00093	EDT Moulded Tray For Paddles	1
0104A00730	Hexa Vessel 1000ml Clear (Set of 6)	1
0104A00242	EDT Merlon Vessel 1000ml - Clear (N)	2
0104BV4114	EDT-Lx V-4 Temperature Controller ETC-15	1
0104A00610	Stirrer for replenishment	2
0104BV4155	EDT V-5 External Probe Without Trim Pot	1
	IntegrOs Lite Version	1

Offline sampling accessories		
Product code	Description	Quantity
0168A00001	Motorised sampling manifold	1
0104A00322	PTFE carrier tubing	6
0104A00334	Syringe manifold	1
0134B00001	Individual vessel temperature sensors with cable and connector	1
0114BSP049	ESP-64 Syringe Pump (with 10 ml syringe)	1
0104A00698	EDT disposable filters (Packet of 50)	1
0121ASC010	ESC-06Dc inert disposable polypropylene vials 20ml (Set of 96)	1
0172ASC017	Sample collector vial tray - 20ml vials	1
0104800453	Auto Dispenser Plate	1
0104B00048	EDT dispenser plate holder (Set of 2)	1
1620200	EDT Sample collector model ESC-6DC with tygon tube SS sleeves	1

Optional accessories		
Product code	Description	Quantity
0104A00477	EDT ClearView™ (230V) Version-2.1	1
0104A00445	Biometrics	1
0909A00008	40 column embedded printer	1
0171A00001	Vibro	1
0104A00785	Connexion+	1
0104A00401	EDT laser marked glass vessel 1000ml - Amber	1
0104A00259	EDT merlon vessel 1000ml amber (N)	1
0104A00440	IntegrOS standard version	1
0104A00438	IntegrOS pro version with 21 CFR	1
1071A00017	TrustE-14 tablet Dispenser plate with provision for manual sampling	1
0104A00179	10 Mesh Customized Stationary Basket (Set Of 6)	1
0171A00005	Auto dispenser plate	1
0104A00054	EDT Basket 10 Mesh (Set of 6)	1
0104A00055	EDT Basket 20 Mesh (Set of 6)	1
0104A00280	EDT Basket 20 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate	1
0104A00252	EDT Laser Marked Basket 20 Mesh	1
0104A00056	EDT Basket 30 Mesh (Set of 6)	1
0104A00057	EDT Basket 40 Mesh (Set of 6)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 | sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 | www.electrolabgroup.com

Tablet Dissolution Tester Trust€-08 with Sample Collector (Dc)

Ordering information:

Optional accessories		
Product code	Description	Quantity
0104A00075	EDT Laser Marked Basket 40 Mesh (Set of 6)	1
0104A00061	EDT Gold Plated Basket 40 Mesh	1
0104A00204	EDT Basket 40 Mesh 24.5 mm (Set of 6) with certificate	1
0104A00058	EDT Basket 60 Mesh (Set of 6)	1
0104A00059	EDT Basket 80 Mesh (Set of 6)	1
0104A00060	EDT Basket 100 Mesh (Set of 6)	1
0104A00069	EDT Basket Holder with 'O' Ring (Set of 6)	1
0104A00206	EDT Holder for Basket 40 Mesh 24.5 mm (Set of 6)	1
0104A00068	EDT Basket Holder with Clip (Set of 6)	1
0104A00299	EDT Basket Holder with 'O' Ring (Set of 6) for 150ml Vessel	1
0104A00067	EDT Basket 'O' Ring (Set of 6)	1
0104A00071	EDT Basket LID 40 Mesh (Set of 6)	1
0104A00070	EDT Basket LID 10 Mesh (Set of 6)	1
0104A00062	EDT Stationary Basket (Set of 6)	1
0104A00772	Stationary Suppository Basket SS (Set of 6)	1
0104A00064	EDT Felodipine Basket with Spacer (Set of 6)	1
0104A00065	EDT Moulded Tray For Basket	1
0104A00097	EDT Sinkers 11 mm ID x 31 mm OD (Set of 6)	1
0104A00096	EDT Sinkers 8 mm ID x 23 mm OD (Set of 6)	1
0104A00094	EDT Sinkers 5 mm ID x 15.5 mm OD (Set of 6)	1
0104A00095	EDT Sinkers 6 mm ID x 18 mm OD (Set of 6)	1
0104A00098	EDT Japanese Sinkers (Set of 6)	1
0104A00245	EDT Customized Japanese Sinkers 15 mm ID X 25-26 mm Height	1
0104A00218	EDT Customized Japanese Sinkers 15 mm ID X 40 mm Height (Set of 6)	1
0104A00220	EDT Japanese Sinkers with 40 Mesh Cover (Set of 6)	1
0104A00286	Customized Japanese Sinkers SS316 Length 25mm x Inner Dia 12mm x Pore Size 3.5mm (Set Of 6)	1
0104A00101	EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big (Set of 6)	1
0104A00102	EDT Spring Type Sinkers 9 mm ID x 30 mm OD - Small (Set of 6)	1
0104A00588	Spring type sinkers 0.8 x 6.8 x 3 (set of 6)	1
0104A00590	Customized Spring Type Sinkers 1 x 10 x 25 (set of 6)	1
0104A00743	Spiral sinkers (4 spiral, SS 316 26 mm length x 10 mm ID (Set of 6)	1
0104A00103	Teflon Coated Sinkers -Spring Type (Big) (Set of 6)	1
0104A00110	EDT Sinkers 'U' Type (Set of 6)	1
0104A00115	EDT Helix Coil Sinkers - CAPWHT-02 (Set of 6)	1
0104A00116	EDT Helix Coil Sinkers - CAPWHT-2S (Set of 6)	1
0104A00112	EDT Helical Coil Sinkers	1
0104A00786	Customize Helix Sinkers 11 mm ID X 31 mm OD With 3 To 4 Helix (Set of 6)	1
0104A00099	EDT Cage Type Sinkers - Big (Set of 6)	1
0104A00100	EDT Cage Type Sinkers - Small (Set of 6)	1
0104A00105	EDT 3 Prong Sinkers - SS (Set of 6)	1
0104A00104	EDT 3 Prong Sinkers - Plastic (Set of 6)	1
0104A00113	EDT Metformin Sinkers - Big (Set of 6)	1
0104A00107	EDT Alfuzosin Sinkers Sinkers 37.5 mm ID x 20 mm OD (Set of 6)	1
0104A00114	EDT Metformin Sinkers - Small (Set of 6)	1
1401A00077	Biphasic Dissolution Paddle (Set of 6)	1
0169A00013	EDT-BTv Laser Marked V-Stir Paddle for 1000mL Glass Vessel (set of 14)	1
0169A00002	EDT-BTv Laser Marked V-Stir Paddle for 1000mL Glass Vessel	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Tablet Dissolution Tester TrustE-08 with Sample Collector (Dc)

Ordering information:

Optional accessories		
Product code	Description	Quantity
0133A00003	EDT Laser Marked Paddle For 4 Ltr Vessel	1
0122A00020	EDT Teflon Paddles for 2 Ltr Vessel (Set of 6)	1
0122A00003	EDT Laser Marked Paddles For 2 Ltr Glass Vessel	1
0104A00789	56mm Paddle Over Disc With 40Mesh (Set Of 6)	1
0104A00783	Paddle Over Disc 80mm Watch glass Assembly with 120 Mesh Screen (Set of 6)	1
0104A00774	EDT Mini Paddles For 150ml Glass Vessel Each	1
0104A00765	EDT 76mm Paddle Over Disc Assembly With Standard Mesh (Set of 6)	1
0104A00704	EDT 55mm Paddle Over Disc Assembly With 20 Mesh (set of 6)	1
0104A00703	EDT 55mm Paddle Over Disc Assembly With 10 Mesh (set of 6)	1
0104A00702	EDT 33mm Paddle Over Disc Assembly With 10 Mesh (set of 6)	1
0104A00699	EDT 33mm Paddle Over Disc Assembly With 20 Mesh (set of 6)	1
0104A00673	Paddle For 500ml vessel	1
0104A00630	EDT Laser Marked Paddles For 1000ml Glass Vessel	1
0114ASP009	ESP-64/84 FEP carrier tube from Disso to Syringe pump (Set of 8)	1
0114ASP001	ESP Glass syringe 25ml	1
0114ASP007	ESP-64/84 FEP carrier tube from syringe pump to Air (Set of 8)	1
0114ASP005	ESP-64/84 FEP carrier tube from syringe pump to Replenishment (Set of 8)	1
0121ASC011	ESC-DC Tygon Tube with SS sleeve (Set of 6)	1
0121ASC012	ESC-DC PTFE Carrier Tube from SC to needle (Set of 6)	1
0104A00516	Vibro for Lx with EKF (Pro)	1
0104A00517	Vibro for Lx with EKF (Standard)	1

Breakdown accessories		
Product code	Description	Quantity
0171B00008	Trust E-08 Molded top cover	1
0171B00030	Trust E-08 Diving Manifold Assembly with 12 Temperature Probe	1
0104BV4229	Trust E-08 Platinum Kit Mechanical Spares	1
0104BV4230	Trust E-08 Platinum Kit Electronic Spares	1
0171B00022	Trust E-08 Diving Manifold Assembly without Temperature Probe and without Sampling Cannula	1
0171B00021	Trust E-8/14 Lift Park Position Opto Sensor	1
0114BSP052	ESP-124 Drive	1
0114BSP005	ESP syringe Opto PCB	1
0114BSP001	ESP-64/84 Processor card	1
0114BSP002	ESP Driver card	1
0114BSP025	ESP Valve Motor Pulley	1
0114BSP019	ESP-124/144 PTFE carrier tube from Disso to Syringe Pump (Set of 14)	1
0114BSP003	ESP Valve Opto PCB	1
0114BSP034	ESP-64/84 Driver card new (SMD)	1
0120BSP001	ESP-124/144 Processor card	1
0114BSP044	ESP Valve (4 way) with bearing housing assembly with pulley for mounting Opto card	1
0114BSP053	ESP-64 Drive	1
0114BSP011	ESP Fuse Set	1
0121BSC048	ESC DC keyboard with panel	1
0114BSP008	ESP Transformer	1
0121BSC045	ESC-DC Processor card 1901003	1
0121BSC060	ESC-DC-Front/Back Opto Card	1
0121BSC059	ESC-DC Proactive Kit	1
0121BSC046	ESC-DC Drive TB 6600	1
0121BSC052	ESC-DC Opto Cutter (Front Back Position)	1
0121BSC051	ESC-DC Side Cover Acrylic	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Tablet Dissolution Tester TrustE-08 with Sample Collector (Dc)

Ordering information:

Breakdown accessories		
Product code	Description	Quantity
0121BSC049	ESC-DC waste tray holder for 8 station	1
0171B00024	Trust E-08 Top cover with LED	1

Assurance Plan		
Product code	Description	Quantity
01104S01	Platinum Assurance Plan for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1
01104S02	Gold Assurance Plan for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1
01104S03	Silver Assurance Plan for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1
01104S20	Bronze Assurance Plan for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

IQ/OQ Document		
Product code	Description	Quantity
01104S15	IQ/OQ Document for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Re-installation Charges		
Product code	Description	Quantity
01104S10	Re-installation Charges with IQ/OQ Document for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Performance Qualification (PQ Documents)		
Product code	Description	Quantity
0168S19	Performance qualification (PQ documents) for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Extended IQ/OQ Document		
Product code	Description	Quantity
	Extended IQ/OQ Document for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Performance Qualification		
Product code	Description	Quantity
	Performance Verification Test (PVT Assist) for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Inspection & Service Charges		
Product code	Description	Quantity
	Inspection & Service Charges for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Upgradation Charges		
Product code	Description	Quantity
	Upgradation Charges for Tablet Dissolution Tester Model TrustE-08 with Dc Sample Collector	1

Single validation Charges		
Product code	Description	Quantity
01104S08	Single validation charges for Tablet Dissolution Tester Model TrustE-14 with Dc Sample Collector	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com