

Air Jet Sieve Shaker Shaktus Air

ELECTROLAB Datasheet



Features:

- 5" touch screen display
- Robust Design
- Complies with USP, pH. Eur, IP
- Integrity operating System Integros 2.0
- Basic Security
- User friendly interface
- 40 Column built in thermal printer (Optional)
- Suitable for dry sieving
- Very low noise and low maintenance
- Weighing balance connectivity with Mettler Toledo, Sartorius, Contech, Ohaus
- Programmable vacuum pressure in multiple units
- Biometric (Optional)

Specifications:

Model	Shaktus Air
Test Speed	10 - 80 RPM
Sieve Diameter	200 mm (fixed)
No. of sieves	1 - 10
Sample Weight	10 - 100 gm / cycle
Pressure / unit	Pa / InH2O / mmHg / InHg / psi
Electrical power	230 VAC / 50 Hz or 110 VAC / 50 Hz $\pm$ 10%
Dimensions (L x W x H)	430mm x 285mm x 264mm
Weight	22 Kg

Ordering Information

Products		
Product code	Description	Quantity
0540500	Shaktus Air (230V)	1
0530500	Shaktus Air (110V)	1
Standard Accessories		
Product code	Description	Quantity
0501A00196	Dry Vacuum Cleaner	1
	Wand	1
Optional Accessories		
Product code	Description	Quantity
0501A00195	High Efficiency Cyclone Model VME 2.1 For Air Jet Sieves Shaker	1
0501A00221	Shaktus Air Sieve Cover	1
	Sieve	1
IQ/OQ Documents		
Product code	Description	Quantity
	IQ/OQ Document for Shaktus Air	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupa Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

License Vs Limits

<u>Shaktus Air</u>	Basic	Network	Advance Reporting	CFR	OPC
Users	50	50	250	250	250
Roles	10	10	10	10	10
Audit Trails	-	-	-	2,00,000	2,00,000
Report Template	-	-	10	10	10
Reports	1 (L.R.)	1 (L.R.)	50,000	50,000	50,000
Products	20	20	1,000	1,000	1,000
Logbook	-	-	2,00,000	2,00,000	2,00,000

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupa Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com