

Electromagnetic Sieve Shaker EMS-8 Monitor



Features:

- Ideal for particle size analysis and separation
- Tri-dimensional sifting motion
- Intermittent and continuous shifting motion with 15 programmable levels of amplitude
- Very low noise operation with low maintenance
- Programmable shake time from 1 min to 99 min
- Suitable for dry and wet sieving
- The display shows the amplitude level
- Instrument is powered by an electromagnetic drive which has no rotating parts, making it maintenance free and quiet in operation

Specifications:

Model	EMS-8 Monitor
Dry sieving	Standard
Mode of operation	Continuous and intermittent
Intermittent operation	At intervals (0.5 sec)
Capacity	Up to 8 sieves of 200 mm dia × 50 mm h Up to 16 sieves of 100 mm dia × 25 mm h Suitable for sieve sizes from 28 µm to 4 mm
Shake time	Programmable from 1 min to 99 min
Power level	Programmable from 5 to 20
Display	16 × 2 character alphanumeric LCD
Product Category	Particle Size Analyzer
Noise level	< 61 dB without sieves at maximum amplitude < 71 dB with sieves and material at maximum amplitude
Amplitude level	Monitor the amplitude through the display
Power	220 / 230 V AC, 50 Hz, 600 VA
Dimensions (W × L × H)	312 mm × 339 mm × 270 mm
Weight	50 kgs without sieves

Ordering information:

Products		
Product code	Description	Quantity
0540300	Electromagnetic Sieve Shaker-EMS 8 Mon(SS)	1
0520300	Electromagnetic Sieve Shaker-EMS 8 Mon	1

Standard Accessories		
Product code	Description	Quantity
0501A00001	EMS-8 Backlight Knob To Fit Locking Arm (Set Of 2)	1
0501A00065	EMS-8 Vertical Rod (Set Of 2)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Electromagnetic Sieve Shaker

EMS-8 Monitor

Optional Accessories		
Product code	Description	Quantity
0501A00220	EMS-8 Lid For 200 mm Dia Test Sieves	1
0501A00217	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 20 - 810 µm) Pharma Grade	1
0501A00216	EMS-8 Test Sieves of 200mm Dia (4.75mm)	1
0501A00215	EMS-8 Test Sieves of 200mm Dia (6.3mm)	1
0501A00214	EMS-8 Test Sieves of 200mm Dia (8mm)	1
0501A00213	EMS-8 Test Sieves of 200mm Dia (10mm)	1
0501A00212	EMS-8 Test Sieves of 200mm Dia (12mm)	1
0501A00211	EMS-8 Test Sieves of 200mm Dia (2mm)	1
0501A00209	EMS-8 Test Sieves Of 200 mm Dia(5.6mm Micron)	1
0501A00208	EMS-8 Test Sieves Of 200 mm Dia(6.3mm Micron)	1
0501A00207	EMS-8 Test Sieves Of 200 mm Dia(6.7mm Micron)	1
0501A00206	EMS-8 Test Sieves Of 200 mm Dia(7.1mm Micron)	1
0501A00205	EMS-8 Test Sieves Of 200 mm Dia(Micron Size 8MM)	1
0501A00204	Kloudface for EMS-8 integrOS Pro Interface Model EKF-06 with life time year License For 110V	1
0501A00203	Kloudface for EMS-8 integrOS Pro Interface Model EKF-06 with (3+1) year License For 110V	1
0501A00202	Kloudface for EMS-8 integrOS Pro Interface Model EKF-06 with 1 year License For 110V	1
0501A002016	EMS-8 Test Sieves Of 200 mm Dia(20Micron)	1
0501A002015	EMS-8 Test Sieves Of 200 mm Dia(25Micron)	1
0501A002014	EMS-8 Test Sieves Of 200 mm Dia(32Micron)	1
0501A002013	EMS-8 Test Sieves Of 200 mm Dia(36Micron)	1
0501A002012	EMS-8 Test Sieves Of 200 mm Dia(4.5mm Micron)	1
0501A002011	EMS-8 Test Sieves Of 200 mm Dia(4.75mm Micron)	1
0501A002010	EMS-8 Test Sieves Of 200 mm Dia(5.0mm Micron)	1
0501A00194	EMS-8 Wet Sieve Attachment For 100 mm Dia	1
0501A00193	EMS-8 Amplitude Sticker	1
0501A00192	EMS-8 Lid For 200 mm Dia Test Sieves	1
0501A00191	EMS-8 Receiver For 200 mm Dia Test Sieves	1
0501A00190	EMS-8 Rubber Leg (Set of 4)	1
0501A00189	EMS-8 Rubber Pad (Set of 2)	1
0501A00188	EMS-8/8 Plus S.S. Body	1
0501A00187	EMS-8 Test Sieves Of 100 mm Dia (28 µm) Pharma Grade	1
0501A00186	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 400 - 38 µm) Pharma Grade	1
0501A00185	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 325 - 45 µm) Pharma Grade	1
0501A00184	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 270 - 53 µm) Pharma Grade	1
0501A00183	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 230 - 63 µm) Pharma Grade	1
0501A00182	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 200 - 75 µm) Pharma Grade	1
0501A00181	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 170 - 90 µm) Pharma Grade	1
0501A00180	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 140 - 106 µm) Pharma Grade	1
0501A00179	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 120 - 125 µm) Pharma Grade	1
0501A00178	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 100 - 150 µm) Pharma Grade	1
0501A00177	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 80 - 180 µm) Pharma Grade	1
0501A00176	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 70 - 212 µm) Pharma Grade	1
0501A00175	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 60 - 250 µm) Pharma Grade	1
0501A00174	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 50 - 300 µm) Pharma Grade	1
0501A00173	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 45 - 355 µm) Pharma Grade	1
0501A00172	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 40 - 425 µm) Pharma Grade	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

☎ +91 - 22 - 4041 3131 ✉ sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

☎ +91 - 22 - 4161 3122 🌐 www.electrolabgroup.com

Electromagnetic Sieve Shaker

EMS-8 Monitor

Optional Accessories		
Product code	Description	Quantity
0501A00171	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 20 - 810 µm) Pharma Grade	1
0501A00170	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 30 - 600 µm) Pharma Grade	1
0501A00169	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 25 - 710 µm) Pharma Grade	1
0501A00168	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 20 - 850 µm) Pharma Grade	1
0501A00167	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 18 - 1.00 mm) Pharma Grade	1
0501A00166	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 16 - 1.18 mm) Pharma Grade	1
0501A00165	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 14 - 1.40 mm) Pharma Grade	1
0501A00164	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 12 - 1.70 mm) Pharma Grade	1
0501A00163	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 10 - 2.00 mm) Pharma Grade	1
0501A00162	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 8 - 2.36 mm) Pharma Grade	1
0501A00161	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 7 - 2.80 mm) Pharma Grade	1
0501A00160	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 6 - 3.35 mm) Pharma Grade	1
0501A00159	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 5 - 4.00 mm) Pharma Grade	1
0501A00158	EMS-8 Test Sieves Of 200 mm Dia (Micron 28 µm) Pharma Grade	1
0501A00157	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 400 - 38 µm) Pharma Grade	1
0501A00156	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 325 - 45 µm) Pharma Grade	1
0501A00155	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 270 - 53 µm) Pharma Grade	1
0501A00154	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 230 - 63 µm) Pharma Grade	1
0501A00153	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 200 - 75 µm) Pharma Grade	1
0501A00152	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 170 - 90 µm) Pharma Grade	1
0501A00151	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 140 - 106 µm) Pharma Grade	1
0501A00150	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 120 - 125 µm) Pharma Grade	1
0501A00149	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 100 - 150 µm) Pharma Grade	1
0501A00148	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 80 - 180 µm) Pharma Grade	1
0501A00147	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 70 - 212 µm) Pharma Grade	1
0501A00146	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 60 - 250 µm) Pharma Grade	1
0501A00145	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 50 - 300 µm) Pharma Grade	1
0501A00144	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 45 - 355 µm) Pharma Grade	1
0501A00143	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 40 - 425 µm) Pharma Grade	1
0501A00142	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 35 - 500 µm) Pharma Grade	1
0501A00141	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 30 - 600 µm) Pharma Grade	1
0501A00140	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 25 - 710 µm) Pharma Grade	1
0501A00139	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 20 - 850 µm) Pharma Grade	1
0501A00138	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 18 - 1.00 mm) Pharma Grade	1
0501A00137	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 16 - 1.18 mm) Pharma Grade	1
0501A00136	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 14 - 1.40 mm) Pharma Grade	1
0501A00135	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 12 - 1.70 mm) Pharma Grade	1
0501A00134	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 10 - 2.00 mm) Pharma Grade	1
0501A00133	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 8 - 2.36 mm) Pharma Grade	1
0501A00132	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 7 - 2.80 mm) Pharma Grade	1
0501A00131	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 6 - 3.35 mm) Pharma Grade	1
0501A00130	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 5 - 4.00 mm) Pharma Grade	1
0501A00129	EMS-8 Test Sieves Of 100 mm Dia (28 µm)	1
0501A00128	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 400 - 38 µm)	1
0501A00127	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 325 - 45 µm)	1
0501A00126	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 270 - 53 µm)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

☎ +91 - 22 - 4041 3131 ✉ sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

☎ +91 - 22 - 4161 3122 🌐 www.electrolabgroup.com

Optional Accessories		
Product code	Description	Quantity
0501A00125	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 230 - 63 µm)	1
0501A00124	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 200 - 75 µm)	1
0501A00123	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 170 - 90 µm)	1
0501A00122	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 140 - 106 µm)	1
0501A00121	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 120 - 125 µm)	1
0501A00120	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 100 - 150 µm)	1
0501A00119	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 80 - 180 µm)	1
0501A00118	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 70 - 212 µm)	1
0501A00117	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 60 - 250 µm)	1
0501A00116	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 50 - 300 µm)	1
0501A00115	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 45 - 355 µm)	1
0501A00114	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 40 - 425 µm)	1
0501A00113	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 35 - 500 µm)	1
0501A00112	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 30 - 600 µm)	1
0501A00111	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 25 - 710 µm)	1
0501A00110	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 20 - 850 µm)	1
0501A00109	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 18 - 1.00 mm)	1
0501A00108	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 16 - 1.18 mm)	1
0501A00107	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 14 - 1.40 mm)	1
0501A00106	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 12 - 1.70 mm)	1
0501A00105	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 10 - 2.00 mm)	1
0501A00104	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 8 - 2.36 mm)	1
0501A00103	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 7 - 2.80 mm)	1
0501A00102	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 6 - 3.35 mm)	1
0501A00101	EMS-8 Test Sieves Of 100 mm Dia (ASTM No. 5 - 4.00 mm)	1
0501A00100	EMS-8 Test Sieves Of 200 mm Dia (Micron 28 Um)	1
0501A00099	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 400 - 38 µm)	1
0501A00098	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 325 - 45 µm)	1
0501A00097	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 270 - 53 µm)	1
0501A00096	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 230 - 63 µm)	1
0501A00095	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 200 - 75 µm)	1
0501A00094	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 170 - 90 µm)	1
0501A00093	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 140 - 106 µm)	1
0501A00092	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 120 - 125 µm)	1
0501A00091	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 100 - 150 µm)	1
0501A00090	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 80 - 180 µm)	1
0501A00089	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 70 - 212 µm)	1
0501A00088	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 60 - 250 µm)	1
0501A00087	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 50 - 300 µm)	1
0501A00086	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 45 - 355 µm)	1
0501A00085	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 40 - 425 µm)	1
0501A00084	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 35 - 500 µm)	1
0501A00083	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 30 - 600 µm)	1
0501A00082	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 25 - 710 µm)	1
0501A00081	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 20 - 850 µm)	1
0501A00080	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 18 - 1.00 mm)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

☎ +91 - 22 - 4041 3131 ✉ sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

☎ +91 - 22 - 4161 3122 🌐 www.electrolabgroup.com

Electromagnetic Sieve Shaker

EMS-8 Monitor

Optional Accessories		
Product code	Description	Quantity
0501A00079	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 16 - 1.18 mm)	1
0501A00078	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 14 - 1.40 mm)	1
0501A00077	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 12 - 1.70 mm)	1
0501A00076	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 10 - 2.00 mm)	1
0501A00075	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 8 - 2.36 mm)	1
0501A00074	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 7 - 2.80 mm)	1
0501A00073	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 6 - 3.35 mm)	1
0501A00072	EMS-8 Test Sieves Of 200 mm Dia (ASTM No. 5 - 4.00 mm)	1
0501A00071	EMS-8 Locking Arm (Horizontal) For 100 mm Dia Test Sieves	1
0501A00070	EMS-8 Allenky Screw - 6 mm	1
0501A00069	EMS-8 Communication Cable	1
0501A00068	EMS-8 Sieve Easy Software Cd	1
0501A00067	EMS-8 Lid For Seives (NEW)	1
0501A00006	EMS-8Wet Sieve Attachment For 200 mm Dia	1
0501A00005	EMS-8 Lid & Receiver For 100 Mm Dia Test Sieves	1
0501A00004	EMS-8 Lid & Receiver For 200 mm Dia Test Sieves	1
0501A00002	EMS-8Locking Arm (Horizontal) For 200mm Dia Test Sieves	1

Breakdown Accessories		
Product code	Description	Quantity
0502B00014	EMS-8 Driver Card For 110V	1
0501B00019	EMS-8 Plus Metal Resistor 600E/100W	1
0501B00018	EMS-8 Platinum Kit	1
0501B00017	EMS-8 Base Cover with Front Panel (SS)	1
0501B00016	EMS-8 - 4Pin Relimate Cable (Display Card To RS232 Card)	1
0501B00015	EMS-8 - 4Pin Relimate Cable (Display Card To Drv. Card)	1
0501B00014	EMS-8 Mains Cord 5 Amp	1
0501B00012	EMS-8 RS-232 Communication Card	1
0501B00011	EMS-8 Casting Pillar (Set Of 4)	1
0501B00010	EMS-8 Base Plate	1
0501B00009	EMS-8 Power Supply On/Off Switch	1
0501B00008	EMS-8 Acrylic Plate (Set of 3)	1
0501B00007	EMS-8 Coil	1
0501B00006	EMS-8 Front Panel with Keypad	1
0501B00005	EMS-8 Display Card	1
0501B00004	EMS-8 LCD Display	1
0501B00003	EMS-8 Driver Card	1
0501B00002	EMS-8 Transformer	1
0501B00001	EMS-8 Power Supply Socket with Filter	1

IQ/OQ Document		
Product code	Description	Quantity
040300S15	IQ/OQ Document For Electrolab Electro Magnetic Sieve Shaker-EMS 8 Mon	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com