

Tablet Friability Tester Fritus-L



Fritus-L

Features:

- Programmable Time & Count Modes
- USP-Compliant 10° Drum Tilt
- Unique Front Loading Drum Design
- Automatic Sample Discharge after Test Completion
- Automatic Drum Positioning for Easy Reloading
- Supports AD Drum™ (Patented) & Abrasion Drum
- Snap-Click Drum Locking Knob
- Maintenance-Free Drive Mechanism
- Gentle Start & Stop Drum Operation
- No Drum Removal Required for Loading/Unloading
- Memory Retains Programmed Settings
- 7-Segment Display with Simple Keypad Operation
- Automatic Drum Detection & Error Indication
- Compact, User-Friendly Design

Specifications:

Model	Fritus-L
Speed	Fixed 60 RPM
Speed Accuracy	±1 RPM
Test Modes	Time Mode / Count Mode
Timer	Hrs:Min 0.01 to 9.59 / Min:Sec 00.01 to 99.59
Counter	1 to 9999 revolutions
Drum Type	Dual Drum
Power Supply	230/110 VAC, 50/60 Hz, Single Phase
Dimensions (L x W x H)	288mm x 289mm x 489mm
Weight	10.7 Kg with Drum

Ordering Information

Products		
Product code	Description	Quantity
0440800	Tablet Friability Tester Model Fritus-L (230V)	1
0430800	Tablet Friability Tester Model Fritus-L (110V)	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com