

Tablet Dissolution Tester Inspire-14 with Sample Collector (Dc)



Inspire-14 with Sample Collector (Dc)

Features:

- LCD display of 40 characters and 4 lines for user friendly operation
- Multilevel security with password
- Alphanumeric data entry of product and user details
- Unique TimeAction™ function which allows media changeover and infinity test run, ideal for sustained and controlled release products
- Complies with USP, IP, Ph. Eur. and JP specifications
- Automatic stirrer height positioning for USP 1, 2, 5, 6, intrinsic test methods
- Easy Snap-Fit™ shaft locking mechanism for positive engagement and wobble free operation
- Clear water bath with metal top-plate
- Precise individual vessel centering system
- Adjustable leveling of instrument
- Vibration free smooth motorized lift movement
- Isolated water circulating pump for precise temperature control of the water bath
- 20 programmable protocols
- 16 programmable sampling intervals
- Easy to empty and clean water bath with bottom drainage
- Designed to minimize routine validation
- Report for online validation of test, monograph, error log during the test
- Programmable wake-up for temperature control
- Programmable sample volume from 1 ml to 20 ml
- Audio-visual indicators for system status
- Power failure detection facility
- Split Rinse to avoid carry over
- Replenishment
- Individual vessel temperature monitoring

Specifications:

Model	Inspire-14 with DC
Speed range	Variable speed from 20 to 300 RPM
Speed accuracy	±1 RPM
Temperature range	20° C to 40° C (Extended range 10°C to 40°C & 20°C to 55°C optional)
Temperature accuracy	±0.1° C
Display	40 x 4 character backlighted LCD
No. of protocols	20
No. of samples	16
Sample interval	Maximum 16 of 1 min to 23 hrs, 59 min
Printer interface	Serial port
Printer supported	80 column Dot Matrix printer
Interface port	RS 485/RS 232
Vessel capacity	150mL, 250mL, 500mL, 1L
Water bath	One piece moulded clear bath
Temperature controller	Heater : 1kW, SS 316 with thermal protection Sensor : RTD
Electrical power	220/230 VAC, 50/60 Hz, 1.3 kW/6 Amp 110 VAC, 50/60 Hz, 1.3 kW/10 Amp
Dimensions (W x L x H)	521 x 1672 x 663 mm
Weight	212 kg

Ordering information:

Products		
Product code	Description	Quantity
01210204	Tablet Dissolution Tester Model Inspire-14 with Dc -230V	1
01110204	Tablet Dissolution Tester Model Inspire-14 with Dc -110V	1

Standard accessories		
Part number	Description	Quantity
0104A00075	EDT laser marked basket - 40 mesh (Set of 6)	2
0104A00066	EDT Basket Rod	12
0104A00068	EDT Basket Holder with clip (Set of 6)	2
0104A00592	Dissolution Inspire Detachable Paddle Blade (Set of 6)	2
0104A00610	Stirrer for replenishment	2
0104A00242	EDT merlon vessel -1000ml clear (N)	14
0104BV4155	EDT V-5 External Probe Without Trim Pot	1
0104BV4114	EDT- Lx V-4 Temperature Controller ETC-15	1

Offline sampling accessories		
Part number	Description	Quantity
0105A00277	Motorised sampling manifold	1
0105A00322	PTFE carrier tubing	12
0105B00101	Syringe filter manifold	1
0134B00001	Individual vessel temperature sensors with cable and connector	12
0120BSP046	ESP-124 syringe pump (With 10mL syringe)	1
0121BSC057	EDT Sample Collector Model – ESC-12DC - with tygon tube ss sleeve	1
0114ASC013	ESC Inert Disposable Polypropylene Vials 20mL (Pack of 100)	2
0104A00698	Disposable filters (Packet of 100)	1

Assurance Plan		
Part number	Description	Quantity
0183S03	Silver Assurance Plan for Tablet Dissolution Tester Model Inspire-14 with DC Sample Collector	1
0183S01	Platinum Assurance Plan for Tablet Dissolution Tester Model Inspire-14 with DC Sample Collector	1
0183S02	Gold Assurance Plan for Tablet Dissolution Tester Model Inspire-14 with DC Sample Collector	1
0183S20	Bronze Assurance Plan for Tablet Dissolution Tester Model Inspire-14 with DC Sample Collector	1

IQ/OQ Document		
Part number	Description	Quantity
0183S15	IQ/OQ Document for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Inspection & Service Charges

Part number	Description	Quantity
01102S15	Inspection & Service Charges for Tablet Dissolution Model Inspire-14 with Dc Sample Collector	1

Performance Qualification

Part number	Description	Quantity
01204S17	Performance Qualification for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Single Validation Charges

Part number	Description	Quantity
01102S04	Single Validation Charges for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Extended Warranty Support Plan

Part number	Description	Quantity
01102S31	Extended Warranty Support Plan for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Performance Verification Test

Part number	Description	Quantity
01102S14	Performance Verification Test (PVT Assist) for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Reinstallation Charges

Part number	Description	Quantity
01102S11	Reinstallation Charges with IQ/OQ Document for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

Dismantling Charges

Part number	Description	Quantity
01102S12	Dismantling Charges for Tablet Dissolution Tester Model Inspire-14 with Dc Sample Collector	1

*The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131  sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122  www.electrolabgroup.com

	Optional accessories	
SR NO	Optional accessories	
1	EDT- 406L S.S. Bath Top Plate	0133A00007
2	EDT Kone Vessel 1000 ml – Amber With Sigma Delta Certificate	0104A00478
3	EDT 20 Mesh Stationary Basket (Set of 6)	0104A00583
4	EDT-Inspire Syringe manifold Tygon Tube with SS Sleeve (Set Of 6)	0104A00706
5	EDT Rubber Holder For External Probe.(Set Of 6)	0104A00693
6	EDT 62mm Paddle Over Disc Assembly With 16 Mesh(set of 6)	0104A00609
7	EDT Moulded Jar LID - Plastic	0104A00625
8	EDT – Adjustable Legs Long (Set of 4)	0104A00631
9	EDT Teflon Paddles (set of 8)	0104A005893
10	EDT Laser Marked Paddle For 4 Ltr Vessel	0133A00003
11	EDT Paddle/Basket Rod	0104A00596
12	EDT 85mm Paddle Over Disc Assembly	0104A00362
13	EDT-14 Bath Top Plate For Hexa Vessel	0104A00781
14	EDT Adjustable Sampling Cannula for low evaporation jar lid (Set of 8)	0104A00048
15	EDT Apex Glass Vessel Clear 2Ltr. With Sigma Delta Certificate	0104A0058945
16	EDT Gasket For Tablet Dispensor Plate	0104A00272
17	EDT Helical Coil Sinkers	0104A00112
18	EDT Adjustable Sampling Cannulas For Tablet Input Device (Set of 6) Length 144 mm	0104A00046
19	EDT Lid for Basket 40 Mesh 24.5 mm (Set of 6)	0104A00205
20	EDT Inside Caliper 6"	0104A00182
21	EDT 90mm Watch Glass Disc (Set of 6)	0104A00782
22	EDT-14L/14Lx Water Bath with Plastic QRC - Clear	0105A00143
23	EDT Laser Marked Glass Vessel 1000 MI - Clear With Sigma Delta Certificate	0104A00400
24	EDT Offline Accessories for 8 Station Tablet Dissolution Tester (Dx)	0104A00488
25	EDT Depth Gauge for 2 Ltr Vessels	0122A00010
26	EDT Basket 40 Mesh 24.5 mm (Set of 6)	0104A00204
27	EDT Intrinsic Attachment	0104A00128
28	EDT Plastic Quick Release Coupling with Shut Off Valve (Set of 2)	0104A00146
29	EDT Glass Vessel 1000 MI - Amber With Sigma Delta Certificate	0104A00399
30	EDT Glass Vessel 1000 MI – Clear With Sigma Delta Certificate	0104A00398
31	EDT Kone Vessel 1000ml - Amber	0104A00011
32	EDT Jar Lock for Glass Vessel Suitable for Composite Top Plate (Set of 16)	0104A00138
33	EDT-14Lx PTFE Bath Top Plate	0104A00614
34	EDT Bath Rod	0104A00262
35	EDT-14L/14Lx Aluminum Top Plate for 250ml Conversion Kit	0105A00395
36	EDT Digital Depth Gauge	0104A00155
37	EDT Customized Cage Sinkers 9.70mm Height X 13.50mm ID (Set of 6)	0104A00771
38	EDT Disposable Filters 35 Micron (Pack Of 100)	0104A00769
39	EDT Customized Sinkers 14mm Internal Diameter X 34 mm Internal Length (Set of 6)	0104A00421
40	EDT 33mm Paddle Over Disc Assembly	0104A00358
41	EDT 50mm Paddle Over Disc Assembly	0104A00359
42	EDT 55mm Paddle Over Disc Assembly	0104A00360
43	EDT 65mm Paddle Over Disc Assembly	0104A00361
44	EDT Laser Marked Mini Paddles for 100 MI Glass Vessel (Set of 6)	0104A00089
45	EDT Laser Marked Paddle for 33mm Disc (Set of 6)	0104A00353
46	EDT Laser Marked Paddle for 49mm Disc (Set of 6)	0104A00354
47	EDT Laser Marked Paddle for 55mm Disc (Set of 6)	0104A00355
48	EDT Laser Marked Paddle for 65mm Disc (Set of 6)	0104A00356
49	EDT Laser Marked Paddle for 85mm Disc (Set of 6)	0104A00357
50	EDT Laser Marked Paddles For 1000 MI Glass Vessel (Set of 6)	0104A00091
51	EDT Laser Marked Paddles For 2 Ltr Glass Vessel	0122A00003
52	EDT Mini Paddles For 100 MI Glass Vessel (Set of 6)	0104A00083
53	EDT Mini Paddles For 150 MI Glass Vessel (Set of 6)	0104A00293
54	EDT Mini Paddles For 250 MI Glass Vessel (Set of 6)	0104A00084
55	EDT Moulded Tray For Paddles	0104A00093
56	EDT Paddle Over Disc 33	0104A00168
57	EDT Paddle Over Disc 50 (Set of 6)	0104A00169
58	EDT Paddle Over Disc 55 (Set of 6)	0104A00170

59	EDT Paddle Over Disc 65 (Set of 6)	0104A00171
60	EDT Paddle Over Disc 85 (Set of 6)	0104A00172
61	EDT Paddles – 1000ml For Stationary Intrinsic Attachment (Set of 6)	0104A00317
62	EDT Paddles For 1000 ML Glass Vessel (Set of 6)	0104A00085
63	EDT Paddles For Enhancer Cell (Set of 6)	0104A00088
64	EDT Teflon Paddles for 2 Ltr Vessel (Set of 6)	0122A00020
65	EDT-208Lx/14Lx paddle without laser marking	0122A00025
66	EDT Paddles For 1000 ml Vessel (Set of 8)	0104A00526
67	EDT Paddles For 1000 ML Glass Vessel	0104A00568
68	EDT 90mm Watch Glass Disc Assembly With 120 Mesh Screen Without Paddle	0104A00576
69	EDT Teflon Paddle	0104A005894
70	EDT Laser Marked Paddles For 1000 ML Glass Vessel	0104A00630
71	EDT 55mm Paddle Over Disc Assembly With 120 Mesh(set of 6)	0104A00637
72	EDT 33mm Paddle Over Disc Assembly With 20 Mesh (set of 6)	0104A00699
73	EDT 33mm Paddle Over Disc Assembly With 10 Mesh (set of 6)	0104A00702
74	EDT 55mm Paddle Over Disc Assembly With 20 Mesh(set of 6)	0104A00704
75	EDT 55mm Paddle Over Disc Assembly With 10 Mesh(set of 6)	0104A00703
76	EDT 76mm Paddle Over Disc Assembly With Standard Mesh(Set Of 6)	0104A00765
77	EDT Digital Thermometer	0104A00124
78	EDT-14L/14Lx Tablet Dispensor Plate	0105A00159
79	EDT I-Disso 6U	0164A00001
80	EDT-14L/14Lx I-Disso 12U (230 V)	0105A00198
81	Laser Marking Charges For EDT Kone Vessel 1000ml	0104A00010S19
82	EDT Hybrid Cannula for Sampling Manifold (Set of 6)	0104A00452
83	EDT Basket LID 10 Mesh (Set of 6)	0104A00070
84	EDT-14L/EDT-14Lx 250 ML Low Volume Conversion Kit	0105A00213
85	EDT Glass Vessel 'V' Type 1000ml - Clear	0104A00016
86	EDT - Offline Accessories for 14 Station Tablet Dissolution Tester (DC)	0171A00003
87	EDT Basket Holder with Clip (Set of 6)	0104A00068
88	EDT Glass Vessel 'D' Type 1000ml - Clear	0104A00014
89	EDT Japanese Sinkers	0104A00622
90	EDT Sampling Cannula Rubber 'O' Ring for Dispensor plate (Set of 10)	0104A00463
91	EDT Sinkers 11 mm ID x 31 mm OD (Set of 6)	0104A00097
92	EDT Basket LID 40 Mesh (Set of 6)	0104A00071
93	EDT Digital Validation Kit	0104A00239
94	EDT-Lx V-4 Lift Disc Opto Card	0104BV4059
95	EDT Cooling Instrument (Temp. From 10° to 40° C)	0104A00330
96	EDT Jar Lock For Merlon Vessel For Composite Top Plate	0104A00139
97	EDT Teflon Coated Sinkers 22 mm Long & 13 cm Length of Gauge Wire	0104A00163
98	EDT Adjustable Sampling Cannulas For 500 ML	0104A00256
99	EDT Sampling Cannulas for 2 Ltr Vessel	0122A00006
100	EDT-14Lx Manual Syringe Sampling Attachment	0104A00640
101	EDT Butterfly Sinkers(Set of 6)	0104A00109
102	EDT KBR Bench Press	0104A00129
103	EDT Low Evaporation Jar Lid for Adjustable Cannulas Suitable for Glass Vessel (Set of 6)	0104A00234
104	EDT Merlon Vessel 1000ml - Clear (N)	0104A00242
105	EDT-14L/14Lx Bath Top Plate – Aluminum	0105A00311
106	EDTV-1234 Bath Temperature Sensor	0104B00013
107	EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 16)	0104A00133
108	EDT Silicon Tube Set For ETC-11	0104A00483
109	EDT Plunger Dial Gauge	0104A00482
110	EDT-14L/14Lx Light with supply Adaptor for I-Disso 12U	0104A00668
111	EDT Pin Sinkers Type - B (Set of 12)	0104A00672
112	EDT Blackout Curtain for I-Disso 06U	0164A00002
113	EDT Vernier Caliper for 250 ml Low Volume Conversion Kit	0104A00221
114	EDT Vibro	0104A00480
115	EDT-14L/14Lx Tablet Dispensor Plate with Felodipine Basket Provision	0105A00289
116	EDT Customized Japanese sinker ID 15mm x IL 32mm (Set of 6)	0104A00678
117	EDT-14L/14Lx IMMERSION (ENHANCER) CELL KIT- 0.5 CM2 For Acrylic Bath	0105A00420
118	EDT-14L/14Lx Aluminum Top Plate for Enhancer Cell Kit Acrylic Bath	0105A00387
119	EDT-14L/14Lx CPU For I-Disso 12U (230 V)	0104A00700
120	Tygon Tube Replacement And Performance Enhancement Kit For EDT-14Lx With D Sample Collector	0104A00682
121	EDT Basket Sinkers 40 Mesh (19mm Height) with Lid (set of 6)	0104A00667

122	EDT Laser Marked Glass Vessel 100 MI - Clear	0104A00020
123	EDT Laser Marked Glass Vessel 250 MI - Clear	0104A00021
124	EDT Laser Marked Kone Vessel 1000 MI - Clear	0104A00029
125	EDT Laser Marked Kone Vessel 1000 MI - Amber	0104A00030
126	EDT Laser Marked Basket 20 Mesh	0104A00252
127	EDT Basket 20 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate	0104A00280
128	EDT Luer Lock FTLLC-6 Male (Set of 8)	0104A00329
129	EDT Glass Vessel 4Ltr - Clear	0133A00001
130	EDT – Intrinsic Plate	0104A00278
131	EDT Merlon Vessel 1000 MI - Clear with Jar Lock Suitable for S.S. Top Plate	0104A00225
132	EDT Merlon Vessel 1 Ltr Amber With Jar Centring Ring For Composite Top Plate	0104A00258
133	EDT Sampling Cannulas Fixed Height 750 ml (Set of 6) Length 108 mm	0104A00157
134	EDT Sampling Cannulas Fixed Height 900 MI (Set of 6) Length 98 mm	0104A00156
135	EDT Sampling Cannulas - Peak (Set of 6)	0104A00049
136	EDT Disposable Filters (PKT of 100)	0104A00050
137	EDT Sinkers 5 mm ID x 15.5 mm OD (Set of 6)	0104A00094
138	EDT Disposable Filters (PKT of 100)	0104A00050
139	EDT Sinkers 6 mm ID x 18 mm OD (Set of 6)	0104A00095
140	EDT Sinkers 8 mm ID x 23 mm OD (Set of 6)	0104A00096
141	EDT Japanese Sinkers (Set of 6)	0104A00098
142	EDT Disposable Filters (PKT of 50)	0104A00698
143	EDT V-1234 Tygon Tube With SS Sleeve (Set of 6)	0104A00250
144	EDT Luer Lock BSFTL6-210-6 Female (Set of 8)	0104A00373
145	EDT-14L/14Lx – Tablet Dispenser Plate with Provision for Stationary Basket	0105A00403
146	EDT Sampling Cannula Fixed Height 250 MI (Set of 6) Length 65mm	0104A00435
147	EDT Clearview (230V) Version-2	0104A00432
148	EDT Kone Vessel 1000ml - Clear	0104A00010
149	EDT Basket 30 Mesh (Set of 6)	0104A00056
150	EDT-14L/14Lx PTFE Carrier Tube For Syringe Manifold (From Sampling Probe To Syringe Manifold) (Set Of 12)	0105A00322
151	Re-Coating For EDT/08LX/08L	0104A00718
152	Re-Coating For EDT/14LX/14L	0104A00715
153	EDT-14LX Lift Cover Set Top & Bottom.	0104A00719
154	EDT Deconex 12 basic (7ltr.)	0104A00723
155	EDT Rubber Holder For Disposable Filter (Set Of 6)	0104A00705
156	EDT Customized Japanese Sinkers 15 mm ID X 25-26 mm Height	0104A00245
157	EDT Rubber Holder For Disposable Filter (Set Of 6)	0104A00705
158	EDT Tygon Tube for Drainage	0104A00309
159	EDT Basket 40 Mesh Teflon Coated (Set of 6)	0104A00674
160	EDT Teflon Coated Basket 20 Mesh (Set of 6)	0104A00381
161	EDT Adjustable Sampling Cannula For Low Evaporation Jar LID (Set of 6)	0104A00047
162	Re-Coating For EDT-14L/14Lx Bath Top Plate – Aluminum	0104A00072
163	EDT Adjustable Sampling Cannulas For 900 ML (Set of 6)	0104A00257
164	EDT Laser Marked Glass 1000 MI - Amber	0104A00024
165	EDT Pin Type Sinkers - A (Set of 12)	0104A00671
166	EDT Deconex 12Basic (Pack of 10 Bottle each of 0.9L/1.3Kg)	0104A00737
167	EDT Tygon Tube for ETC-11	0104A00150
168	EDT-14L/14Lx Aluminum Top Plate for 150ml Conversion Kit	0105A00383
169	EDT Jar Lid for 100ml Vessels	0104A00249
170	EDT Cooling Instrument (Temp. From 10° to 40° C) - 3 litre	0104A00570
171	EDT Jar Lid for 250ml Vessels	0104A00389
172	EDT Jar Lock for Glass Vessel Suitable for Aluminum Top Plate (Set of 16)	0104A00385
173	EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big (Set of 6)	0104A00101
174	Offline Sampling Accessories for EDT-14Lx Dissolution Tester with Dx Sample Collector	0104A00711
175	EDT Laser Marked Glass Vessel 1000 MI - Amber With Sigma Delta Certificate	0104A00401
176	EDT Glass Vessel 1000ml - Amber	0104A00005
177	EDT Glass Vessel 1000ml - Clear	0104XA00004
178	EDT Laser Marked Glass Vessel 1000 MI - Clear	0104A00023
179	EDT Glass Vessel 1000ml - Clear	0104A00004
180	EDT Laser Marked Basket 40 Mesh	0104A00629
181	EDT Basket 40 Mesh	0104A00567
182	EDT Basket 40 Mesh for 150ml Vessel (Set of 6)	0104A00300
183	EDT Laser Marked Basket 80 Mesh (Set of 6)	0104A00077
184	EDT Basket Holder with 'O' Ring (Set of 6) for 150ml Vessel	0104A00299

185	EDT Digital Vernier Caliper for 1000 MI Vessel (150mm)	0104A00365
186	EDT Sampling Manifold Cannula S.S for 2Ltr (Set of 6)	0122A00019
187	EDT Validation Kit - Analog	0104A00118
188	EDT Basket Holder with 'O' Ring (Set of 6)	0104A00069
189	EDT Digital Dial Gauge	0104A00451
190	EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 8)	0104A00515
191	EDT Glass Vessel 100ml - Clear	0104A00001
192	EDT Felodipine Basket with Spacer (Set of 6)	0104A00064
193	EDT Laser Marked Basket 60 Mesh (Set of 6)	0104A00076
194	EDT Laser Marked Basket 100 Mesh (Set of 6)	0104A00078
195	EDT Sampling Cannula Fixed Height 500 MI (Set of 6) Length 108 mm	0104A00231
196	EDT Low Evaporation Jar LID (Set of 6)	0104A00042
197	Edt Hybrid Sampling Cannula For 2 Ltr.(Set Of 6)	0122A00011
198	EDT Tube Connecting Coupling With Nozzle (Set of 2)(Replacement Of QRC)	0104A00180
199	EDT-14L/LX Side Cir. Water Bath-Clear with Asmbly	0105A00307
200	EDT Legs for 2Ltr Water Bath (Set of 4)	0122A00021
201	EDT Transdermal Rod For Cylinder - Height 107.9 mm (Set of 8)	0104A00766
202	EDT Tube Connecting Coupling (Set of 2)(Replacement Of QRC)	0104A00181
203	EDT 33mm Watchglass Disc Assembly With 16 Mesh Screen	0104A00404
204	EDT 55mm Watchglass Disc Assembly With 16 Mesh Screen	0104A00405
205	EDT 33mm Watchglass Disc Assembly With 16 Mesh Screen	0104A00404
206	EDT 90mm Watchglass Disc Assembly With 16 Mesh Screen	0104A00406
207	EDT Kone Vessel 250ml - Clear	0104A00450
208	EDT Clearview (110V) Version-2.1	0104A00491
209	EDT 6U/12U Camera with USB Cable	0164A00003
210	EDT Glass Vessel 4Ltr - Amber	0133A00002
211	EDT Glass Vessel 2Ltr – Amber	0122A00018
212	EDT Laser Marked Glass Vessel 2 Ltr - Amber	0122A00017
213	EDT Magnetic Stand of Dial Gauge for 2 Ltr Vessel	0122A00016
214	EDT Digital Depth Gauge for 2 Ltr Vessels	0122A00005
215	EDT Glass Vessel 2Ltr - Clear	0122A00002
216	EDT Laser Marked Glass Vessel 2 Ltr Clear	0122A00001
217	EDT-14L/14Lx – Auto Dispenser Plate	0105A00454
218	EDT-14L/14Lx Bath Guard (Set of 2)	0105A00436
219	EDT-14L/14Lx Water Bath For Enhancer Cell with Complete Assembly (For ETC-11)	0105A00431
220	EDT-14L/14Lx Water Bath For Enhancer Cell with Complete Assembly (For ETC-15)	0105A00430
221	EDT-14L/14Lx – Auto Dispenser Plate	0105A00454
222	EDT-14Lx Top Cover	0105A00423
223	EDT-14L/14Lx – Water Circulation Silicon Tube Set	0105A00396
224	EDT-14L/14Lx 150 MI Low Volume Conversion Kit (Aluminum)	0105A00380
225	EDT-14L/14Lx PTFE Carrier Tube From Disso to Ismatec Pump (Set of 14)	0105A00338
226	EDT-14L/14Lx Blackout Curtain for I-Disso 12U	0105A00332
227	EDT-14L/14Lx -Motorised Sampling Manifold (For Peristaltic Offline)	0105A00326
228	EDT-14L/Lx -Immersion (Enhancer) Cell Kit 2cm ²	0105A00319
229	EDT-14L/14Lx S.S. Strip for Composite Top Plate (Set of 8)	0105A00313
230	EDT-14L/14Lx 150 MI Low Volume Conversion Kit (Composite)	0105A00304
231	EDT-14L/14Lx 150 MI Low Volume Conversion Kit (Composite)	0105A00304
232	EDT-14L/14Lx -Motorised Sampling Manifold (For Syringe Offline)	0105A00277
233	EDT-14L/14Lx Push Button for Tablet Dispenser Plate	0105A00270
234	EDT-14L/EDT-14Lx Composite Top Plate for 100ml Low Volume Conversion Kit	0105A00223
235	EDT-14L/EDT-14Lx Composite Top Plate for 250 MI Low Volume Conversion Kit	0105A00222
236	EDT-14L/14Lx Water Bath Clear with Quick Release Coupling & Nozzle Connector	0105A00161
237	EDT-14L/14Lx Composite Top Plate for Merlon Vessel	0105A00137
238	EDT-14L/14Lx Composite Top Plate for Merlon Vessel	0105A00137
239	EDT-14L/14Lx Composite Top Plate for Glass Vessel	0105A00136
240	EDT Jar Lock for Merlon Vessel Suitable for S.S. Top Plate (Set of 28)	0105A00133
241	EDT-14L/14Lx S.S. Top Plate for Merlon Vessel	0105A00131
242	EDT-14L/14Lx 1 Ltr. Conversion Kit for 250 MI Vessels	0104A00580
243	EDT-14L/14Lx 1 Ltr. Conversion Kit for 150 MI Vessels	0104A00578
244	EDT Cooling Instrument (Temp. From 10° to 40° C) - 6 litre	0104A00571
245	EDT Basket Holder with 'O' Ring	0104A00566
246	EDT 33mm Watch glass Disc.Assembly with 120 Mesh Screen	0104A00525
247	EDT 90mm Watch glass Disc.Assembly with 120 Mesh Screen	0104A00524

248	EDT 55mm Watch glass Disc.Assembly with 120 Mesh Screen	0104A00523
249	EDT Plastic 100ml Beakers For cleaning(Set of 8)	0104A00522
250	EDT-8/14LX Tablet Locater (set of 8)	0104A00520
251	EDT-8/14LX Sampling port seal (set of 8)	0104A00519
252	EDT-V 1234 Merlon Vessels set of 6 OLD type	0104A00506
253	EDT Customized Cage Type Sinkers Size 16mm ID x 13 mm HT With Spacer (Set of 6)	0104A00476
254	EDT Customized Cage Type Sinkers Size 12mm ID X 11mm Length (Set of 6)	0104A00475
255	EDT Replenishment Cannula (Set of 6)	0104A00473
256	EDT Laser Marked Kone Vessel 1000 ml - Clear With Sigma Delta Certificate	0104A00472
257	EDT Kone Vessel 1000 ml – Clear With Sigma Delta Certificate	0104A00471
258	EDT Customized Cage Type Sinkers Size 10mm ID Length x 10mm HT With Lid (Set of 6)	0104A00465
259	EDT Tablet Dispensing Derlin 'O' Ring for Dispenser plate (Set of 6)	0104A00464
260	EDT Manifold Rubber 'O' Ring for Dispenser plate (Set of 2)	0104A00462
261	EDT Customized Cage Type Sinkers Size 14mm ID x 13 mm HT With Spacer (Set of 6)	0104A00461
262	EDT Cage Type Sinkers - Small With Spacer (Set of 6)	0104A00460
263	EDT Cage Type Sinkers - Big With Spacer (Set of 6)	0104A00459
264	EDT Customized Granules Dispensing Adapter For Dispenser Plate (Set of 12)	0104A00458
265	EDT Spacer For Cage Type Sinkers Size 22ID X 22mm HT with lid (Set of 6)	0104A00457
266	EDT Spacer For Small Cage Type Sinkers (Set of 6)	0104A00456
267	EDT Spacer For Big Cage Type Sinkers (Set of 6)	0104A00455
268	EDT Deconex 12 Basic - 1000 ml	0104A00449
269	EDT 40 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)	0104A00447
270	EDT Customized Cage Type Sinkers – Medium With Lid (Set of 6)	0104A00446
271	EDT Teflon Coated Basket 10 Mesh (Set of 6)	0104A00444
272	EDT Basket LID 20 Mesh (Set Of 6)	0104A00441
273	EDT Bath Nozzle Rubber Gasket (Set of 2)	0104A00437
274	EDT Silicon Tube Set For ETC-15	0104A00433
275	EDT Basket Rod for 250ml Vessel	0104A00428
276	EDT Basket 40 Mesh for 250ml Vessel (Set of 6)	0104A00426
277	EDT Basket Sinkers 40 Mesh (9.00mm) with Lid (Set of 6)	0104A00425
278	EDT Disposable Filter 3.00mm (Pkts of 100)	0104A00424
279	EDT Big Retaining Ring for for Enhancer Cell Kit – 4cm (Set of 8)	0104A00415
280	EDT Top Cover with Base Plate for ETC-11LX	0104A00414
281	EDT Sampling Cannula Fixed Height 250 ML (Set of 6) Length 71 mm	0104A00413
282	EDT Luer Lock PMLP-21 (Male) (Set of 8)	0104A00412
283	EDT Silicon Tube for Drainage	0104A00411
284	EDT 120 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)	0104A00409
285	EDT 120 Mesh Screen For 55mm Watchglass Disc Assembly (Set of 6)	0104A00408
286	EDT 120 Mesh Screen For 33mm Watchglass Disc Assembly (Set of 6)	0104A00407
287	EDT Luer Lock FSLR-1 (Female) (Set of 8)	0104A00402
288	EDT Washer for Stationary Basket (Set of 2)	0104A00390
289	EDT Analog Pressure Gauge For KBR Press	0104A00388
290	EDT Analog Pressure Gauge For KBR Press	0104A00388
291	EDT Luer Lock FTLL-6 Female (Set of 8)	0104A00384
292	EDT Glass Vessel Trolley	0104A00379
293	EDT Customized Transparent Moulded Jar Lid - Plastic (Set of 6)	0104A00378
294	EDT Customized PTFE Carrier Tubing (1mm) with Ferrule Connection (Set Of 6)	0104A00377
295	EDT Sampling Cannula Fixed Height 500 ML (Set of 6) Length 127 mm	0104A00376
296	EDT Customized Helix Coil Sinkers 9.5mm ID X 20mm Height (Set of 6)	0104A00374
297	EDT Hardened Polished Surface Plate (for Intrinsic)	0104A00372
298	EDT Bolts for Fixing the Die and Surface Plate (for Intrinsic)(Set of 3)	0104A00371
299	EDT Stand For 1000ml Vessel (Composite)	0104A00370
300	EDT Deconex CIP Power-X (14 Kg)	0104A00369
301	EDT Online Connector (Set of 8)	0104A00368
302	EDT Glass Vessel 'D' Type 1000 ML - Amber With Ring	0104A00364
303	EDT 15ml Replenishment Bottle With Rubber Cap	0104A00351
304	EDT Luer Lock MLRL210-6 Male (Set of 8)	0104A00350
305	EDT Customized 10 Mesh Basket Sinkers with Lid (12.7mm ID – 24.5mm Height) (Set of 6)	0104A00348
306	EDT Customized 10 Mesh Basket Sinkers with Lid (9.1mm ID – 18.8mm Height) (Set of 6)	0104A00347
307	EDT Prednisone Tablet (Pkt of 30)	0104A00345
308	EDT Prednisone Standard USP-1559006 Reference 250 Mg	0104A00344
309	EDT Disc-Filter-Millex-HN, 33mm, 0.45um, Nylon Membrane, Non-Sterile (Pkt of 250)	0104A00343
310	EDT Disc Filter-Millex -HN, 33mm, 0.45um PES Membrane, Non-Sterile (Pkt of 250)	0104A00342

311	EDT Probe for Vibration Meter	0104A00336
312	EDT Sampling Cannula Fixed Height For Enhancer Cell Kit -200ml (Set of 6)	0104A00335
313	EDT Cylindrical Flat Bottom Glass Vessel 200ml – Clear for Enhancer Cell (Set of 6)	0104A00333
314	EDT Enhancer Cell for Enhancer Cell Kit – 2 CM (Set of 6)	0104A00327
315	EDT Sampling Manifold Cannula S.S. (SET Of 6)	0104A00324
316	EDT Luer Lock BSMLRL-6 Male (Set of 8)	0104A00321
317	EDT Enhancer Cell for Enhancer Cell Kit - 4 CM (Set of 6)	0104A00320
318	EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6)	0104A00316
319	EDT Sinkers 20 Mesh with Lid (Set of 6)	0104A00314
320	EDT Flat Bottom Glass Vessel 1000 MI - Clear for Stationary Intrinsic Attachment (Set of 6)	0104A00316
321	EDT Vernier Caliper for 150 ml Low Volume Conversion Kit	0104A00302
322	EDT Basket LID 40 Mesh for 150ml Vessel (Set of 6)	0104A00301
323	EDT Basket Rod for 150ml Vessel	0104A00298
324	EDT Jar Lid for 150ml Vessels	0104A00297
325	EDT Centering Device for 150 ml Low Volume Conversion Kit	0104A00296
326	EDT Sampling Cannula - Fixed Height for 150 MI Vessel (Set of 6)	0104A00295
327	EDT Sampling Cannula - Fixed Height for 150 MI Vessel (Set of 6)	0104A00295
328	EDT Laser Marked Glass Vessel 150 MI - Clear	0104A00292
329	EDT Glass Vessel 150 MI - Clear	0104A00291
330	EDT Basket Sinkers 40 Mesh (11.7mm Height) with Lid (Set of 6)	0104A00288
331	EDT Sinkers 10 Mesh with Lid (Set of 6)	0104A00284
332	EDT Basket 10 Mesh 24.5 mm (Set of 6) With Laser Marking & Certificate	0104A00279
333	EDT Deconex-CIP Fresh	0104A00276
334	EDT Deconex-CIP Alpha- X	0104A00275
335	EDT U Glass Tube	0104A00274
336	EDT Teflon Coated Basket Shaft With O Ring & Holder	0104A00273
337	EDT Rubber Ring For Low Evaporation Jar Lid	0104A00271
338	EDT Spirit Level	0104A00269
339	EDT Top Cover with Base Plate for ETC-11	0104A00268
340	EDT Basket 100 Mesh "V" Type	0104A00266
341	EDT Basket 40 Mesh "H" Type	0104A00264
342	EDT Basket 100 Mesh "H" Type	0104A00263
343	EDT S.S Knob for Bath Rod	0104A00261
344	EDT Merlon Vessel 1000ml Amber (N)	0104A00259
345	EDT SS Filters with Holders – 80 Micron (Set of 6)	0104A00255
346	EDT SS Filters with Holders – 40 Micron (Set of 6)	0104A00254
347	EDT SS Filters with Holders – 80 Micron (Set of 6)	0104A00255
348	EDT Rubber Holder For SS Filter	0104A00253
349	EDT Peak Glass Lense to Measure Basket Mesh	0104A00248
350	EDT Polished Punch	0104A00246
351	EDT Deconex 24 Liq (7 Kg)	0104A00241
352	EDT Deconex 25 Organ Acid (5 Kg)	0104A00240
353	EDT Vernier Caliper with Stopper for 1000ml Vessels	0104A00237
354	EDT Low Evaporation Jar Lid for Sampling Cannulas Fixed Height Suitable for Glass Vessel (Set of 6)	0104A00235
355	EDT Sampling Cannula Fixed Height 100 MI (Set of 6) Length 73.5 mm	0104A00232
356	EDT Dispenser Jar Ring Black (Set of 8)	0104A00229
357	EDT Merlon Vessel 1000 MI - Amber with Jar Lock Suitable for Composite Top Plate	0104A00228
358	EDT Merlon Vessel 1000 MI - Amber with Jar Lock Suitable for Composite Top Plate	0104A00228
359	EDT Vernier Caliper with Stopper for 1000ml Vessels	0104A00237
360	EDT Dispenser Jar Ring Black (Set of 8)	0104A00229
361	EDT Merlon Vessel 1000 MI - Amber with Jar Lock Suitable for Composite Top Plate	0104A00228
362	EDT Merlon Vessel 1000 MI - Amber with Jar Lock Suitable for S.S. Top Plate	0104A00227
363	EDT Merlon 1000 MI - Clear with Jar Centering Ring Suitable for Composite Top Plate	0104A00226
364	EDT Spindle Latch with Spring (Set of 8)	0104A00224
365	EDT Centering Device for 250 MI Low Volume Conversion Kit	0104A00223
366	EDT Japanese Sinkers with 40 Mesh Cover (Set of 6)	0104A00220
367	EDT Sinkers 40 Mesh with Lid (Set of 6)	0104A00219
368	EDT Customized Japanese Sinkers 15 mm ID X 40 mm Height (Set of 6)	0104A00218
369	EDT Holder for Basket 40 Mesh 24.5 mm (Set of 6)	0104A00206
370	EDT Deconex 12 Basic (500ml)	0104A00203
371	EDT Stop Watch	0104A00185
372	EDT Transdermal Cylinder - Height 50.8 mm(Set of 6)	0104A00184
373	EDT 'D' Type Basket 40 Mesh (Set of 6)	0104A00177

374	EDT 'D' Type Basket 10 Mesh (Set of 6)	0104A00174
375	EDT Luer Lock FTLL210-6 Female (Set of 8)	0104A00167
376	EDT Luer Lock MTLL220-6 Male (Set of 8)	0104A00166
377	EDT Gauge Block 25 mm	0104A00160
378	EDT Basket Validation Gauge with Certificate	0104A00154
379	EDT Moulded Tray For Basket, Basket Holder & Rod	0104A00152
380	EDT Plastic Quick Release Coupling with Nozzle (Set of 2)	0104A00145
381	EDT Vessel Centering Ring for Merlon Vessel Suitable for Composite Top Plate (Set of 8)	0104A00141
382	EDT Vessel Centering Ring for Glass Vessel Suitable for Composite Top Plate (Set of 8)	0104A00140
383	EDT Vessel Centering Plate for Merlon Vessel Suitable for S.S. Top Plate	0104A00135
384	EDT TNT Washer for Merlon Vessel Suitable for S.S.> Top Plate (Set of 16)	0104A00134
385	EDT Jar Lock for Glass Vessel Suitable for S.S. Top Plate (Set of 16)	0104A00132
386	EDT Digital Verticality Meter	0104A00125
387	EDT Thermometer - Mercury	0104A00123
388	EDT Digital Tachometer	0104A00122
389	EDT Depth Gauge for 1000 MI Vessel	0104A00121
390	EDT Centering Gauge	0104A00120
391	EDT Centering Device for 1000 MI Vessel	0104A00119
392	EDT Helix Coil Sinkers - CAPWHT-2S (Set of 6)	0104A00116
393	EDT Centering Device for 1000 MI Vessel	0104A00119
394	EDT Helix Coil Sinkers - CAPWHT-02 (Set of 6)	0104A00115
395	EDT Metformin Sinkers - Small (Set of 6)	0104A00114
396	EDT Metformin Sinkers - Big (Set of 6)	0104A00113
397	EDT Sinkers 'U' Type(Set of 6)	0104A00110
398	EDT Alfuzosin Sinkers Sinkers 37.5 mm ID x 20 mm OD (Set of 6)	0104A00107
399	EDT 3 Prong Sinkers - SS (Set of 6)	0104A00105
400	EDT 3 Prong Sinkers - Plastic (Set of 6)	0104A00104
401	EDT Spring Type Sinkers 9 mm ID x 30 mm OD - Small (Set of 6)	0104A00102
402	EDT Cage Type Sinkers - Small (Set of 6)	0104A00100
403	EDT Laser Marked Felodipine Basket (Set of 6)	0104A00081
404	EDT Laser Marked Stationary Basket (Set of 6)	0104A00080
405	EDT Laser Marked Basket 40 Mesh (Set of 6)	0104A00075
406	EDT Laser Marked Basket 10 Mesh (Set of 6)	0104A00072
407	EDT Basket 'O' Ring (Set of 6)	0104A00067
408	EDT Basket Rod	0104A00066
409	EDT Moulded Tray For Basket	0104A00065
410	EDT Stationary Basket (Set of 6)	0104A00062
411	EDT Basket 100 Mesh (Set of 6)	0104A00060
412	EDT Basket 80 Mesh (Set of 6)	0104A00059
413	EDT Basket 60 Mesh (Set of 6)	0104A00058
414	EDT Basket 20 Mesh (Set of 6)	0104A00055
415	EDT Basket 20 Mesh (Set of 6)	0104A00055
416	EDT Inline Filters - 0.45 Micron (Set of 6)	0104A00052
417	EDT SS Filters with Holders (Set of 6)	0104A00051
418	EDT Sampling Cannulas Fixed Height (Set of 6)	0104A00045
419	EDT Jar LID For Staggered Start - White (Set of 6)	0104A00044
420	EDT Moulded Jar LID - Plastic (Set of 6)	0104A00041
421	EDT Rubber Stand For Glass Vessel	0104A00040
422	EDT Laser Marked Glass Vessel With Ring 1000ml - Clear	0104A00037
423	EDT Laser Marked Glass Vessel 'V' Type 1000ml - Clear	0104A00035
424	EDT Laser Marked Glass Vessel 'D' Type 1000 MI - Clear	0104A00033
425	EDT Laser Marked Flat Bottom 1000ml - Amber	0104A00032
426	EDT Laser Marked Flat Bottom 1000ml - Amber	0104A00032
427	EDT Laser Marked Glass Vessel 250 MI - Amber	0104A00022
428	EDT Laser Marked Glass Vessel 250 MI - Amber	0104A00022
429	EDT Glass Vessel With Ring 1000ml - Clear	0104A00018
430	EDT Glass Vessel V Type 1000 ml – Amber	0104A00017
431	EDT Flat Bottom 1000 MI - Clear	0104A00012
432	EDT Glass Vessel 250ml - Amber	0104A00003
433	EDT Glass Vessel 250ml - Clear	0104A00002
434	EDT-14L/14Lx S.S. Top Plate for Glass Vessel	0103A00130
435	EDT Vibration Meter	0104A00127
436	EDT Spring Assembly for Tablet Dispenser Plate (Set of 2)	0104A00484

437	EDT Deconex 12 Basic (10 Bottles of 1000 ml)	0104A00442
438	EDT Spring Type Sinkers 9.5 mm ID x 30 mm OD - Big	0104A00623
439	EDT Felodipine Basket (Set of 6)	0104A00063
440	EDT 40 Mesh Stationary Basket (Set of 6)	0104A00584
441	EDT Basket 40 Mesh (Set of 6)	0104A00057
442	EDT-BTv Basket rod (1000 ml)	0169A00008
443	EDT Clearview (230V)	0104A00410
444	EDT-08/14Tx Manifold Motor Driv Card	0171B00006
445	EDT Basket 20 Mesh "V" Type	0104A00267
446	EDT 'D' Type Basket 20 Mesh (Set of 6)	0104A00175
447	EDT Basket 20 Mesh "H" Type	0104A00265
448	EDT Glass Vessel 500 ml - Clear	0104A00608
449	EDT 16 Mesh Screen For 90mm Watchglass Disc Assembly (Set of 6)	0104A00581
450	EDT- 406L Water Bath	0133A00005
451	Inspire 14 Top Cover	0104A00550
452	Connexion	0104A00785

	Breakdown Parts	
1	EDT-14Lx V-5 Individual Vessel Temp. Sensor PCB with Cable & Connector	0105BV4152
2	EDT 14L V-1234 Diving Manifold Bearing's Housing (Set of 2)	0105BV4148
3	EDT-14Lx V-4 – 10 Pin FRC Cable (Processor Card to Stirrer Opto Card)	0105BV4116
4	EDT-14Lx V-4 Front Panel With Keypad	0105BV4105
5	EDT-14L V-1234 Double Pully Belt Tensioner	0105B00087
6	EDT-14L V-1234 Belt Tightener	0105B00084
7	EDT-14L V-1234 Water Circulation Tygon Tube Set	0105B00058
8	EDT-14L V-1234 Big Rubber Leg for Water Bath 44mm (Set of 8)	0105B00056
9	EDT-14L V-1234 Small Rubber Leg for Water Bath 22mm (Set of 8)	0105B00052
10	EDT-14L/LX V-1234 Ind. Ves. Temp. Sensor(Set of 14)	0105B00012
11	EDT 14L V-1234 Manifold Bearing (Set of 2)	0105B00008
12	EDT-Lx V-5 Driver Card	0104BV5157
13	EDT D-Shell 9-Pin RS232/ Printer Port with 10 pin Frc Cable	0104BV4187
14	EDT-Lx V5 Heater Card for ETC-15 Without Trim Pot	0104BV4154
15	EDT-Lx V5 Compatch Card Without Trim Pot	0104BV4153
16	EDT-Lx V-45 Mains Supply Filter for ETC-15 (230V)	0104BV4141
17	EDT - Lx V-4 SMPS - 12V	0104BV4140
18	EDT-LX V-4 Mains Cord	0104BV4137
19	EDT-Lx V-4 Inlet Connector for Side Circulation Waterbath	0104BV4135
20	EDT-Lx V-4 EMIS Filter for ETC-15 (110V)	0104BV4132
21	EDT-Lx V-4 EMIS Filter for ETC-15 (230V)	0104BV4131
22	EDT-Lx V-4 Hi-Pot Filter for ETC-15	0104BV4111
23	EDT-Lx V-4 Temperature Controller ETC-15	0104BV4114
24	EDT-Lx V-4 RTC Component Kit	0104BV4094
25	EDT Lx V-4 Communication Cable (Syringe Pump to Sample Collector)	0104BV4091
26	EDT Lx V4 Communication Cable (Disso to Syringe Pump)	0104BV4090
27	EDT-Lx V-4 Serial Printer Cable	0104BV4087
28	EDT-Lx V-4 SMPS Card	0104BV4062
29	EDT-Lx V-4 BLDC Drive	0104BV4060
30	EDT-Lx V-4 Lift Opto Bottom	0104BV4048
31	EDT-Lx V4 Main Cable 15 AMP	0104BV4030
32	EDT V-1234 Bath Coupling Mounting Plate (Set of 2)	0104B00096
33	EDT V-1234 Inlet & Outlet Connector for Waterbath	0104B00087
34	EDT V-1234 Drain Rubber for Water Bath	0104B00050
35	EDT Dispensor Plate Holder (Set of 2)	0104B00048
36	EDT V-1234 Inline Fuse Holders	0104B00047
37	EDT V1234 Mains Switch (Yellow)	0104B00046
38	EDT V-1234 Rocker Switch Mains Red	0104B00045
39	EDT V-1234 Lift Motor 230V-2A	0104B00041
40	EDT V-1234 Lift Motor 110V-3A	0104B00037
41	EDT V-1234 Heater 110V-3A	0104B00035
42	EDT V-1234 Links Cable 5 Amp	0104B00032
43	EDT V-1234 Mains Fuse Holder	0104B00029
44	EDT V-1234 Thermal Cut-Off	0104B00028
45	EDT V-1234 Timing Pulley Set	0104B00026
46	EDT V-1234 Drainage Valve	0104B00025
47	EDT V-1234 Manifold Opto	0104B00009
48	EDT-08L/14L V-234 Proxy Card	0104B00004
49	EDT-08L/14L V-34 Manifold Bush Highspeed	0104B00003
50	EDT V -1234 LED with Cable (Set of 2)	0104BV4173
51	EDT-BTv Reader for Single/Double Pulley Spindle	0169B00008
52	EDT-Lx V-4 Heater- 230V-2A (For ETC-15) Bend Nozzle & Thermal Cutoff is Separate	0104BV4192
53	EDT-Lx V-4 Heater- 230V-2A (For ETC-15) Old Straight Nozzle	0104BV4191
54	EDT-Inspire BLDC Drive 8 Pole	0104BV4188
55	EDT V-1234 Cut spacer 0.2mm	0104BV4177
56	EDT V-1234 Cut Spacer 0.5mm	0104BV4176
57	EDT V-1234 Individual Vessel Temperature Sensor (Set of 8)	0104B00012
58	EDT V-1234 Single Pulley Spindle	0104B00014
59	EDT V-1234 Double Pulley Spindle	0104B00015

60	EDT-Lx V-4 Fuse Set	0104BV4027
61	EDT-8/14Lx V-4-Inspire Driver Card	0104BV4163
62	EDT-14L/14Lx Water Bath with Plastic QRC - Clear	0105A00143
63	EDT-Lx V-4 Heater- 110V (For ETC-15) Old Straight Nozzle	0104B00101
64	EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 230V	0104BV4128
65	EDT Wear & Tear Kit (For D)	0105BV4165
66	EDT-LX V4 Inspire Processor Card	0104BV4183
67	EDT-Lx V4 Compatch Card	0104BV4080
68	EDT-Lx V4 Gear for Lift Motor	0104BV4085
69	EDT-Lx V-45 Advance Filter for ETC-15 (230V)	0104BV4158
70	EDTV-1234 Heater 230V-2A	0104B00039
71	EDT- V-4-Inspire Stirrer opto card	0104BV4185
72	EDT-Lx V-4 Lift Disc Opto Card	0104BV4059
73	EDT V-1234 MF Connector for Manifold	0104B00083
74	EDT V-1234 External Probe	0104B00018
75	EDTV-1234 Bath Temperature Sensor	0104B00013
76	EDT-Lx V4 Stirrer Opto Card	0104BV4001
77	EDT V-1234 Individual Vessel Temperature Sensor	0104B00049
78	EDT-14Lx V-5 Individual Vessel Temperature Sensor PCB with Cable & Connector with out trim pot	0104BV5158
79	EDT- Wear & Tear Kit (For Dx)	0105BV4164
80	EDT V-1234 Height adjustment spacer 0.5mm	0104BV4178
81	EDT-14L V-1234 Timing Belt Set	0105B00022
82	EDT-Lx V-4 BLDC Motor with Gear	0104BV4118
83	EDT-Lx V-4 BLDC Motor with Gear And Drive(8 Pole)	0104BV4112
84	EDT-Lx V-4 LAN Adaptor	0104BV4055
85	EDT Lift lead screw Nut	0104BV4182
86	EDT-Lx V-4 Circulation Pump Kit (ETC-15) – 110V	0104BV4127
87	EDT V-1234 Circulation Pump 110V-3A	0104B00034
88	Inspire Processor card for 10 to 55 degC	0104BV4237
89	Inspire Lift Park Position Opto Sensor	0104BV4201
90	Inspire Program Patch	0104BV4238
91	Inspire 14 Diving Manifold Assembly without Temperature Probe and without Sampling Cannula	0104BV4200
92	Inspire Lift Park Position Opto Sensor	0104BV4201
93	Inspire -14 Serial Printer Cable	0104BV4204
94	Inspire-14 Platinum Kit Mechanical Spares	0104BV4227
95	Inspire-14 Platinum Kit Electronic Spares	0104BV4228