

Diffusion Cell Apparatus EDC-07



Features:

- USP/EP/IP Compliance
- 7 station apparatus
- Ideal for transdermal and topical products
- 10 programmable protocols with 12 sample intervals
- Ergonomically designed to save bench space
- Vertical diffusion cells available in 5, 7.5, 12, 15, and 20 ml sizes
- Specially designed for academic instruction
- Simultaneous timer function allows cumulative sampling intervals
- Equipped with seven diffusion cells
- 30–1000 RPM, settable
- External temperature sensor for validating media temperature

Specifications:

Model	EDC-07
Number of Cells	7
No. of protocol	10 maximum
Sample interval	12
No. of motor	1
Range	30.0°C to 40.0°C
Control accuracy	± 0.5°C
Resolution	0.1°C
Heater sensor	RTD (PT 100)
Heater wattage	100 watts
External probe sensor	RTD (PT 100)
Type	Magnetic Coupling
Flow rate	15 litres/min
No. of Timer	1
Range	1 min to 999 hrs 59 min
Display	20 x 4 LCD
Supply	230 VAC, 50 Hz OR 110 VAC, 60 Hz
Fuse	2 Amp/250 V (For Input 230 V AC,50Hz)
Power Consumption	1 K VA
Dimensions (L x W x H)	455 x 450 x 392 mm

Ordering information:

Products		
Product code	Description	Quantity
1020100	EDC -07-230V	1
1010100	EDC -07-110V	1

* The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Diffusion Cell Apparatus

Standard accessories		
Product code	Description	Quantity
1001A00060	EDC Diffusion Cell Without Purging – 7.5mL	1
1001A00071	EDC Thumb Nut (Set of 8)	1
1001A00017	EDC-07 Magnetic Candy 3×6 mm (Set of 7)	1
1001A00067	EDC Dummy Cell	1
1001A00069	EDC Acrylic Hood	1
0104B00018	EDC V-1234 External Probe	1
1001A00150	EDC Plastic Syringe	1
1001A00039	EDC-07 Sample Withdrawal Syringe – 1 ml (Set of 7)	1
1001B00005	EDC Temp calibration Assembly with probe holder	1
1001A00023	EDC-07 Clamp Plate (Set of 7)	1

Optional Accessories		
Product code	Description	Quantity
1001A00056	EDC Diffusion Cell Without Purging – 15mL (20mm ID × 34mm OD)	7
1001A00054	EDC Diffusion Cell Without Purging – 12.5mL (15mm ID × 34mm OD)	7
1001A00148	EDC Diffusion Cell Without Purging – 12mL (15mm ID × 34mm OD)	7
1001A00143	EDC Diffusion Cell Amber colour Without Purging – 12mL (15mm ID × 34mm OD)	7
1001A00052	EDC Diffusion Cell Without Purging – 7.5mL (15mm ID × 34mm OD)	7
1001A00050	EDC Diffusion Cell Without Purging – 5mL (9mm ID × 34mm OD)	7
1001A00012	EDC Glass Diffusion Cell Without Purging – 12ml – 34mm Orifice	7
1001A00109	EDC Glass Diffusion Cell Without Purging – 20ml – 34mm Orifice	7
1001A00004	EDC Glass Diffusion Cell Without Purging – 5ml – 34mm Orifice	7
1001A00037	EDC Glass Diffusion Cell Without Purging – 7.5ml – 34mm Orifice	7
1001A00112	EDC-07 Clamp Plate 20 ml (Set of 7)	7
1001A00068	EDC-07 Acrylic Cell Holder with SS Plate – 30 mm Orifice (Set of 7)	7
1001A00142	EDC-07 Sample Withdrawal Syringe -10ml (Set of 7)	1
1001A00023	EDC-07 Clamp Plate 7.5ml/12ml (Set of 7)	1
1001A00141	EDC Diffusion Cell Without Purging – 12ml (15mm ID X 30mm OD)	7
1001A00060	EDC Diffusion Cell Without Purging – 7.5mL (15mm ID X 30mm OD)	7
1001A00058	EDC Diffusion Cell Without Purging – 5mL (9mm ID X 30mm OD)	7
1101A00074	EDC-07 Top plate (Cell holder mounting plate)	1
1001A00010	EDC Glass Diffusion Cell Without Purging - 12ml – 30mm Orifice	7
1001A00006	EDC Glass Diffusion Cell Without Purging - 7.5ml – 30mm Orifice	7
1001A00136	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 30mm Orifice With Donar cap.	7
1001A00137	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 30mm Orifice With Donar cap.	7
1001A00040	EDC Glass Diffusion Cell Without Purging - 5 ML – 15mm Orifice (Old)	7
1001A00138	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 15mm Orifice With Donar cap.	7
1001A00139	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 15mm Orifice With Donar cap.	7
1001A00025	EDC-07 Donor Cap Close 7.5ml/12ml (Set of 7)	1
1001A00027	EDC-07 Tygon Tube (Set of 14)	1
1001A00026	EDC-07 Donor Cap Open 7.5 & 12 ML (Set of 7)	1
1001A00015	EDC-07 Holder Sealing Rubber (Set of 7)	1
1001A00044	EDC-07 Donor Cap Close 15ml (Set of 7)	1
1001A00024	EDC-07 Donor Washer 7.5ml/12ml (Set of 7)	1
1001A00140	EDC-07 Silicon Tube 4 x 7 (Set of 14) – 1	1
1001A00033	EDC-07 Glass Donor Cap Open 7.5 & 12 ML (Set of 7) – 1	1

* The information contained in this document is believed to be correct but **ELECTROLAB** accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Diffusion Cell Apparatus

Optional Accessories		
Product code	Description	Quantity
1001A00018	EDC-07 Magnetic Candy 3x10 mm (Set of 7) – 1	1
1001A00019	EDC-07 Clamp Plate 5ml (Set of 7) – 1	1
1001A00014	EDC-07 Acrylic Cell Holder with SS Plate – 34 mm Orifice (Set of 7) – 1	1
1001A00108	EDC-07 Acrylic Cell Holder – 34 mm Orifice (Set of 7) – 1	1
1001A00004	EDC-07 Certified Glass Cell (15 mm Orifice & 7.5ml Volume) Without Cap – 1	1
1001A00032	EDC-07 Acrylic Donor Cap Open 7.5 & 12 ML (Set of 7) – 1	1
1001A00012	EDC Glass Diffusion Cell Without Purging – 12ml – 34mm Orifice	1
1001A00109	EDC Glass Diffusion Cell Without Purging – 20ml – 34mm Orifice	1
1001A00004	EDC Glass Diffusion Cell Without Purging – 5ml – 34mm Orifice	1
1001A00037	EDC Glass Diffusion Cell Without Purging – 7.5ml – 34mm Orifice	1
1001A00112	EDC-07 Clamp Plate 20 ml (Set of 7)	1
1001A00068	EDC-07 Acrylic Cell Holder with SS Plate – 30 mm Orifice (Set of 7)	1
1001A00142	EDC-07 Sample Withdrawal Syringe -10ml (Set of 7)	1
1001A00023	EDC-07 Clamp Plate 7.5ml/12ml (Set of 7)	1
1001A00141	EDC Diffusion Cell Without Purging – 12ml (15mm ID X 30mm OD)	1
1001A00060	EDC Diffusion Cell Without Purging – 7.5mL (15mm ID X 30mm OD)	1
1001A00058	EDC Diffusion Cell Without Purging – 5mL (9mm ID X 30mm OD)	1
1101A00074	EDC-07 Top plate (Cell holder mounting plate)	1
1001A00010	EDC Glass Diffusion Cell Without Purging - 12ml – 30mm Orifice	1
1001A00006	EDC Glass Diffusion Cell Without Purging - 7.5ml – 30mm Orifice	1
1001A00136	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 30mm Orifice With Donar cap.	1
1001A00137	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 30mm Orifice With Donar cap.	1
1001A00040	EDC Glass Diffusion Cell Without Purging - 5 ML – 15mm Orifice (Old)	1
1001A00138	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 15mm Orifice With Donar cap.	1
1001A00139	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 15mm Orifice With Donar cap.	1
1001A00025	EDC-07 Donor Cap Close 7.5ml/12ml (Set of 7)	1
1001A00027	EDC-07 Tygon Tube (Set of 14)	1
1001A00026	EDC-07 Donor Cap Open 7.5 & 12 ML (Set of 7)	1
1001A00015	EDC-07 Holder Sealing Rubber (Set of 7)	1
1001A00044	EDC-07 Donor Cap Close 15ml (Set of 7)	1
1001A00024	EDC-07 Donor Washer 7.5ml/12ml (Set of 7)	1
1001A00140	EDC-07 Silicon Tube 4 x 7 (Set of 14)	1
1001A00033	EDC-07 Glass Donor Cap Open 7.5 & 12 ML (Set of 7) – 1	1
1001A00006	EDC Glass Diffusion Cell Without Purging - 7.5ml – 30mm Orifice	1
1001A00136	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 30mm Orifice With Donar cap.	1
1001A00137	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 30mm Orifice With Donar cap.	1
1001A00040	EDC Glass Diffusion Cell Without Purging - 5 ML – 15mm Orifice (Old)	1
1001A00138	EDC Glass Diffusion Cell Amber colour Without Purging - 7.5ml – 15mm Orifice With Donar cap.	1
1001A00139	EDC Glass Diffusion Cell Amber colour Without Purging - 12ml – 15mm Orifice With Donar cap.	1
1001A00025	EDC-07 Donor Cap Close 7.5ml/12ml (Set of 7)	1
1001A00027	EDC-07 Tygon Tube (Set of 14)	1
1001A00026	EDC-07 Donor Cap Open 7.5 & 12 ML (Set of 7)	1
1001A00015	EDC-07 Holder Sealing Rubber (Set of 7)	1
1001A00044	EDC-07 Donor Cap Close 15ml (Set of 7)	1
1001A00024	EDC-07 Donor Washer 7.5ml/12ml (Set of 7)	1
1001A00140	EDC-07 Silicon Tube 4 x 7 (Set of 14)	1
1001A00033	EDC-07 Glass Donor Cap Open 7.5 & 12 ML (Set of 7)	1

* The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Diffusion Cell Apparatus

Optional Accessories		
Product code	Description	Quantity
1001A00018	EDC-07 Magnetic Candy 3x10 mm (Set of 7)	1
1001A00019	EDC-07 Clamp Plate 5ml (Set of 7)	1
1001A00014	EDC-07 Acrylic Cell Holder with SS Plate – 34 mm Orifice (Set of 7)	1
1001A00108	EDC-07 Acrylic Cell Holder – 34 mm Orifice (Set of 7)	1
1001A00004	EDC-07 Certified Glass Cell (15 mm Orifice & 7.5ml Volume) Without Cap	1
1001A00032	EDC-07 Acrylic Donor Cap Open 7.5 & 12 ML (Set of 7)	1
1001A00031	EDC-07 Glass Donor Cap Close 7.5ml/12ml (Set of 7)	1
1001A00030	EDC-07 Acrylic Donor Cap Close 7.5ml/12ml (Set of 7)	1
1001A00029	EDC-07 Glass Donor Cap Open 5ml (Set of 7)	1
1001A00028	EDC-07 Donor Cap Close 5ml (Set of 7)	1
1001A00022	EDC-07 Donor Cap Open 5ml (Set of 7)	1
1001A00020	EDC-07 Donor Washer 5ml (Set of 7)	1
1001A00017	EDC-07 Magnetic Candy 3x6 mm (Set of 7)	1
1001A00016	EDC-07 Cell Cap (Set of 7)	1
1001A00003	EDC-07 Certified Glass Cell (15 mm Orifice & 12.5ml Volume) Without Cap	1
1001A00001	EDC Adjustable Leg with Check Nut (Set of 4)	1
1001A00009	EDC - 'O' Ring for Lid of 12ml Cell	1
1001A00044	EDC-07 Donor Cap Close 15ml (Set of 7)	1
1001A00065	EDC Teflon Donor Cap – Close 7.5ml (Set of 7)	1
1001A00066	EDC Teflon Donor Cap – Close 12.5ml (Set of 7)	1
1001A00070	EDC Teflon Donor Cap – Open 7.5ml (Set of 7)	1
1001A00151	EDC Diffusion Cell Without Purging – 5ml (15mm Orifice X 30mm OD)	1
1001A00008	EDC – Acrylic Lid of 12ml Cell	1
1001A00145	EDC Validation Kit	1
1001A00146	EDC Platinum Kit	1
1001A00075	EDC Diffusion Cell Without Purging – 10ml (34)	1
1001A00111	EDC Teflon Donor Cap – Close 20ml (Set of 7)	1
1001A00108	EDC-07 Acrylic Cell Holder – 34 mm Orifice (Set of 7)	1
1001A00005	EDC Glass Lid Of 12ml Cell	1

Assurance Plan		
Product code	Description	Quantity
1001S01	Platinum Assurance Plan for Diffusion Cell Apparatus Model EDC-07	1
1001S02	Gold Assurance Plan for Diffusion Cell Apparatus Model EDC-07	1
1001S03	Silver Assurance Plan for Diffusion Cell Apparatus Model EDC-07	1
1001S20	Bronze Assurance Plan for Diffusion Cell Apparatus Model EDC-07	1

* The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

+91 - 22 - 4041 3131 sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

+91 - 22 - 4161 3122 www.electrolabgroup.com

Diffusion Cell Apparatus

Breakdown Accessories		
Product code	Description	Quantity
1001B0035	EDC-07 Heater Card New	1
1001B0020	EDC-07 Circulation Pipe with Barb Assy.-1 & 2	1
1001B0029	EDC-07 Compactah Card	1
1001B0021	EDC-07 DCV Valve with Luer Lock FTLL 230-1	1
1001B0018	EDC-07 Driven Pulley Assembly (Set of 7)	1
1001B0027	EDC-07 Driver Card	1
1001B0030	EDC-07 Heater Card	1
1001B0002	EDC-07 Luer Lock Female FTLL 240-1 (Set of 7)	1
1001B0023	EDC-07 Luer Lock Male MTLL 240-1 (Set of 7)	1
1001B0028	EDC-07 OPTO Card	1
1001B0024	EDC-07 Processor Card with Patch	1
1001B0019	EDC-07 Timer Belt 17XXL (Set of 3)	1
1001B0034	EDC-07 Front Panel With Keypad	1
1001B0013	EDC – Circulation Pump	1
1001B0032	EDC Bath Temperature Probe	1
1001B0031	EDC Clamps for holding Heater (Set of 4)	1
1001B0001	EDC Heater Assembly	1
1001B0026	EDC Luer Lock FTLL 230-1 (Set of 2)	1
1001B0016	EDC Motor	1
1001B0005	EDC Luer Lock FTLL 230-1 (Set of 2)	1
1001B0033	EDC Temperature Probe	1

Service Charges		
Product code	Description	Quantity
1001S05	Service Charges for Diffusion Cell Apparatus Model EDC-07	1

IQ/OQ Document		
Product code	Description	Quantity
1001S15	IQ/OQ Document for Diffusion Cell Apparatus Model EDC-07	1

Reinstallation Charges		
Product code	Description	Quantity
1001S10	Reinstallation Charges without IQ/OQ Document for Diffusion Cell Apparatus Model EDC-07	1
1001S11	Reinstallation Charges with IQ/OQ Document for Diffusion Cell Apparatus Model EDC-07	1

Dismantling Charges		
Product code	Description	Quantity
1001S12	Inspection Charges for Diffusion Cell Apparatus Model EDC-07	1

Single Validation Charges		
Product code	Description	Quantity
1001S08	Single Validation Charges for Diffusion Cell Apparatus Model EDC-07	1

Inspection & Service Charges		
Product code	Description	Quantity
1001S08	Inspection & Service Charges for Diffusion Cell Apparatus Model EDC-07	1

In-House Repairing		
Product code	Description	Quantity
1001S07	In-House Repairing for Diffusion Cell Apparatus Model EDC-07	1

Upgradation		
Product code	Description	Quantity
1001S08	Upgradation for Diffusion Cell Apparatus Model EDC-07	1

Inspection Charges for Diffusion Cell Apparatus Model EDC-07		
Product code	Description	Quantity
1001S04	Inspection Charges for Diffusion Cell Apparatus Model EDC-07	1

* [The information contained in this document is believed to be correct but ELECTROLAB accepts no liability for any errors and reserves the right to alter specifications without notice.](#)

ELECTROLAB (India) PVT. LTD. (CIN NO : U33110MH2002PTC136735)

Head Office: 401, Tirupati Udyog, I. B. Patel Road, Off. Western Express Highway, Goregaon (East), Mumbai - 400 063, India.

☎ +91 - 22 - 4041 3131 ✉ sales@electrolabgroup.com

Factory: EL 23/24, T. T. C., Electronic Zone, M. I. D. C, Mahape, Navi Mumbai - 400 710, India.

☎ +91 - 22 - 4161 3122 🌐 www.electrolabgroup.com